biotical Entamoeba Card

Rapid Test to detect
Entamoeba in card format

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biotical health
ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE:
QUALITY MANAGEMENT SYSTEM

biotical

Entamoeba card

ENGLISH

WHAT IS BIOTICAL ENTAMOEBA CARD USED FOR?

biotical Entamoeba card is a rapid immunochromatography test for the qualitative detection of *Entamoeba histolytica* and *Entamoeba dispar* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional can make a presumptive diagnosis of amebiasis.

WHAT IS AMEBIASIS?

Amebiasis is an infection of the gastrointestinal tract due to *Entamoeba histolytica*, a protozoan capable of invading the intestinal mucosa and can spread to other organs. *Entamoeba dispar* (amoeba morphologically similar to *E. histolytica*), is the species without invasive potential and is the cause of most asymptomatic infections.

Amebiasis uses to occur in areas where sanitary conditions are inadequate. On the entire planet, each year about 40 to 50 million people develop amebiasis, and between 40,000 and 70,000 of these people die as a result of the disease, making amebiasis the fourth leading cause of death due to protozoan infection.

Entamoeba exists in two forms: trophozoite (active parasite) or cyst (latent parasite). Many people have no symptoms or have few symptoms; however, they excrete cysts in stool and, therefore, they can spread the infection.

The disease can manifest itself as acute, chronic or as an asymptomatic infection.

Possible symptoms of amebiasis include flatulence, painful abdominal cramps, intermittent diarrhea or constipation and, in most severe cases, dysentery (loose stools with mucus and blood), formation of amebomas (large lumps in the colon) and abscesses in the liver. The infection occurs in contaminated water, fresh vegetables and fruits or other undercooked or washed food. It is also possible that feces cysts are transported to food through flies or cockroaches.

The diagnosis is usually analyzing the stool and, if necessary, performing a colonoscopy, ultrasound or blood test.

biotical Entamoeba card is a qualitative immunochromatographic test for the detection of and Entamoeba histolytic and dispar antigens in stool samples. The result is obtained in 10 minutes.

HOW DOES BIOTICAL ENTAMOEBA CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- **In the test line (T):** with mouse monoclonal antibodies against *Entamoeba*.
- **In the control line (C):** with rabbit polyclonal antibodies against a specific protein.

Sample area (S): the pad where the sample is placed has been dispensed with a test solution (mouse monoclonal antibodies against Entamoeba) conjugated to red polystyrene latex and a control solution (specific binding protein) conjugated to green polystyrene latex, forming two conjugate colored complexes.

TEST LINE (T):

- **If the sample is positive:** The antigens of the patient's sample will react with the red colored conjugate complex (anti-Entamoeba monoclonal antibodies - red latex microspheres) and advance by capillarity through the membrane. The anti-Entamoeba antibodies present in the membrane of the test line (T) will capture the colored test complex resulting in the red line.
- **If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red colored complex. The anti-Entamoeba antibodies present in the membrane of the test line (T) will not capture the red colored antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the *biotical Entamoeba card*. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the colored test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- *biotical Entamoeba card* test
- Instructions for use
- Collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

VERY IMPORTANT: Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.


TEST PROCEDURE

Allow test, samples and controls to reach room temperature (15-30 °C) prior to testing. Do not open pouches until the performance of the assay.


- Proceed to shake the stool collection tube in order to assure good sample dispersion.
- Remove *biotal Entamoeba card* test from its sealed bag just before using it.
- Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid
- **Read the results at 10 minutes. Do not read the test result later than 10 minutes.**

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(please refer to the illustration below)

	<i>biological test</i>	Interpretation of the results
1	C = GREEN T = Ø	There is no <i>Entamoeba histolytica</i> or <i>E. dispar</i> presence. No infection caused by <i>Entamoeba histolytica</i> or <i>E. dispar</i> .
2	C = GREEN T = RED	There is <i>Entamoeba histolytica</i> or <i>E. dispar</i> presence. Infection caused by <i>Entamoeba histolytica</i> or <i>E. dispar</i> .
3-4	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

- **NEGATIVE:** Only the GREEN line appears across the results window in the site marked with the letter C (control line).
- **POSITIVE:** In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.
- **INVALID:** Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test should be performed no later than 2 hours after opening the bag containing it.
- The sample must be recent or frozen without fixatives or preservatives. The samples must be taken during the first week of appearance of the symptoms to ensure a greatest presence of pathogens in the sample.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.

biotical Entamoeba card

Rapid test to detect *Entamoeba* in card format

- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Entamoeba histolytica* and/or *Entamoeba dispar* in the sample. To confirm the presence of infection by *Entamoeba histolytica* or *Entamoeba dispar* it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques (biochemical or microscopy methods).
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. A negative result does not indicate absence of infection by *Entamoeba histolytica* or *Entamoeba dispar*. It is possible to appear false negatives by not reaching the detection limit of the test.

EXPECTED VALUES

Entamoeba histolytica is the third main cause of morbidity and mortality due to parasitic disease in humans after malaria and schistosomiasis. It is estimated that *E. histolytica* and *E. dispar* parasitize approximately 10% of the world population, of which 90% are asymptomatic infections. About between 50000 and 100000 deaths every year are due to this disease.

In developed countries, amoebiasis tends to be more common in older and immunodepressed individuals. However, in tropical regions, the epidemiology of amoebiasis is completely different and is more common among the general population and particularly among patients attending health care centres with diarrhoea.

Out of 10% of the world's population infected by *E. histolytica*, only 1% develops invasive form of the disease. Only approximately 10% (pathogenic infections) need to be treated. In developed countries, *E. histolytica* infection is not endemic and occurs mostly after traveling to areas of endemicity.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation was performed, using faecal samples, comparing the results obtained by an immunochromatographic test (biotical Entamoeba, biotical) and other commercial qPCR tests (VIASURE *Entamoeba histolytica* Real Time Detection Kit and VIASURE *Entamoeba dispar* Real Time Detection Kit, biotical).

The results were as follows:

		qPCR test: VIASURE <i>Entamoeba histolytica</i> Real Time Detection Kit		
		+	-	Total
IC test: <i>biotical Entamoeba card</i>	+	4	1	5
	-	1	107	108
	Total	5	108	113

	Mean Value <i>E. histolytica</i>	95% confidence interval
Sensitivity	80%	28.4 - 99.5%
Specificity	99%	94.9 - 100%
PPV	80%	28.4 - 99.5%
NPV	97%	94.9 - 100%

		qPCR test: VIASURE <i>Entamoeba dispar</i> Real Time Detection Kit		
		+	-	Total
IC test: <i>biotical Entamoeba card</i>	+	15	1	16
	-	6	107	113
	Total	21	108	129

	Mean Value <i>E. dispar</i>	95% confidence interval
Sensitivity	71%	47.8 - 88.7%
Specificity	99%	94.9 - 100%
PPV	94%	69.8 - 99.8%
NPV	95%	88.8 - 98.0%

		qPCR test: VIASURE <i>Entamoeba histolytica</i> Real Time Detection Kit and VIASURE <i>Entamoeba dispar</i> Real Time Detection Kit		
		+	-	Total
IC test: <i>biotical Entamoeba card</i>	+	19	2	21
	-	7	214	221
	Total	26	216	242

	Mean Value <i>E. histolytica</i> + <i>E. dispar</i>	95% confidence interval
Sensitivity	73%	52.2 - 88.4%
Specificity	99%	96.7 - 100%
PPV	90%	69.6 - 98.8%
NPV	97%	93.6 - 98.7%

The results showed a high sensitivity and specificity to detect *Entamoeba histolytica* and *E. dispar* using biotical Entamoeba

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Entamoeba; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Campylobacter jejuni</i>	<i>Helicobacter pylori</i>	<i>Shigella boydii</i>
<i>Campylobacter coli</i>	<i>Listeria monocytogenes</i>	<i>Shigella dysenteriae</i>
<i>Clostridium difficile</i>	<i>Salmonella enteritidis</i>	<i>Shigella flexneri</i>
<i>Escherichia coli</i> O157:H7	<i>Salmonella paratyphi</i>	<i>Shigella sonnei</i>
<i>Cryptosporidium parvum</i>	<i>Salmonella typhi</i>	<i>Staphylococcus aureus</i>
<i>Giardia lamblia</i>	<i>Salmonella typhimurium</i>	

REFERENCES

1. BLESSMANN J., et al. "Longitudinal Study of Intestinal *Entamoeba histolytica* Infections in Asymptomatic Adult Carriers". *Journal of clinical microbiology*, Oct. 2003, p. 4745–4750 Vol. 41, No. 10.
2. GONIN P., et al. "Detection and Differentiation of *Entamoeba histolytica* and *Entamoeba dispar* Isolates in Clinical Samples by PCR and Enzyme-Linked Immunosorbent Assay". *Journal of clinical microbiology*, Jan. 2003, p. 237–241 Vol. 41, No. 1

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