biotical Giardia Card

Rapid Test to detect
Giardia in card format

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biotical health
ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE:
QUALITY MANAGEMENT SYSTEM

biotical Giardia card

ENGLISH

WHAT IS BIOTICAL GIARDIA CARD?

biotical Giardia card is a colored chromatographic immunoassay test for the qualitative detection of *Giardia* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of giardiasis.

WHAT IS GIARDIASIS?

Giardia is one of the most frequent parasites that cause water-borne diseases in humans: cryptosporidiosis and giardiasis.

Giardiasis is an intestinal disease caused by *Giardia intestinalis* (also called *Giardia lamblia* and *Giardia duodenalis*), a single-celled microscopic parasite that is found in the intestine of people and animals and is transmitted through the feces of an infected person or animal.

This parasite is protected by an outer covering (cyst) that allows them to survive outside the body and in the environment for long periods of time. Its infective form is usually by fecal-oral transmission. Among the most common symptoms of both diseases are: watery diarrhea, cramps and stomach disorders. In immunocompromised people or in elderly people, complications can appear that lead to serious illness and even death.

The reference method for the diagnosis of Giardiasis diseases is the identification of cysts in an examination with optical microscopy. With this complex technique, non-immediate results are obtained. Considering that Giardiasis has treatment, it is important to do an immediate diagnosis.

biotical Giardia card, is a qualitative immunochromatographic test for the detection of *Giardia* antigens in stool samples. The result is obtained in 10 minutes.

HOW DOES BIOTICAL GIARDIA CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- **In the test line (T):** with mouse monoclonal antibodies against *Giardia*.
- **In the control line (C):** with rabbit polyclonal antibodies against a specific protein.

Sample area (S): the pad where the sample is placed has been dispensed with a test solution (mouse monoclonal antibodies against *Giardia* for the T line) conjugated to red polystyrene latex and a control solution (specific binding protein) conjugated to green polystyrene latex, forming two conjugate colored complexes.

TEST LINE (T):

- **If the sample is positive:** The antigens of the patient's sample will react with the red colored conjugate complex (anti-*Giardia* monoclonal antibodies - red latex microspheres) and advance by capillarity through the membrane. The anti-*Giardia* antibodies present in the membrane of the T-line will capture the colored test complex resulting in the red line.
- **If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red colored complex. The anti-*Giardia* antibodies present in the membrane of the T-line will not capture the red colored antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the *biotical Giardia card*. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the colored test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- *biotical Giardia card* test
- Instructions for use
- Collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

VERY IMPORTANT: Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.


TEST PROCEDURE

Allow test, samples and controls to reach room temperature (15-30 °C) prior to testing. Do not open pouches until the performance of the assay.


- Proceed to shake the stool collection tube in order to assure good sample dispersion.
- Remove *biotal Giardia card* test from its sealed bag just before using it.
- Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid
- **Read the results at 10 minutes. Do not read the test result later than 10 minutes.**

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(please refer to the illustration below)

	<i>biotical test</i>	Interpretation of the results
1	C = GREEN T = Ø	There is no <i>Giardia</i> presence. No infection caused by <i>Giardia</i> .
2	C = GREEN T = RED	There is <i>Giardia</i> presence. Infection caused by <i>Giardia</i> .
3-4	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

- **NEGATIVE:** Only the GREEN line appears across the results window in the site marked with the letter C (control line).
- **POSITIVE:** In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.
- **INVALID:** Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test should be performed no later than 2 hours after opening the bag containing it.
- The sample must be recent or frozen without fixatives or preservatives. The samples must be taken during the first week of appearance of the symptoms to ensure the greatest presence of pathogens in the sample.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.

biotical Giardia card

Rapid test to detect *Giardia* in card format

- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Giardia* in the sample. To confirm the presence of infection by *Giardia* it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques (biochemical methods or microscopy).
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory for confirmation. A negative result does not indicate absence of infection by *Giardia*. It is possible to appear false negatives by not reaching the detection limit of the test.

EXPECTED VALUES

Giardia infection has become recognized as one of the most common causes of waterborne disease (found in both drinking and recreational water) in humans. *Giardia* is found worldwide.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using an immunochromatographic test (biotical *Giardia*, biotical) and these results were confirmed using a microscopy technique and PCR (positive results).

The results were as follows:

		Microscopy technique/PCR		
		+	-	Total
IC test: <i>biotical Giardia card</i>	+	61	0	61
	-	2	191	193
	Total	63	191	254

	Sensitivity	Specificity	PPV	NPV
<i>biotical Giardia</i> vs Microscopy technique	97%	> 99%	> 99%	99%

The results showed a high sensitivity and specificity to detect *Giardia* using *biotical Giardia card*.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical *Giardia*; no cross reactivity against gastrointestinal organism occasionally present in faeces:

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