biotical
H. pylori
Card

*Rapid Test to detect
Helicobacter pylori in card
format*

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EN - ES

Rev. 01
04-2021

ISO 13485 CERTIFIED
MEDICAL DEVICE:
QUALITY MANAGEMENT SYSTEM

biotical *H. pylori* card

ENGLISH

WHAT IS BIOTICAL *H. PYLORI* CARD USED FOR?

biotical *H. pylori* card is a rapid immunochromatography test for the qualitative detection of *Helicobacter pylori* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional can make a presumptive diagnosis to detect *Helicobacter pylori* on the patient.

WHAT IS HELICOBACTER PYLORI (H.PYLORI)?

Helicobacter pylori (HP) is a spiral-shaped bacterium that has bacterial flagella at one end, which gives it the mobility to colonize the stomach. It is one of the most prevalent human infections, affecting approximately 60% of the population, although in populations with lower socioeconomic status and worse health conditions, it reaches 80%.

The infection can develop in an inflammatory response in the stomach and chronic gastritis. Although only a small percentage of patients will develop symptoms despite suffering from gastritis, it seems to depend on both host factors and the strain of HP which causes the infection.

In addition to causing chronic gastritis, this bacterium is related to gastric and duodenal ulcer, lymphoma and gastric cancer. *H. pylori* can be transmitted from one person to another by direct contact with saliva, vomit or stool. It can also be transmitted through contaminated food or water.

biotical *H. pylori* card is a qualitative immunochromatographic test for the detection of HP antigens in stool samples. The result is obtained in 10 minutes.

HOW DOES BIOTICAL HELICOBACTER PYLORI CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- **In the test line (T):** with mouse monoclonal antibodies against *H. Pylori*.
- **In the control line (C):** with rabbit polyclonal antibodies against a specific protein.
- **Sample area (S):** the pad where the sample is placed has been dispensed with a test solution (mouse monoclonal antibodies against *Helicobacter pylori*) conjugated to red polystyrene latex and a control solution (specific binding protein) conjugated to red polystyrene latex, forming two conjugate colored complexe.

TEST LINE (T):

- **If the sample is positive:** The antigens of the patient's sample will react with the red colored conjugate complex (anti-*Helicobacter pylori* monoclonal antibodies-red latex microspheres) and advance by capillarity through the membrane. The anti-*Helicobacter pylori* antibodies present in the membrane of the test line (T) will capture the colored test complex resulting in the red line.
- **If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red colored complex. The anti-*Helicobacter pylori* antibodies present in the membrane of the test line (T) will not capture the red colored antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

- **Regardless of whether the sample is positive or not,** to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the red control line will always appear.

EXTERNAL POSITIVE CONTROL

Packs of 30 units of positive control swabs are available to the customer.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C (35,6-86°F). The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

biotical *H. pylori* Card

Rapid Test to detect *Helicobacter pylori* in card format

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CONDITIONS OF TRANSPORT

The preservation of the kit has been evaluated under normal storage conditions (24 months at room temperature), accelerated stability (5 days at 45°C or 113°F) and in a climatic chamber to simulate transport conditions (18 cycles from 2 to 40°C or 35,6-104°F, with 90% relative humidity), keeping the characteristics of the kit unchanged.

PRECAUTIONS

- For professional in vitro diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.

- Used tests should be managed as sanitary waste (sanitary waste container).
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the biotical *Helicobacter pylori* test. Do not use any other commercial kit component, if compatibility is not specified.

- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

- Pale yellow or orange bands may appear in the test strip results window, at the positions where the colored test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- biotical *H. pylori* card test
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

The kit contains materials for analysis of up to 25 samples.


SPECIMEN COLLECTION AND PREPARATION

Stool samples.

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8 °C or 35,6-46,4°F) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C (-4°F). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- Take out the cap of the collection tube (1) and use the stick to pick up sufficient stool sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2) making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample (approx. 50mg) and add it to the collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

TEST PROCEDURE

Allow test, samples and controls to reach room temperature (15-30°C or 59-86°F) prior to testing. do not open pouches until the performance of the assay.

- Proceed to shake the stool collection tube in order to assure good sample dispersion.
- Remove the biotical *H. pylori* card test from its sealed bag just before using it.
- Take the stool collection tube, cut the end of the cap (4) and dispense 3 drops in the circular window marked with the letter S (5), avoid adding solid particles.
- **Read the results at 10 minutes. Do not read the test result later than 10 minutes.**

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

biotal *H. pylori* Card


Rapid Test to detect *Helicobacter pylori* in card format

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
INTERPRETATION OF RESULTS


(please refer to the illustration below)


1 - NEGATIVE

2 - POSITIVE

3 - NO VALID RESULT

4 - NO VALID RESULT

	biotal test	Interpretation of the results
1	C = RED T = NO LINE	There is no <i>Helicobacter pylori</i> presence. No infection caused by <i>Helicobacter pylori</i> .
2	C = RED T = RED	There is <i>Helicobacter pylori</i> presence. <i>Helicobacter pylori</i> infection, which might mean gastrointestinal diseases (stomach and duodenum) like gastritis, peptic ulcer disease or gastric carcinoma.
3-4	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

- **NEGATIVE:** Only the RED line appears across the results window in the site marked with the letter C (control line).
- **POSITIVE:** In addition to the RED control line, a RED line also appears in the site marked with the letter T (test line) in the results window.
- **NO VALID:** Total absence of the control coloured line (RED) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A red line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Helicobacter pylori* in the sample. To confirm the presence of infection by *Helicobacter pylori* it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other invasive techniques (endoscopy).
- False positives may occur in the test with bloody samples or mucous stools, so the use of the test in this class of samples is not recommended.
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory or invasive techniques for confirmation. A negative result does not indicate absence of infection by *Helicobacter pylori*. It is possible to appear false negatives by not reaching the detection limit of the test.

biotical *H. pylori* Card

Rapid Test to detect *Helicobacter pylori* in card format

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EXPECTED VALUES

Studies have found that more than 90% of patients with duodenal ulcer and 80% of patients with gastric ulcer are infected with *Helicobacter pylori*. biotical *H. pylori* has been compared with different methods: cultures, Urea Breath Test and Urease Test, demonstrating an overall accuracy of >92%.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit):

Detection limit range is: 0.78ng/mL – 0.09ng/mL of *H. pylori* recombinant outer membrane protein.

Clinical sensitivity and specificity:

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical *H. pylori*, biotical) and qPCR technique (VIASURE *Helicobacter pylori* Real Time Detection kit, CerTest). The samples were taken from patients with the same as *Helicobacter pylori* infection symptoms arrived in Gastroenterology area in a Spanish hospital.

The results were as follows:

		qPCR test: Viasure <i>Helicobacter pylori</i> Real Time Detection kit		
		+	-	Total
IC test: biotical <i>H. pylori</i> card	+	54	1	55
	-	1	60	61
	Total	55	61	116

	biotical <i>H. pylori</i> vs Viasure <i>Helicobacter pylori</i> Real Time Detection Kit	
	Mean Value	95% confidence interval
Sensitivity	98,2%	90,3 – 100,0%
Specificity	98,4%	91,2 – 100,0%
PPV	98,2%	90,3 – 100,0%
NPV	98,4%	91,2 – 100,0%
Accuracy	98,3%	93,9 – 99,8%

The results showed a high sensitivity and specificity to detect *Helicobacter pylori* using biotical *H. pylori* card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical *H. pylori*; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Campylobacter coli</i>	<i>Salmonella enteritidis</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	<i>Salmonella paratyphi</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i>	<i>Salmonella typhi</i>	<i>CShigella sonnei</i>
<i>Escherichia coli O157:H7</i>	<i>Salmonella typhimurium</i>	<i>Staphylococcus aureus</i>
<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>	<i>Yersinia enterocolitica</i>

REFERENCES

- 1 Bruce E. Dunn, Hartley Cohen & Martin J. Blaser. *Helicobacter pylori*. Clin. Microbiol. Rev. 10 (4), 720-741, Oct. (1997).
- 2 Martin J. Blaser. *Helicobacter pylori* and gastric diseases. BMJ; 316: 1507-1510 (1998).
- 3 John L. Telford, Antonello Covacci, Rino Rappuoli & Paolo Ghiara. Immunobiology of *Helicobacter pylori* infections. Current Opinion in Immunology, 9; 498-503 (1997).

biotical H. pylori Card

Rapid Test to detect Helicobacter pylori in card format

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