biological Rotavirus Card

Rapid Test to detect
Rotavirus in card format

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biological health
ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE:
QUALITY MANAGEMENT SYSTEM

biotical Rotavirus card

ENGLISH

WHAT IS BIOTICAL ROTAVIRUS CARD?

biotical Rotavirus card is a colored chromatographic immunoassay test for the qualitative detection of *Rotavirus* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of *Rotavirus* infection.

GASTROENTERITIS AND ROTAVIRUS

Rotavirus and *Adenovirus* are the main cause of severe diarrhea in children under 5 years, causing high morbidity and mortality worldwide.

Acute gastroenteritis is a very common disease that affects children and young people around the world. Enteric *Rotavirus* is the first cause of severe diarrhea in children under 5 years of age, followed by *Adenovirus*.

Rotavirus is transmitted by feco-oral contact, causing watery diarrhea, vomiting, headache, fever and stomachache. Symptoms usually begin 1 or 2 days after becoming infected and may last 1 to 3-4 days.

The most common procedures to identify a gastroenteritis caused by *Rotavirus* involve the isolation and identification of the virus through techniques that require more than 1 day. However, *biotical Rotavirus card*, is an immunochromatographic qualitative test, easy to use and with almost immediate results, for the detection of *Rotavirus* directly from the patient's stool sample, providing results in 10 minutes.

HOW DOES BIOTICAL ROTAVIRUS CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- **In the test line (T):** with mouse monoclonal antibodies against *Rotavirus*.
- **In the control line (C):** with rabbit polyclonal antibodies against a specific protein.

Sample area (S): the pad where the sample is placed has been dispensed with a preparation of reagents of the T line (mouse monoclonal antibodies against *Rotavirus*) conjugated to red polystyrene latex; and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex, forming two conjugate colored complexes.

TEST LINE (T):

- **If the sample is positive:** The antigens of the patient's sample will react with the red colored conjugate complex (anti-*Rotavirus* monoclonal antibodies - red latex microspheres for the T line) and advance by capillarity through the membrane. The anti-*Rotavirus* antibodies present in the T line will capture the colored test complex resulting in the red line.
- **If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red colored complex. The anti-*Rotavirus* antibodies present in the T line will not capture the red colored antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30 °C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the *biotical Rotavirus test*. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the colored test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- *biotical Rotavirus card* test
- Instructions for use
- Collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20 °C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

biotal Rotavirus card


Rapid test to detect *Rotavirus* in card format

TEST PROCEDURE

Allow test, samples and controls to reach room temperature (15-30 °C) prior to testing. Do not open pouches until the performance of the assay.


- Proceed to shake the stool collection tube in order to assure good sample dispersion.
- Remove *biotal Rotavirus card* test from its sealed bag just before using it.
- Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid
- **Read the results at 10 minutes. Do not read the test result later than 10 minutes.**

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(please refer to the illustration below)

	<i>biotal</i> test	Interpretation of the results
1	C = GREEN T = Ø	There is no <i>Rotavirus</i> presence. No viral infection caused by <i>Rotavirus</i> .
2	C = GREEN T = RED	There is <i>Rotavirus</i> presence. Viral infection caused by <i>Rotavirus</i> .
3	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

- **NEGATIVE:** Only the GREEN line appears across the results window in the site marked with the letter C (control line).
- **POSITIVE:** In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.
- **INVALID:** Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test must be performed within 2 hours after opening the sealed pouch.
- Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
- *biotical Rotavirus* card will only indicate the presence of *Rotavirus* in the specimen (qualitative detection) and should be used for the detection of *Rotavirus* antigens in human faeces specimens only. Neither the quantitative value nor the rate of increase in *Rotavirus* antigens concentration can be determined by this test.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- Samples should be collected as soon as possible once the symptoms appear. After one week of infection the number of viruses decrease, making the sample less reactive.
- Some stool samples can decrease the intensity of the control line.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods are recommended. A negative result does not at any time preclude the possibility of *Rotavirus* infection s it may be due to a very low antigen concentration.
- This test provides a presumptive diagnosis of *Rotavirus* infections. All results must be interpreted together with other clinical information and laboratory findings available to the physician.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using an immunochromatographic test (*biotal Rotavirus, biotal*) and these results were confirmed with a commercial available test (Ridascreen® *Rotavirus* ELISA Test, r-Biopharm).

The results were as follows:

		Ridascreen® <i>Rotavirus</i> ELISA Test		
		+	-	Total
IC test: <i>biotal Rotavirus</i> card	+	18	1	19
	-	0	43	43
	Total	18	44	62

	Sensitivity	Specificity	PPV	NPV
<i>biotal Rotavirus</i> vs Ridascreen® <i>Rotavirus</i> ELISA Test	>99%	98%	>94%	>99%

The results showed a high sensitivity and specificity to detect *Rotavirus* using *biotal Rotavirus*.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of *biotal Rotavirus*; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

- | | | |
|-------------------------------|---------------------------------|--------------------------------|
| <i>Adenovirus</i> | <i>Escherichia coli</i> O157:H7 | <i>Salmonella typhi</i> |
| <i>Astrovirus</i> | <i>Giardia lamblia</i> | <i>Shigella boydii</i> |
| <i>Campylobacter coli</i> | <i>Helicobacter pylori</i> | <i>Shigella dysenteriae</i> |
| <i>Campylobacter jejuni</i> | <i>Listeria monocytogenes</i> | <i>Shigella flexneri</i> |
| <i>Clostridium difficile</i> | <i>Norovirus</i> | <i>Shigella sonnei</i> |
| <i>Cryptosporidium parvum</i> | <i>Salmonella enteritidis</i> | <i>Staphylococcus aureus</i> |
| <i>Entamoeba histolytica</i> | <i>Salmonella paratyphi</i> | <i>Yersinia enterocolitica</i> |
| <i>Enterovirus</i> | <i>Salmonella typhimurium</i> | |

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