ENGLISH - ESPAÑOL

Ref: RTB25F
RTB25FC
RTB25FCN

biotical FOB Card

Rapid Test to detect
human FOB in card
format

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biotical health
ready to test

IU RTB25F-RTB25FC-RTB25FCN Rev.03 / 2024-03

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical FOB Card

ENGLISH

WHAT IS BIOTICAL FOB CARD?

biotical FOB card is a coloured chromatographic immunoassay test (qualitative cut off test) for the detection of Human Haemoglobin (hHb) in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to detect faecal occult blood (gastrointestinal bleeding) in the lower digestive tract.

COLORECTAL CANCER AND FOB

Colorectal cancer is the second leading cause of disease and death in the western world. Screening, when using faecal occult blood tests (FOB), is based on the detection of bleeding from the digestive tract in stool samples: important symptom of neoplasms of the colon, cancer in its early stage and adenomatous polyps of the intestine.

Haemoglobin is a protein found in the red blood cells of vertebrates, which carries oxygen to the body's organs and tissues and carbon dioxide from the organs and tissues to the lungs. It can be released in the gastrointestinal tract and eliminated through the faeces, in diseases associated with gastrointestinal bleeding.

The presence of blood in the stool may be the only symptom of colorectal cancer, although not all the presence of blood in the stool is caused by cancer.

Colonoscopy, endoscopy and biopsy are the reference techniques used to detect colon and stomach cancer. These methods are invasive and require very laborious work to obtain a result. However, biotical FOB card, is a semi-quantitative immunochromatographic test, easy to use and with almost immediate results, being a very specific technique to detect blood loss in upper and lower digestive tract.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25F	biotical FOB Card	25 test without positive or negative controls.
RTB25FC	biotical FOB Card (+PC)	25 test with positive control, without negative control.
RTB25FCN	biotical FOB Card (+PC +NC)	25 test with positive and negative control.

HOW DOES BIOTICAL FOB CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In the **GREEN line (C)**: with rabbit polyclonal antibodies against a specific protein.
- In the **RED line (T)**: with mouse monoclonal antibodies against human Haemoglobin.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against human Haemoglobin) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex.

TEST LINES:

If the sample is positive: the antigens of the patient's sample will react with the red coloured conjugate complex (anti- human Haemoglobin monoclonal antibodies red latex microspheres for the line T) and advance by capillarity through the membrane. The anti-human Haemoglobin antibodies present in the T line will capture the coloured conjugate complex resulting in a red line.

If the sample is negative: if the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with red coloured complex. The anti-human Haemoglobin antibodies present in the T line will capture the coloured antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.

- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25FC & RTB25FCN references).
- Negative Control (only RTB25FCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.


Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. Use only fresh samples or fresh frozen samples without preservatives or fixatives. Once diluted, the sample can be stored for up to 2 days at room temperature or 5 days refrigerated.

Collection tube


1 Stool sample


2 Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

3 Cut the end of the cap.


4 Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 3 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID RESULT

INVALID RESULT

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>FOB</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>FOB</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25FC** & **RTB25FCN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical FOB Positive Control is an external quality control for biotical FOB card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical FOB Positive Control is manufactured with human haemoglobin extract dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by **rotating the positive control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes.**

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

NEGATIVE CONTROL: INTENDED USE & COMPOSITION


biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.


Results window

GREEN C line (Control)
RED T line (Test FOB)

Sample window (S)

- 1 Add 15 drops of diluent and introduce the negative control swab into the testing tube.
- 2 Rotating 30 seconds and extract the liquid.
- 3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of human Haemoglobin in the sample, without determining its origin (colorectal cancer, haemorrhoids, blood in urine, stomach irritations, etc.). To confirm its origin, it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques.
- False positives may occur in the test with bloody samples or mucous stools, if they have been taken during menstrual period, if the patient has haemorrhoids, blood in urine or constipation, so the use of the test in this class of samples is not recommended.

- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory or invasive techniques for confirmation. A negative result does not indicate absence of blood in the stool, as bleeding from the stomach can occur intermittently or the blood may not have been distributed throughout the stool sample. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity.

Cut-off value of biotical FOB Card is 50ng/mL (5.1µg hHb/g faeces) for human haemoglobin.

Clinical sensitivity and specificity.

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical FOB Card, biotical) and an ELISA test (Human Haemoglobin ELISA kit, ICL).

The results were as follows:

Data for cut-off value 5 µg hHb/g.

		Human Haemoglobin ELISA kit, ICL (ELISA test)		
		+	-	Total
Biotical FOB Card	+	38	11	49
	-	0	49	49
	Total	38	60	98

<i>Biotical FOB Card vs Human Haemoglobin ELISA kit, ICL</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99,9%	90,8% - 100 %
SPECIFICITY	81,7%	69,6% - 90,5%
PPV	77,6%	66,9% - 85,5%
NPV	>99,9%	-
ACCURACY	88,8%	80,8% - 94,3%

Specificity for a cut-off value of 1 µg hHb/g: 93,5% (mean value); 82,1-98,6% (95% confidence interval).

The results showed a high sensitivity and specificity to detect human haemoglobin using biotical FOB Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical FOB Card; no cross reactivity against gastrointestinal organism occasionally present in faeces.

<i>Adenovirus</i>	<i>Escherichia coli</i> O157:H7	<i>Salmonella enteritidis</i>
<i>Astrovirus</i>	<i>Giardia</i>	<i>Salmonella paratyphi</i> A
Human calprotectin	<i>Helicobacter pylori</i>	<i>Salmonella typhi</i>
<i>Campylobacter coli</i>	<i>Bovine Haemoglobin</i>	<i>Salmonella typhimurium</i>
<i>Campylobacter jejuni</i>	Pig Haemoglobin	<i>Shigella boydii</i>
<i>Clostridium difficile</i> antigen GDH	Bovine lactoferrin	<i>Shigella dysenteriae</i>
<i>Clostridium perfringens</i>	Human lactoferrin	<i>Shigella flexneri</i>
<i>Cryptosporidium</i>	<i>Legionella</i>	<i>Shigella sonnei</i>
<i>Entamoeba dispar</i>	<i>Listeria monocytogenes</i>	<i>Streptococcus pyogenes</i>
<i>Entamoeba histolytica</i>	Norovirus G1	Bovine transferrin
<i>Escherichia coli</i> O111:H8	Norovirus GII	Human transferrin
<i>Escherichia coli</i> O26:H11	<i>Rotavirus</i>	

SYMBOLS FOR IVD COMPONENTS & REAGENTS

	In vitro diagnostic device		Manufacturer
	Consult instructions for use		Contains sufficient for <n> test
	Keep dry	DIL	Sample diluent
	Use by	LOT	Batch code
	Temperature limitation	REF	Catalogue number

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REF

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REFERENCES / BIBLIOGRAFÍA /

1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.

Sierra de Guadarrama, 1. Ed. Dublín.
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es