ENGLISH - ESPAÑOL

Ref: RTB25TF
RTB25TFC
RTB25TFCN

biotical Transferrin Card

Rapid Test to detect
human transferrin in
card format

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biotical health
ready to test

IU RTB25TF-RTB25TFC-RTB25TFCN Rev. 02 / 2024-04

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical Transferrin Card

ENGLISH

WHAT IS BIOTICAL TRANSFERRIN CARD?

biotical Transferrin Card is a coloured chromatographic immunoassay test (qualitative cut-off test) for the determination of human transferrin (hTf) in stool samples. biotical Transferrin Card offers a simple, a highly sensitive and non-invasive screening assay to make a presumptive diagnosis of gastrointestinal bleeding.

COLORECTAL CANCER AND TRANSFERRIN

Transferrin is a blood-derived component that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases. Transferrin is stable in faeces and a good marker to detect loss of blood from the upper and lower intestine (gastrointestinal bleeding).

When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme and heme-derived porphyrins in amounts that depend on the site, the amount of bleeding and the transit time through the gut. Immunochemical tests detect intact or nearly intact human haemoglobin, being a very specific technique to detect the loss of blood from the lower intestine, because blood from lower sites is less degraded during transit. The haemoglobin is unstable in faeces, leading to false negative results.

The detection of faecal transferrin, which is more stable than haemoglobin, provides an alternative way of diagnosing the disease in the upper digestive tract.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25TF	biotical Transferrin Card	25 tests without positive or negative controls.
RTB25TFC	Biotical Transferrin Card (+PC)	25 tests with positive control, without negative control.
RTB25TFCN	biotical Transferrin Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL TRANSFERRIN CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against human Transferrin.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against human Transferrin) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to red polystyrene latex.

TEST LINES:

If the sample is positive: The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-human Transferrin monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-human Transferrin present in the T line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-human Transferrin antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25TFC & RTB25TFCN references).
- Negative Control (only RTB25TFCN reference).

The kit contains materials for the analysis of a maximum of 25 samples or controls.

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

SPECIMEN COLLECTION AND PREPARATION


Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):


- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

Collection tube

1 Stool sample

2 Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.


- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 3 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of Transferrin is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. Transferrin presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25TFC** & **RTB25TFCN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical Transferrin Positive Control is an external quality control for biotical TRANSFERRIN card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical Transferrin Positive Control is manufactured with human transferrin extract dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL - INTERPRETATION OF RESULTS

	biotal Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube **1**.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S **3**.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotal Test	Interpretation of the results
1	<p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>

2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>
---	--	---

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of human transferrin in the sample, nevertheless, it can be due to several causes, besides colorectal bleeding, such as haemorrhoids, blood in urine or stomach irritations. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source of the blood.
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory or invasive techniques for confirmation. A negative result does not indicate absence of blood in the stool, as bleeding from the stomach can occur intermittently or the blood may not have been distributed throughout the stool sample. It is possible to appear false negatives by not reaching the detection limit of the test.
- False positives may occur in the test with bloody samples or mucous stools, if they have been taken during menstrual period, if the patient has haemorrhoids, blood in urine or constipation, so the use of the test in this class of samples is not recommended.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity.

Cut-off value of biotical Transferrin is 4ng/mL (0.4µg hTf/g faeces) for human transferrin.

Clinical sensitivity and specificity.

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical transferrin, biotical) and an ELISA test (Human transferrin ELISA kit, ICL).

The results were as follows: Data for cut-off value 0.4µg hTf /g.

		Human Transferrin ELISA kit, ICL (ELISA test)		
		+	-	Total
Biotical Transferrin Card	+	38	3	41
	-	12	47	59
	Total	50	50	100

<i>Biotical Transferrin Card vs Human transferrin ELISA kit, ICL</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	76.0%	61.8 - 86.9%
SPECIFICITY	94.0%	83.5 - 98.7%
PPV	92.7%	80.1 - 98.5%
NPV	79.7%	67.2 - 98.0%
ACCURACY	85.0%	76.5 - 91.4%

Sensitivity for a cut-off value of 2 µg hTf/g: 81.1% (mean value); 64.8-92.0% (95% confidence interval).
 Specificity for a cut-off value of 0.2 µg hTf/g: 96.3% (mean value); 81.0-99.9% (95% confidence interval).

The results showed a high sensitivity and specificity to detect Transferrin using biotical Transferrin Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Transferrin Card; no cross reactivity against other gastrointestinal pathogens, substances and/or faecal markers present in faeces.

<i>Adenovirus</i>	<i>Entamoeba dispar</i>	Human lactoferrin	<i>Shigella boydii</i>
<i>Astrovirus</i>	<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Shigella dysenteriae</i>
Human calprotectin	<i>Escherichia coli</i> O111	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
<i>Campylobacter coli</i>	<i>Escherichia coli</i> O26	Norovirus G1	<i>Shigella sonnei</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli</i> O157	Norovirus GII	<i>Streptococcus pneumococcal</i>
<i>C. difficile</i> antigen GDH	<i>Giardia</i>	Rotavirus	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile</i> Toxin A	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>	Bovine transferrin
<i>Clostridium difficile</i> Toxin B	Bovine haemoglobin	<i>Salmonella paratyphi</i> A	<i>Yersinia enterocolitica</i> O:3
<i>Clostridium perfringens</i>	Pig haemoglobin	<i>Salmonella typhi</i>	<i>Yersinia enterocolitica</i> O:9
<i>Cryptosporidium</i>	Human haemoglobin	<i>Salmonella typhimurium</i>	

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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		ELISA test: Human Transferrin ELISA kit		
		+	-	Total
IC Test: Biotal Transferrin Card	+	38	3	41
	-	12	47	59
	Total	50	50	100

<i>Biotal Transferrin Card vs Human transferrin ELISA kit</i>	VALOR MEDIO	95% INTERVALO DE CONFIANZA
SENSIBILIDAD	76.0%	61.8 – 86.9%
ESPECIFICIDAD	94.0%	83.5 – 98.7%
PPV	92.7%	80.1 – 98.5%
NPV	79.7%	67.2 – 98.0%
PRECISIÓN	85.0%	76.5 – 91.4%

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REF

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REFERENCES / BIBLIOGRAFÍA /

1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. "Screening for colorectal cancer using the faecal occult blood test, Hemoccult". Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. "Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med". 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. "Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening". Ann Intern Med. 1997; 126: 808-810.
4. Chien-Hua Chiang et al., "A Comparative Study of Three Fecal Occult Blood Test in Upper Gastrointestinal Bleeding". Kaohsiung J. Med. Sci. May 2006, 22; 223-228.
5. Virtanen et al., "Higher concentrations of serum transferrin receptor in children than in adults", Am. J. Clin. Nutr. 1999; 69: 256-60.

Sierra de Guadarrama, 1. Ed. Dublín.
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es