ENGLISH - ESPAÑOL

Ref: RTB25FT
RTB25FTC
RTB25FTCN

biotical FOB + Transferrin Card

Rapid Test to detect
FOB + Transferrin in
card format

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biotical FOB+Transferrin Card

ENGLISH

WHAT IS BIOTICAL FOB+TRANSFERRIN CARD?

biotical FOB + Transferrin card is a colored chromatographic immunoassay test (qualitative cut off test) for the simultaneous detection of Human Hemoglobin (hHb) and Human Transferrin (hTf) in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to detect fecal occult blood (gastrointestinal bleeding) in lower and/or upper digestive tract.

FOB AND TRANSFERRIN

Colorectal cancer is the second leading cause of disease and death in the western world. Screening, when using fecal occult blood tests (FOB), is based on the detection of bleeding from the digestive tract in stool samples: important symptom of neoplasms of the colon, cancer in its early stage and adenomatous polyps of the intestine.

Hemoglobin is a protein found in the red blood cells of vertebrates, which carries oxygen to the body's organs and tissues and carbon dioxide from the organs and tissues to the lungs. It can be released in the gastrointestinal tract and eliminated through the faeces, in diseases associated with gastrointestinal bleeding.

Transferrin is a blood-derived component that can reach the gastrointestinal tract and be eliminated in the faeces due to some gastrointestinal diseases associated with blood loss. Transferrin is stable in the faeces and a good marker for detecting blood loss in the upper and lower part of the intestine (gastrointestinal bleeding). Hemoglobin is unstable in feces and can cause false negatives. Detection of "faecal" transferrin, more stable than hemoglobin, and is an alternative to diagnose gastrointestinal bleeding.

The presence of blood in the stool may be the only symptom of colorectal cancer, although not all the presence of blood in the stool is caused by cancer.

Colonoscopy, endoscopy and biopsy are the reference techniques used to detect colon and stomach cancer. These methods are invasive and require very laborious work to obtain a result. However, biotical FOB + Transferrin card, is a semi-quantitative immunochromatographic test, easy to use and with almost immediate results, being a very specific technique to detect blood loss in upper and lower digestive tract.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25FT	biotical FOB+Transferrin Card	25 tests without positive or negative controls.
RTB25FTC	biotical FOB+Transferrin Card (+PC)	25 tests with positive control, without negative control.
RTB25FTCN	biotical FOB+Transferrin Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL FOB+TRANSFERRIN CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In **GREEN** line (Control): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (Test): with mouse monoclonal antibodies against human Transferrin (hTf).
- In **BLUE** line (Test): with mouse monoclonal antibodies against human Hemoglobin (hHb).
- **Sample area (S)**: the pad where the sample is placed has been dispensed with: a preparation of reagents of the **RED** line (mouse monoclonal antibodies against hTf) conjugated to red polystyrene latex; with another preparation of reagents of the **BLUE** line (mouse monoclonal antibodies against hHb) conjugated to blue polystyrene latex; and with another preparation for the control solution (specific binding protein) conjugated to **GREEN** polystyrene latex; forming three conjugate colored complexes.

TEST LINES:

If the sample is positive: the antigens of the patient's stool sample will react with the red colored conjugate complex (anti-human Transferrin monoclonal antibodies red latex microspheres) and / or blue colored conjugate complex (anti-human Hemoglobin monoclonal antibodies blue latex microspheres) and advance by capillarity through the membrane. The anti-human Transferrin antibodies present in the RED line will capture the colored conjugate complex resulting in a red line, and the anti-human Hemoglobin antibodies present in the BLUE line will capture the colored conjugate complex resulting in a blue line.

If the sample is negative: if the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red- or blue-colored complexes. The antihuman Transferrin and anti-human Hemoglobin antibodies present in the RED and BLUE line will not capture the colored antigen-complex (not formed) and the red and blue lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.

- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Collection tubes with diluent.
- Positives Controls (only RTB25FTC & RTB25FTCN references).
- Negative Control (only RTB25FTCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION.

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen -20°C. In that case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.


Patients should not collect samples during their menstrual period if they have bleeding haemorrhoids, blood in urine or if they have strained during bowel movement.

Specimen preparation (see illustration):


- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. The stool collection tube with dilute sample can be stored until 2 days at room temperature or until 5 days in the refrigerator (2-8°C) prior to testing.

Collection tube


1 Stool sample


2 Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

3 Cut the end of the cap.

4 Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

TRANSFERRIN NEGATIVE
FOB NEGATIVE

TRANSFERRIN POSITIVE
FOB NEGATIVE


TRANSFERRIN NEGATIVE
FOB POSITIVE


TRANSFERRIN POSITIVE
FOB POSITIVE

INVALID

INVALID

INVALID

INVALID

	biotical Test	Interpretation of the results
1	GREEN LINE NO LINE NO LINE	NEGATIVE FOR TRANSFERRIN & HEMOGLOBIN Only the GREEN line appears. No presence of human Transferrin (hTF) or FOB/human Hemoglobin (hHB) is detected in the sample.
2	GREEN LINE RED LINE NO LINE	POSITIVE FOR TRANSFERRIN – NEGATIVE FOR HEMOGLOBIN In addition to the GREEN control line, a RED line also appears. Human Transferrin (hTF) presence is detected in the sample. No presence of FOB/human Haemoglobin (hHB) is detected.
3	GREEN LINE NO LINE BLUE LINE	POSITIVE FOR HEMOGLOBIN – NEGATIVE FOR TRANSFERRIN In addition to the GREEN control line, a BLUE line also appears. FOB/human Haemoglobin (hHB) presence is detected in the sample. No presence of Human Transferrin (hTF) is detected.
4	GREEN LINE RED LINE BLUE LINE	POSITIVE FOR TRANSFERRIN & HEMOGLOBIN In addition to the GREEN control line, a RED and BLUE lines also appear. Human Transferrin (hTF) and FOB/human Haemoglobin (hHB) presence are detected in the sample.

	biotal Test	Interpretation of the results
5-8	Any other result	<p>INVALID</p> <p>Total absence of the control line (GREEN) regardless the appearance or not of any RED – BLUE lines cause an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.</p>

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the coloured RED-BLUE bands in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (GREEN) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25FTC** & **RTB25FTCN**).

POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotal Transferrin Positive Control and biotal FOB Positive Control are two external quality controls for biotal FOB + Transferrin card. Use these controls swabs to check that the sample reagents and the test card are working properly.

biotal Transferrin Positive Control: transferrin extract dried in swab containing a preservative. biotal FOB Positive Control: hemoglobin extract dried in swab containing a preservative.

POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- 1 Add 15 drops of diluent and introduce the positive control swab into the testing tube, 2 rotating 30 seconds and extract the liquid.
 - 3 Add 4 drops in the circular window (S).
- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
 - Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
 - Remove the biotical device from its sealed bag just before using it.
 - Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
 - Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

biotical Test	Interpretation of the results
<p>GREEN LINE RED LINE NO BLUE LINE</p>	<p>RESULT: human Transferrin (hTF) POSITIVE CONTROL.</p> <p>A GREEN line (Control) and a RED line appear across the results window during the test performance due to presence of the positive control swab.</p>
<p>GREEN LINE NO RED LINE BLUE LINE</p>	<p>RESULT: human Haemoglobin (hHB) POSITIVE CONTROL.</p> <p>A GREEN line (Control) and a BLUE line appear across the results window during the test performance due to presence of the positive control swab.</p>
	<p>INVALID RESULTS FOR POSITIVE CONTROL.</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of any of the test lines (RED - BLUE), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.


1 Add 15 drops of diluent and introduce the positive control swab into the testing tube, 2 rotating 30 seconds and extract the liquid.

3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube 1.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the Biotal test Card device from its sealed bag just before using it.
- Dispense 4 drops from the testing tube in the corresponding circular window 3. Use a separate pipette and test for each control.
- Read the result at 10 minutes. Do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

biotical Test		Interpretation of the results
	GREEN LINE NO RED LINE NO BLUE LINE	RESULT: CONTROL NEGATIVE. Only the control line (C) appears across the results window during the test performance.
		RESULT: INVALID FOR NEGATIVE CONTROL. Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new negative control swab.

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of human hemoglobin in the sample, without determining its origin (hemorrhoids, blood in urine, stomach irritations, colorectal cancer, etc.). To confirm its origin, it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques.
- False positives may occur in the test with bloody samples or mucous stools, if they have been taken during menstrual period, if the patient has haemorrhoids, blood in urine or constipation, so the use of the test in this class of samples is not recommended.
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory or invasive techniques for confirmation. A negative result does not indicate absence of blood in the stool, as bleeding from the stomach can occur intermittently or the blood may not have been distributed throughout the stool sample. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit).

Cut-off values of biotical FOB + Transferrin are 50ng/mL (5.1µg hHb/g faeces) for human haemoglobin and 4ng/mL (0.4µg hTf/g faeces) for human transferrin.

Clinical sensitivity and specificity.

An evaluation with faecal samples was performed comparing the results obtained by an immunochromatographic test (biotical FOB + Transferrin, biotical) and others commercial ELISA tests (Human Hemoglobin ELISA kit, ICL) and (Human Transferrin ELISA Kit, ICL).

The results for hTF were as follows:(data for cut-off value 0.4 µg hTF/g).

		Human Transferrin ELISA Kit, ICL (ELISA test)		
		+	-	Total
Biotical FOB+TRANSFERRINA Card (RED LINE: TRANSFERRIN)	+	38	3	41
	-	12	47	59
	Total	50	50	100

Biotical FOB+Transferrina Card vs Human Transferrin ELISA Kit, ICL (TRANSFERRIN)	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	76,0%	61,8% - 86,9%
SPECIFICITY	94,0%	83,5% - 98,8%
PPV	92,7%	80,7% - 97,5%
NPV	79,7%	70,4% - 86,6%
ACCURACY	85,0%	76,5% - 91,4%

Sensitivity for a cut-off value of 2 µg hTf/g: 81.1% (mean value); 64.8-92.0% (95% confidence interval).
Specificity for a cut-off value of 0.2 µg hTf/g: 96.3% (mean value); 81.0-99.9% (95% confidence interval).

The results for hHG were as follows: (data for cut-off value 5 µg hHb/g).

		Human Hemoglobin ELISA kit, ICL (ELISA test)		
		+	-	Total
Biotical FOB+TRANSFERRINA Card (BLUE HEMOGLOBIN)	+	38	11	49
	-	0	49	49
	Total	38	60	98

Biotical FOB+Transferrina Card vs Evaluation Giardia (Giardia)	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99,9%	90,8% - 100,0%
SPECIFICITY	81,7%	69,6% - 90,5%
PPV	77,6%	66,9% - 85,5%
NPV	>99,9%	-
ACCURACY	88,8%	80,8% - 94,3%

Specificity for a cut-off value of 1 µg hHb/g: 93.5% (mean value); 82.1-98.6% (95% confidence interval).

The results showed a high sensitivity and specificity to detect human haemoglobin and human transferrin using biotical FOB + Transferrin.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical FOB + Transferrin; no cross reactivity against other gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces.

For RED LINE: Transferrin detection line:

<i>Adenovirus</i>	<i>Entamoeba dispar</i>	<i>Human lactoferrin</i>	<i>Shigella boydii</i>
<i>Astrovirus</i>	<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Shigella dysenteriae</i>
<i>Human calprotectin</i>	<i>Escherichia coli</i> O111	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
<i>Campylobacter coli</i>	<i>Escherichia coli</i> O26	<i>Norovirus GI</i>	<i>Shigella sonnei</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli</i> O157	<i>Norovirus GII</i>	<i>Streptococcus pneumococcal</i>
<i>Clostridium difficile antigen GDH</i>	<i>Giardia</i>	<i>Rotavirus</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile Toxin A</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>	<i>Bovine transferrin</i>
<i>Clostridium difficile Toxin B</i>	<i>Bovine hemoglobin</i>	<i>Salmonella paratyphi A</i>	<i>Yersinia enterocolitica</i> O:3
<i>Clostridium perfringens</i>	<i>Pig hemoglobin</i>	<i>Salmonella typhi</i>	<i>Yersinia enterocolitica</i> O:9
<i>Cryptosporidium</i>	<i>Human hemoglobin</i>	<i>Salmonella typhimurium</i>	

For BLUE LINE: FOB detection line:

<i>Adenovirus</i>	<i>Entamoeba dispar</i>	<i>Human lactoferrin</i>	<i>Shigella boydii</i>
<i>Astrovirus</i>	<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Shigella dysenteriae</i>
<i>Human calprotectin</i>	<i>Escherichia coli</i> O111	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
<i>Campylobacter coli</i>	<i>Escherichia coli</i> O26	<i>Norovirus GI</i>	<i>Shigella sonnei</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli</i> O157	<i>Norovirus GII</i>	<i>Streptococcus pneumococcal</i>
<i>Clostridium difficile antigen GDH</i>	<i>Giardia</i>	<i>Rotavirus</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile Toxin A</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>	<i>Bovine transferrin</i>
<i>Clostridium difficile Toxin B</i>	<i>Bovine hemoglobin</i>	<i>Salmonella paratyphi A</i>	<i>Human transferrin</i>
<i>Clostridium perfringens</i>	<i>Pig hemoglobin</i>	<i>Salmonella typhi</i>	<i>Yersinia enterocolitica</i> O:3
<i>Cryptosporidium</i>	<i>Bovine lactoferrin</i>	<i>Salmonella typhimurium</i>	<i>Yersinia enterocolitica</i> O:9

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. "Screening for colorectal cancer using the faecal occult blood test, Hemoccult". Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. "Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med". 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. "Suggested technique for Faecal Occult Blood testing and interpretation in colo-rectal cancer screening". Ann Intern Med. 1997; 126: 808-810.
4. Chien-Hua Chiang et al., "A Comparative Study of Three Fecal Occult Blood Test in Upper Gastrointestinal Bleeding". Kaohsiung J. Med. Sci. May 2006, 22; 223-228.
5. Virtanen et al., "Higher concentrations of serum transferrin receptor in children than in adults", Am. J. Clin. Nutr. 1999; 69: 256-60.

Sierra de Guadarrama, 1. Ed. Dublín
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es