ENGLISH - ESPAÑOL

Ref: RTB25LT
RTB25LTC
RTB25LTCN

biotical Lactoferrin Card

Rapid Test to detect
human lactoferrin in
card format

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biotical health
ready to test

IU RTB25LT-RTB25LTC-RTB25LTCN Rev. 02 / 2024-04

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical Lactoferrin Card

ENGLISH

WHAT IS BIOTICAL LACTOFERRIN CARD?

biotical Lactoferrin card is a coloured chromatographic immunoassay test (qualitative cut off test) for the detection of human Lactoferrin (hLf) in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of Inflammatory Bowel Diseases (IBD) as well as a follow up of the treatment and to predict the risk of relapse.

GASTROINTESTINAL INFLAMMATION AND LACTOFERRIN

Lactoferrin is a glycoprotein component of neutrophil secondary granules, a primary component of the acute inflammatory response and is released from faecal leukocytes. This protein is resistant to proteolysis in the faeces and may serve as a marker of inflammation in the intestine. The major cause of faecal neutrophils in patients with chronic diarrhoea is chronic inflammatory bowel disease of the colon (i.e., Crohn's Disease and ulcerative colitis). Lactoferrin has been also studied as a predictor of infection with invasive enteropathogens in children with diarrhoea. Bacterial inflammatory diarrhoea may be caused by Shigella, Salmonella, Campylobacter and Clostridium difficile.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25LT	biotical Lactoferrin Card	25 tests without positive or negative controls.
RTB25LTC	biotical Lactoferrin Card (+PC)	25 tests with positive control, without negative control.
RTB25LTCN	biotical Lactoferrin Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL LACTOFERRIN CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against human Lactoferrin (hLf).
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents of the T line (mouse monoclonal antibodies against human Lactoferrin) conjugated with latex of red polystyrene; and a control solution (specific binding protein) conjugated to green polystyrene latex, forming two conjugate-coloured complexes.

TEST LINE:

If the sample is positive: The antigens of the patient's sample will react with the red coloured conjugate complex (anti-human Lactoferrin monoclonal antibodies red latex microspheres for the line T) and advance by capillarity through the membrane. The anti-human Lactoferrin antibodies present in the T line will capture the coloured test complex resulting in the red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-human Lactoferrin antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25LTC & RTB25LTCN references).
- Negative Control (only RTB25LTCN reference)

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION


Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. The stool collection vial with diluted sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

Collection tube


1 Stool sample


2 Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

3 Cut the end of the cap.


4 Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID

INVALID

	biotal Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of human lactoferrin is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. Human lactoferrin presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references ended with **RTB25LTC** & **RTB25LTCN**).


POSITIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Lactoferrin Positive Control is an external quality control for biotal Lactoferrin card. Use this control swab to check that the sample reagents and the test card are working properly.

biotal Lactoferrin Positive Control is manufactured with inactivated Lactoferrin extract dried in swab containing a preservative.

POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by **rotating the positive control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes**. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotal Test	Interpretation of the results
1	 <p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. Do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test must be carried out within 2 hours after opening the sealed bag.
- An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- The intensity of test line may vary from very strong at high human lactoferrin concentration to faint when the human lactoferrin concentration is close to the cut-off value of the test.
- biotical Lactoferrin Card should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of human lactoferrin in faecal samples; nevertheless, it can be also due to several causes besides inflammatory bowel disease (IBD). A positive result should be followed up with additional diagnostic procedures. Endoscopy and histology on biopsy specimens are the methods for detecting and quantifying bowel inflammation.
- A negative result is not meaningful because of it is possible the human lactoferrin concentration in the stool sample is lower than the cut-off value. Negative results do not exclude inflammation, some diseases such as celiac sprue and microscopic colitis polyps that involve mainly mononuclear inflammation.

- Lactoferrin is a component of breast milk; the test will be positive in breast fed children and should not be used to evaluate neonates receiving breast milk.
- Stool from patients with active inflammatory bowel diseases that usually involve significant neutrophilic inflammation of the intestine, such as Crohn’s disease and ulcerative colitis, would be positive for faecal lactoferrin. biotical Lactoferrin Card could be sensitive for this diagnosis in patients with chronic diarrhoea.
- Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
- Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity:

Cut-off value of biotical Lactoferrin Card is 100ng/mL (10µg hLf/g faeces) for human lactoferrin.

Clinical sensitivity and specificity:

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical Lactoferrin, biotical) and other commercial ELISA tests (Human Lactoferrin ELISA kit, ICL).

The results were as follows:

		ELISA test: Human Lactoferrin ELISA kit		
		+	-	Total
Biotical Lactoferrin Card	+	27	13	40
	-	3	47	50
	Total	30	60	90

<i>Biotical Lactoferrin Card vs Human Lactoferrin ELISA kit</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	90.00%	73.5% – 97.9%
SPECIFICITY	78.3%	65.8% – 87.9%
PPV	67.5%	50.9% – 81.4%
NPV	94.0%	83.5% – 98.0%
ACCURACY	82.2%	72.7% – 89.4%

Specificity for a cut-off value of 2.5 µg hLf/g: 100.0% (mean value); 83.9-100.0% (95% confidence interval).

The results showed a high sensitivity and specificity to detect human lactoferrin using biotical Lactoferrin Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Lactoferrin Card; no cross reactivity against other faecal markers occasionally present in faeces:

<i>Adenovirus</i>	<i>Entamoeba dispar</i>	Bovine lactoferrin	<i>Shigella boydii</i>
<i>Astrovirus</i>	<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Shigella dysenteriae</i>
Human calprotectin	<i>Escherichia coli</i> O111	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
<i>Campylobacter coli</i>	<i>Escherichia coli</i> O26	Norovirus G1	<i>Shigella sonnei</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli</i> O157	Norovirus GII	<i>Streptococcus pneumococcal</i>
<i>C. difficile</i> antigen GDH	<i>Giardia</i>	Rotavirus	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile</i> Toxin A	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>	Bovine transferrin
<i>Clostridium difficile</i> Toxin B	Bovine haemoglobin	<i>Salmonella paratyphi</i> A	Human transferrin
<i>Clostridium perfringens</i>	Pig haemoglobin	<i>Salmonella typhi</i>	<i>Yersinia enterocolitica</i> O:3
<i>Cryptosporidium</i>	Human haemoglobin	<i>Salmonella typhimurium</i>	<i>Yersinia enterocolitica</i> O:9

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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REF

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REFERENCES / BIBLIOGRAFÍA /

1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. Clinica Chimica Acta 381 Feb 2007, p. 63-68.
2. Quail, M.A. et al. Fecal Lactoferrin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. Inflamm Bowel Dis, Vol 15 No 5; May 2009, p. 756-759.
3. Gaya D.R., et al. Faecal lactoferrin in the assessment of Crohn's disease activity. Q J Med 2005, Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Lactoferrin and PMN-Elastase, CRP, and Clinical Indices. Am. J. Gastroenterol. 2008; Vol 103, p. 162-169.

Sierra de Guadarrama, 1. Ed. Dublín.
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es

biotical health S.L.U.