ENGLISH - ESPAÑOL

Ref: RTB25CL
RTB25CLC
RTB25CLCN

biotical

Calprotectin + Lactoferrin Card

Rapid Test to detect
Calprotectin +
Lactoferrin in card
format

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biotical health
ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical

ENGLISH

Calprotectin + Lactoferrin Card

WHAT IS BIOTICAL CALPROTECTIN + LACTOFERRIN CARD?

biotical Calprotectin + Lactoferrin card is a coloured chromatographic immunoassay test (qualitative cut off test) for the simultaneous detection of Human Calprotectin (hCp) and Human Lactoferrin (hLf) in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis and detect gastrointestinal inflammation present in various pathologies.

GASTROINTESTINAL INFLAMMATION, CALPROTECTIN AND LACTOFERRIN

Human Calprotectin (hCp) is a calcium-associated protein that has bacteriostatic and fungicidal properties. In inflammatory processes their plasma levels rise considerably, with high levels of hCp being found in the faeces of patients with various inflammatory bowel diseases (IBD), such as Crohn's disease and ulcerative colitis. This faecal hCp assay is useful for differentiating an Inflammatory Bowel Disease (IBD) from functional Irritable Bowel Syndrome (IBS). hCp, once released until its degradation, remains stable at ambient temperature for 7 days in the faeces, thus being an ideal analyte. This is because it is extremely resistant to proteolytic degradation (proteolysis).

Human Lactoferrin (hLf) is a glycoprotein that is present in secreted fluids such as saliva, tears and breast milk, and acts as a non-specific defense protein. When inflammation occurs in the gastrointestinal tract, hLf is released to prevent bacterial growth by limiting the availability of iron and causing direct damage to the cell membrane in cooperation with lysozymes. Its effect is increased by the presence of specific IgA antibodies secreted directly against bacteria.

hLf is elevated in faeces of patients with Inflammatory Bowel Diseases (IBD), such as ulcerative colitis, Crohn's disease and other pathologies characterized by idiopathic and chronic inflammation that affects the gastrointestinal tract. Paediatric patients and adults with IBD may have a wide variety of clinical symptoms (including abdominal pain and diarrhoea) that may be non-specific.

biotical Calprotectin + Lactoferrin card is a simple, non-invasive test, easy to use and with almost immediate results, based on the detection of Human Calprotectin and Human Lactoferrin in stool samples. This serves to distinguish patients with inflammation (invasive bacterial infection, IBD, etc.) from those without inflammation (viral, toxigenic, etc.). In addition, the determination of faecal Human Calprotectin and Human Lactoferrin concentrations may be useful for predicting imminent clinical recurrence in various IBD. Considering that children need general anaesthesia for colonoscopies, the benefits will be even greater.

TABLE OF REFERENCES


In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25CL	biotical Calprotectin + Lactoferrin Card	25 tests without positive or negative controls.
RTB25CLC	biotical Calprotectin + Lactoferrin Card (+PC)	25 tests with positive control, without negative control.
RTB25CLCN	biotical Calprotectin + Lactoferrin Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL CALPROTECTIN+LACTOFERRIN CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In **RED T₁ line (Calprotectin)**: with mouse monoclonal antibodies against human Calprotectin (hCp).
- In **GREEN C line (Control)**: with rabbit polyclonal antibodies against a specific protein.
- In **RED T₂ line (Lactoferrin)**: with mouse monoclonal antibodies against human Lactoferrin (hLf).
- **Sample area (S)**: The pad where the sample is placed has been dispensed with a preparation of reagents of the T1 line (mouse monoclonal antibodies against human Calprotectin) conjugated with latex of red polystyrene; with another preparation of reagents of the T2 line (mouse monoclonal antibodies against hLf) conjugated to red polystyrene latex; and with another preparation for the control solution (specific binding protein) conjugated to green polystyrene latex, forming three conjugate coloured complexes.

TEST LINES:

If the sample is positive: the antigens of the patient's stool sample will react with one (or both) of the red coloured conjugate complexes (anti-human Calprotectin monoclonal antibodies-latex red microspheres for the line T1 and anti- human Lactoferrin monoclonal antibodies-red latex microspheres for the line T2) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-human Calprotectin present in the T1 line and / or the antibodies anti-human Lactoferrin present in the T2 line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complexes. The anti-human Calprotectin and anti-human Lactoferrin antibodies present in the T1 and T2 line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Collection tubes with diluent.
- Positives Controls (only RTB25CLC & RTB25CLCN references).
- Negative Control (only RTB25CLCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION


Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):


- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample ①, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion ②. The stool collection tube with dilute sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

Collection tube



① Stool sample


② Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

③ Cut the end of the cap.


④ Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)


CALPROTECTIN NEGATIVE
LACTOFERRIN NEGATIVE

CALPROTECTIN POSITIVE
LACTOFERRIN NEGATIVE

CALPROTECTIN NEGATIVE
LACTOFERRIN POSITIVE

CALPROTECTIN POSITIVE
LACTOFERRIN POSITIVE

INVALID

INVALID

INVALID

INVALID

	biotal Test	Interpretation of the results
1	T1 = NO LINE C = GREEN LINE T2 = NO LINE	NEGATIVE FOR CALPROTECTIN & LACTOFERRIN Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of Calprotectin (hCp) or Lactoferrin (hLf) are detected in the sample.
2	T1 = RED LINE C = GREEN LINE T2 = NO LINE	POSITIVE FOR CALPROTECTIN – NEGATIVE FOR LACTOFERRIN In addition to the GREEN control line, a RED line also appears in the site marked with the letter T1 (hCp detection line) in the results window. Calprotectin (hCp) presence is detected in the sample. No presence of Lactoferrin (hLf) is detected.

	biotical Test	Interpretation of the results
3	T1 = NO LINE C = GREEN LINE T2 = RED LINE	POSITIVE FOR LACTOFERRIN – NEGATIVE FOR CALPROTECTIN In addition to the GREEN control line, a RED line also appears in the site marked with the letter T2 (hLf detection line) in the results window. Lactoferrin (hLf) presence is detected in the sample. No presence of Calprotectin (hCp) is detected.
4	T1 = RED LINE C = GREEN LINE T2 = RED LINE	POSITIVE FOR CALPROTECTIN & LACTOFERRIN In addition to the GREEN control line, a RED line also appears in the site marked with the letter T1 (hCp detection line) and a RED line appears in the site marked with the letter T2 (hLf detection line) in the results window. Calprotectin (hCp) and Lactoferrin (hLf) presence are detected in the sample.
5-8	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of any test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T1 – T2) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25CLC** & **RTB25CLCN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical Calprotectin Positive Control and biotical Lactoferrin Positive Control are two external quality controls for biotical Calprotectin + Lactoferrin card. Use these control swabs to check that the sample reagents and the test card are working properly.

biotical Calprotectin Positive Control: calprotectin extract dried in swab containing a preservative. biotical Lactoferrin Positive Control: lactoferrin extract dried in swab containing a preservative.

POSITIVE CONTROL: TESTING PROCEDURE


Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biological device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

biological Test	Interpretation of the results
 <p>T1 = RED LINE C = GREEN LINE T2 = NO LINE</p>	<p>RESULT: CALPROTECTIN POSITIVE CONTROL.</p> <p>A GREEN line (C) and a RED line (T1) appear across the results window during the test performance due to presence of the positive control swab.</p>
 <p>T1 = NO LINE C = GREEN LINE T2 = RED LINE</p>	<p>RESULT: LACTOFERRIN POSITIVE CONTROL.</p> <p>A GREEN line (C) and a RED line (T2) appear across the results window during the test performance due to presence of the positive control swab.</p>

INVALID RESULTS FOR POSITIVE CONTROL.

Total absence of the control line (GREEN), regardless the appearance or not of any of the test lines (RED), or only appears the control line.

Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.


NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the negative control swab into the testing tube **1**.
- Mix the solution by **rotating the negative control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes**. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

biotical Test	Interpretation of the results
 <p>T1 = NO LINE C = GREEN LINE T2 = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE.</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
	<p>RESULT: INVALID FOR NEGATIVE CONTROL.</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test must be carried out up to 2 hours after opening the bag that contains it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. In case of occur, the test should be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with samples from another source (non-faecal or non-human) has not been tested so is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of Calprotectin and/or Lactoferrin in the sample, without determining its origin (IBD, colorectal cancer, enteropathies, etc.). To confirm inflammatory bowel disease, it is necessary for a specialist to evaluate the results together with the patient's clinical evidence. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests by invasive diagnostic techniques (endoscopy or biopsy).
- False positives may occur in the test in samples from patients previously treated with non-inflammatory steroids (may cause positive Calprotectin results) or infants and new-borns (may cause positive Lactoferrin results), bloody samples or mucous stools (could cause Positive Calprotectin and Lactoferrin results from nonspecific reactions). The use of the test in this class of samples is not recommended.

- If a negative test result is obtained and symptoms persist, it is advisable to repeat the test after a few days. Negative test results may be inconclusive and invasive techniques are required for confirmation. A negative result does not indicate absence of intestinal inflammation. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit).

Cut-off values of biotical Calprotectin + Lactoferrin are 500ng/mL (50µg hCp/g faeces) for human calprotectin and 100ng/mL (10µg hLf /g faeces) for human lactoferrin.

Clinical sensitivity and specificity.

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical Calprotectin + Lactoferrin Card, biotical) and others commercial immunoassays (Calprest®, Eurospital, and IBD EZ VUE®, TechLab®).

The results for hCp were as follows:

		Calprest®, Eurospital (immunochromatographic)		
		+	-	Total
Biotical Calprotectin + Lactoferrin Card (T1: CALPROTECTIN)	+	34	2	36
	-	2	26	28
	Total	36	28	64

<i>Biotical Calprotectin + Lactoferrin Card vs Calprest®, Eurospital (CALPROTECTIN)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	94,4%	81,3% - 99,3%
SPECIFICITY	92,9%	76,5% - 99,1%
PPV	94,4%	81,7% - 98,5%
NPV	92,9%	77,1% - 98,1%
ACCURACY	93,8%	84,8% - 98,3%

The results for hLf were as follows:

		IBD EZ VUE®, TechLab® (immunochromatographic)		
		+	-	Total
Biotal Calprotectin + Lactoferrin Card (T2: LACTOFERRIN)	+	20	0	20
	-	0	44	44
	Total	20	44	64

<i>Biotal Calprotectin + Lactoferrin Card vs IBD EZ VUE®, TechLab® (LACTOFERRIN)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99,9%	83,16% - 100,00%
SPECIFICITY	>99,9%	91,96% - 100,00%
PPV	>99,9%	-
NPV	>99,9%	-
ACCURACY	>99,9%	94,40% - 100,00%

The results showed a high sensitivity and specificity to detect human calprotectin and human lactoferrin using biotal Calprotectin + Lactoferrin Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotal Calprotectin + Lactoferrin Card; no cross reactivity against other faecal markers occasionally present in faeces:

Bovine and pig haemoglobin
Human haemoglobin

Bovine and pig transferrin
Human lactoferrin (T1 Line)

Human calprotectin (T2 Line)
Human transferrin

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. Feb 2007, p. 63-68 2.
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. , Vol 15 No 5; May 2009, p. 756-759
3. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. , Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. 2008; Vol 103, p. 162-169.
5. Gisbert, J.P. et al. Role of Biological Markers in Inflammatory Bowel Disease. Gastroenterol. Hepatol. May 2007,30(3), p. 117-29.

Sierra de Guadarrama, 1.
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es

biotical health S.L.U.