ENGLISH - ESPAÑOL

Ref: RTB25R  
RTB25RC  
RTB25RCN

# biotical Rotavirus Card

Rapid Test to detect  
Rotavirus in card format

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biotical health  
ready to test

IU RTB25R-RTB25RC-RTB25RCN Rev. 03 / 2024-04

ISO 13485 CERTIFIED  
MEDICAL DEVICE  
QUALITY MANAGEMENT SYSTEM

# biotical Rotavirus Card

ENGLISH

## WHAT IS BIOTICAL ROTAVIRUS CARD?

biotical Rotavirus card is a coloured chromatographic immunoassay test for the qualitative detection of Rotavirus antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of Rotavirus infection.

## GASTROENTERITIS AND ROTAVIRUS

Rotavirus and Adenovirus are the main cause of severe diarrhoea in children under 5 years, causing high morbidity and mortality worldwide.

Acute gastroenteritis is a very common disease that affects children and young people around the world. Enteric Rotavirus is the first cause of severe diarrhoea in children under 5 years of age, followed by Adenovirus.

Rotavirus is transmitted by feco-oral contact, causing watery diarrhoea, vomiting, headache, fever and stomachache. Symptoms usually begin 1 or 2 days after becoming infected and may last 1 to 3-4 days.

The most common procedures to identify a gastroenteritis caused by Rotavirus involve the isolation and identification of the virus through techniques that require more than 1 day. However, biotical Rotavirus card, is an immunochromatographic qualitative test, easy to use and with almost immediate results, for the detection of Rotavirus directly from the patient's stool sample, providing results in 10 minutes.


## TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25R	biotical Rotavirus Card	25 tests without positive or negative controls.
RTB25RC	biotical Rotavirus Card (+PC)	25 tests with positive control, without negative control.
RTB25RCN	biotical Rotavirus Card (+PC +NC)	25 tests with positive and negative control.

## HOW DOES BIOTICAL ROTAVIRUS CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against *Rotavirus*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against *Rotavirus*) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex.

### TEST LINES:

**If the sample is positive:** The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-*Rotavirus* monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-*Rotavirus* present in the T line will capture the complex giving rise in a red line.

**If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-*Rotavirus* antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

### CONTROL LINE (C):

**Regardless of whether the sample is positive or not**, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

## STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

## GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.

- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles, and mask. Do not eat, drink, or smoke in the working area.

## MATERIALS

### MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25RC & RTB25RCN references).
- Negative Control (only RTB25RCN reference).

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

## SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.

### Specimen preparation (see illustration):

- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample (aprox. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**.

Collection tube


**1** Stool sample

**2** Shake vigorously

## TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

3 Cut the end of the cap.


4 Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	<b>NEGATIVE</b> Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>Rotavirus</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	<b>POSITIVE</b> In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>Rotavirus</i> presence is detected in the sample.
3-4	Any other result	<b>INVALID</b> Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25RC** & **RTB25RCN**).


## POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical Rotavirus Positive Control is an external quality control for biotical ROTAVIRUS card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical Rotavirus Positive Control: inactivated *Rotavirus* antigens dried in swab containing a preservative.

## POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

## POSITIVE CONTROL - INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p><b>RESULT: CONTROL POSITIVE</b></p> <p>A GREEN line (C) and a RED line (T) appear across the results window (C) during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p><b>INVALID RESULT</b></p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

## NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube 1.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotal Test	Interpretation of the results
1	<p>C = GREEN LINE T = NO LINE</p>	<p><b>RESULT: CONTROL NEGATIVE</b></p> <p>Only the control line (C) appears across the results window during the test performance.</p>

2		<p><b>RESULT: INVALID FOR NEGATIVE CONTROL</b></p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>
---	--	---

## LIMITATIONS OF THIS TEST

- The test must be performed within 2 hours after opening the sealed pouch.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- After one week of infection, the number of viruses in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
- This Card should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of Rotavirus in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods or microscopy) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation still persist, a Rotavirus determination should be carried out with another technique (for example microscopy).

## PERFORMANCE CHARACTERISTICS

### Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using an immunochromatographic test (biotical Rotavirus, biotical) and these results were confirmed with a commercially available test (Ridascreen® Rotavirus ELISA Test, r-Biopharm).

The results were as follows:

		Ridascreen® Rotavirus ELISA Test		
		+	-	Total
Biotical Rotavirus Card	+	18	1	19
	-	0	43	43
	Total	18	44	62

<i>Biotal Rotavirus Card vs Ridascreen® Rotavirus ELISA Test</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99%	81.5% - 100%
SPECIFICITY	97.7%	88.0% - 99.9%
PPV	94.7%	72.2% - 99.2%
NPV	>99%	-
ACCURACY	98.4%	91.3% - 99.9%

The results showed a high sensitivity and specificity to detect Rotavirus using biotal Rotavirus Card.

## CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotal Rotavirus Card; no cross reactivity against gastrointestinal organism occasionally present in faeces.

- |                               |                                 |                                |
|-------------------------------|---------------------------------|--------------------------------|
| <i>Adenovirus</i>             | <i>Escherichia coli O157:H7</i> | <i>Salmonella typhi</i>        |
| <i>Astrovirus</i>             | <i>Giardia lamblia</i>          | <i>Shigella boydii</i>         |
| <i>Campylobacter coli</i>     | <i>Helicobacter pylori</i>      | <i>Shigella dysenteriae</i>    |
| <i>Campylobacter jejuni</i>   | <i>Listeria monocytogenes</i>   | <i>Shigella flexneri</i>       |
| <i>Clostridium difficile</i>  | <i>Norovirus</i>                | <i>Shigella sonnei</i>         |
| <i>Cryptosporidium parvum</i> | <i>Salmonella enteritidis</i>   | <i>Staphylococcus aureus</i>   |
| <i>Entamoeba histolytica</i>  | <i>Salmonella paratyphi</i>     | <i>Yersinia enterocolitica</i> |
| <i>Enterovirus</i>            | <i>Salmonella typhimurium</i>   |                                |

## SYMBOLS FOR IVD COMPONENTS & REAGENTS

- |  |                              |  |                                  |
|--|------------------------------|--|----------------------------------|
|  | In vitro diagnostic device   |  | Manufacturer                     |
|  | Consult instructions for use |  | Contains sufficient for <n> test |
|  | Keep dry                     |  | Sample diluent                   |
|  | Use by                       |  | Batch code                       |
|  | Temperature limitation       |  | Catalogue number                 |

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