ENGLISH - ESPAÑOL

Ref: RTB25RA  
RTB25RAC  
RTB25RACN

biotical

# Rotavirus + Adenovirus Card

Rapid Test to detect  
Rotavirus + Adenovirus  
in card format

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biotical health  
ready to test

ISO 13485 CERTIFIED  
MEDICAL DEVICE  
QUALITY MANAGEMENT SYSTEM

# biotical Rotavirus + Adenovirus Card

ENGLISH

## WHAT IS BIOTICAL ROTAVIRUS + ADENOVIRUS CARD?

biotical Rotavirus + Adenovirus card is a colored chromatographic immunoassay test for the simultaneous qualitative detection of Rotavirus and Adenovirus antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of Rotavirus and/or Adenovirus infection.

## GASTROENTERITIS, ROTAVIRUS AND ADENOVIRUS

Rotavirus and Adenovirus are the main cause of severe diarrhea in children under 5 years, causing high morbidity and mortality worldwide.

These viruses are transmitted by feco-oral contact, causing watery diarrhea, vomiting, headache, fever and stomachache. Symptoms usually begin 1 or 2 days after becoming infected and may last 1 to 10 days, depending on the virus (3 days for Rotavirus; 5-8 days for Adenovirus).

The most common procedures to identify a gastroenteritis caused by Rotavirus or Adenovirus involve the isolation and identification of the virus through techniques that require more than 1 day. However, biotical Rotavirus + Adenovirus card, is an immunochromatographic qualitative test, easy to use and with almost immediate results, for the detection of Rotavirus and/or Adenovirus directly from the patient's stool sample, providing results in 10 minutes.


## TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25RA	biotical Rotavirus + Adenovirus Card	25 tests without positive or negative control.
RTB25RAC	biotical Rotavirus + Adenovirus Card (+PC)	25 tests with positive control, without negative control.
RTB25RACN	biotical Rotavirus + Adenovirus Card (+PC +NC)	25 tests with positive and negative control.

## HOW DOES BIOTICAL ROTAVIRUS + ADENOVIRUS CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In the T1 **RED** test line (**Adenovirus**): with mouse monoclonal antibodies against *Adenovirus*.
- In the C **GREEN** line (**Control**): with rabbit polyclonal antibodies against a specific protein.
- In the T2 **RED** test line (**Rotavirus**): with mouse monoclonal antibodies against *Rotavirus*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents of the T1 line (mouse monoclonal antibodies against *Adenovirus*) conjugated with latex of **RED** polystyrene; with another preparation of reagents of the T2 line (mouse monoclonal antibodies against *Rotavirus*) conjugated to **RED** polystyrene latex; and with another preparation for the control solution (specific binding protein) conjugated to **GREEN** polystyrene latex, forming three conjugate colored complexes.

#### TEST LINE (T1/T2):

- **If the sample is positive:** The antigens of the patient's stool sample will react with one (or both) of the red colored conjugate complexes (anti-*Adenovirus* monoclonal antibodies latex red microspheres for the line T1 and anti-*Rotavirus* monoclonal antibodies - red latex microspheres for the line T2) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-*Adenovirus* present in the T1 line and / or the antibodies anti-*Rotavirus* present in the T2 line will capture the complex giving rise in a red line.
- **If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red colored complexes. The anti-*Adenovirus* and anti-*Rotavirus* antibodies present in the T1 and T2 line will not capture the red colored antigen-complex (not formed) and the red lines will not appear.

#### CONTROL LINE (C):

- **Regardless of whether the sample is positive or not**, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

---

## STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

---

## GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.

- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- Used tests should be managed as sanitary waste (sanitary waste container).
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

## MATERIALS

---

### MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Tubes for the sample dilution with sample diluent.
- Positive Control (only RTB25RAC & RTB25RACN references).
- Negative Control (only RTB25RACN reference).

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.


## SPECIMEN COLLECTION AND PREPARATION

---

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20 °C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

### Specimen preparation (see illustration) :

- Take out the cap of the stool collection tube **1** and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **2**, to collect faecal sample (approx. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **3**.

## TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 4.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 5. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)


ADENOVIRUS NEGATIVE  
ROTAVIRUS NEGATIVE

ADENOVIRUS POSITIVE  
ROTAVIRUS NEGATIVE

ADENOVIRUS NEGATIVE  
ROTAVIRUS POSITIVE

ADENOVIRUS POSITIVE  
ROTAVIRUS POSITIVE

INVALID

INVALID

INVALID

INVALID

	biotical Test	Interpretation of the results
1	T1 = NO LINE C = GREEN T2 = NO LINE	<b>NEGATIVE FOR ROTAVIRUS + ADENOVIRUS</b> Only the <b>GREEN</b> line appears across the results window in the site marked with the letter C (control line). No presence of ROTAVIRUS or ADENOVIRUS is detected in the sample.
2	T1 = RED C = GREEN T2 = NO LINE	<b>POSITIVE FOR ADENOVIRUS</b> In addition to the <b>GREEN</b> control line, a <b>RED</b> line also appears in the site marked with the letter T1 (test line 1) in the results window. ADENOVIRUS presence is detected in the sample.
3	T1 = NO LINE C = GREEN T2 = RED	<b>POSITIVE FOR ROTAVIRUS</b> In addition to the <b>GREEN</b> control line, a <b>RED</b> line also appears in the site marked with the letter T2 (test line 2) in the results window. ROTAVIRUS presence is detected in the sample.

	biotical Test	Interpretation of the results
4	T1 = RED C = GREEN T2 = RED	<b>POSITIVE FOR ADENOVIRUS + ROTAVIRUS</b> In addition to the GREEN control line, two RED lines also appears in the site marked with the letter T1 and T2 in the results window. ROTAVIRUS and ADENOVIRUS are detected in the sample.
5-8	Any other result	<b>INVALID</b> Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the tests lines (T1 and T2) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25RAC** & **RTB25RACN**).



## POSITIVE CONTROL: INTENDED USE & COMPOSITION

biotical Rotavirus Positive Control and Adenovirus Positive Control are two external quality controls for Biotical Rotavirus + Adenovirus card test. These control swabs must be used to evaluate that the reagents and the operation of the test are suitable.

biotical Rotavirus Positive Control is manufactured with inactivated Rotavirus antigens dried in swab, containing a preservative. biotical Adenovirus Positive Control is manufactured with inactivated Adenovirus antigens dried in swab, containing a preservative.

## POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1 Add 15 drops of diluent and introduce the positive control swab into the testing tube, 2 rotating 30 seconds and extract the liquid.

3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

## POSITIVE CONTROL – INTERPRETATION OF RESULTS.

biotical Test	Results analysis .
<p>T1 = RED LINE C = GREEN LINE T2 = NO LINE</p>	<p><b>RESULT: Adenovirus POSITIVE CONTROL</b></p> <p>A GREEN(C) line and a RED(T1) line appear in the result window during testing due to the presence of adenovirus antigen.</p>
<p>T1 = NO LINE C = GREEN LINE T2 = RED LINE</p>	<p><b>RESULT: POSITIVE CONTROL of Rotavirus</b></p> <p>A GREEN(C) line and a RED(T2) line appear in the result window during testing due to the presence of rotavirus antigen.</p>
	<p><b>INVALID RESULT</b></p> <p>Total absence of the control line (GREEN), regardless of the appearance or not of any of the T1-T2 test lines (RED), or only the control line appears.</p> <p>Incorrect procedural techniques or reagent deterioration are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.


## NEGATIVE CONTROL: TESTING PROCEDURE.

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- 1 Add 15 drops of diluent and introduce the positive control swab into the testing tube, 2 rotating 30 seconds and extract the liquid.
  - 3 Add 4 drops in the circular window (S).
- **Add 15 drops of diluent** and immediately introduce the negative control swab into the testing tube 1.
  - Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least **30 seconds** 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
  - Remove the Biotical test Card device from its sealed bag just before using it.
  - **Dispense 4 drops** from the testing tube in the corresponding circular window 3. Use a separate pipette and test for each control.
  - **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

biotical Test	Interpretation of the results
 <p>T1 = NO LINE C = <b>GREEN LINE</b> T2 = NO LINE</p>	<p><b>RESULT: CONTROL NEGATIVE.</b></p> <p>Only the control line (C) appears across the results window during the test performance.</p>
	<p><b>RESULT: INVALID FOR NEGATIVE CONTROL.</b></p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

## LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The samples must be taken during the first week of appearance of the symptoms to ensure the greatest presence of pathogens in the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of Rotavirus/Adenovirus in the sample. To confirm the presence of infection by Rotavirus/Adenovirus it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques (biochemical methods and microscopy).
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory (biochemical methods and microscopy). A negative result does not indicate absence of infection by Rotavirus/Adenovirus. It is possible to appear false negatives by not reaching the detection limit of the test.

## PERFORMANCE CHARACTERISTICS

### Clinical sensitivity and specificity:

An evaluation, with faecal samples, was performed using an immunochromatographic test (biotical Rotavirus+Adenovirus, biotical) and these results were compared with a commercial available test (Ridascreen® Rotavirus ELISA Test, r-Biopharm) for T2 line (rotavirus) and T1 line (adenovirus) were confirmed using PCR.

The results were as follows:

		Ridascreen® Rotavirus ELISA Test (ELISA)		
		+	-	Total
biotical Rotavirus + Adenovirus card (T2: ROTAVIRUS)	+	18	1	19
	-	0	43	43
	Total	18	44	62

Biotal Rotavirus + Adenovirus Card vs Ridascreen® Rotavirus ELISA Test (ROTA VIRUS)	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99%	81.5% - 100%
SPECIFICITY	97.7%	87.9% - 99.9%
PPV	94.7%	72.2% - 99.2%
NPV	>99%	-
ACCURACY	98.4%	91.3% - 99.9%

		Adenovirus PCR (PCR)		
		+	-	Total
biotal Rotavirus + Adenovirus card (T1: ADENOVIRUS)	+	7	0	7
	-	0	52	52
	Total	7	52	59

Biotal Rotavirus + Adenovirus Card vs Adenovirus PCR (ADENOVIRUS)	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99%	59.0% - 100%
SPECIFICITY	>99%	93.2% - 100%
PPV	>99%	-
NPV	>99%	-
ACCURACY	>99%	93.9% - 100%

The results showed a high sensitivity and specificity to detect Rotavirus + Adenovirus using biotal Rotavirus + Adenovirus Card.

## CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotal Rotavirus + Adenovirus, no cross reactivity against gastrointestinal pathogens occasionally present in faeces.

<i>Adenovirus (for T2)</i>	<i>Escherichia coli O157:H7</i>	<i>Salmonella typhimurium</i>
<i>Astrovirus</i>	<i>Giardia Lambia</i>	<i>Salmonella typhi</i>
<i>Campylobacter coli</i>	<i>Helicobacter pylori</i>	<i>Shigella boydii</i>
<i>Campylobacter jejuni</i>	<i>Listeria Monocytogens</i>	<i>Shigella dysenteriae</i>
<i>Clostridium difficile</i>	<i>Norovirus</i>	<i>Shigella flexneri</i>
<i>Cryptosporidium parvum</i>	<i>Rotavirus (for T1)</i>	<i>Shigella sonnei</i>
<i>Enterovirus</i>	<i>Salmonella enteritidis</i>	<i>Streptococcus aureus</i>
<i>Entamoeba Hystolitica</i>	<i>Salmonella paratyphi</i>	<i>Yersinia Enterocolitica</i>

## SYMBOLS FOR IVD COMPONENTS & REAGENTS

---

In vitro diagnostic device

Consult instructions for use

Keep dry

Use by

Temperature limitation

Manufacturer

Contains sufficient for <n> test

DIL

Sample diluent

LOT

Batch code

REF

Catalogue number

\*\*\*

# biotical Rotavirus + Adenovirus Card

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