ENGLISH - ESPAÑOL

Ref: RTB25AT
RTB25ATC
RTB25ATCN

biotical Astrovirus Card

Rapid Test to detect
Astrovirus in card
format

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biotical health

ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical Astrovirus Card

ENGLISH

WHAT IS BIOTICAL ASTROVIRUS CARD?

biotical Astrovirus card is a coloured chromatographic immunoassay test for the qualitative detection of *Astrovirus* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of *Astrovirus* infection.

GASTROENTERITIS AND ASTROVIRUS

Acute gastroenteritis is a very common disease that affects children and young people worldwide, with *Astrovirus* infection being one of the main causal agents in children under 5 years of age, causing high morbidity and mortality worldwide.

Astrovirus is transmitted by faecal-oral contact, causing watery diarrhoea, vomiting, headache, fever and stomach-ache. Symptoms usually begin 1 or 2 days after becoming infected and may last 1 to 3 days.

The most common procedures to identify a gastroenteritis caused by *Astrovirus* involve the isolation and identification of the virus through techniques that require more than 1 day. However, biotical Astrovirus card, is an immunochromatographic qualitative test, easy to use and with almost immediate results, for the detection of Astrovirus directly from the patient’s stool sample, providing results in 10 minutes.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25AT	biotical Astrovirus Card	25 tests without positive or negative controls.
RTB25ATC	biotical Astrovirus Card (+PC)	25 tests with positive control, without negative control.
RTB25ATCN	biotical Astrovirus Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL ASTROVIRUS CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against *Astrovirus*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against *Astrovirus*) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to red polystyrene latex.

TEST LINES (T):

If the sample is positive: The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-*Astrovirus* monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-*Astrovirus* present in the T line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-*Astrovirus* antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.

- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25ATC & RTB25ATCN references).
- Negative Control (only RTB25ATCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.


The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.


Specimen preparation (see illustration):

- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample (aprox. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.


- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>Astrovirus</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>Astrovirus</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25ATC** & **RTB25ATCN**).

POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical Astrovirus Positive Control is an external quality control for biotical ASTROVIRUS card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical Astrovirus Positive Control: inactivated *Astrovirus* antigens dried in swab containing a preservative.

2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>
---	--	---

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The samples must be taken during the first week of appearance of the symptoms to ensure the presence of pathogens in the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of Astrovirus in the sample. To confirm the presence of infection by Astrovirus it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques (biochemical methods and microscopy).
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory techniques (biochemical methods and microscopy). A negative result does not indicate absence of infection by Astrovirus. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotal Astrovirus Card, biotal) and another commercial ELISA Astrovirus assay (Ridascreen® Astrovirus ELISA Test, r-Biopharm). The samples were taken from patients with the same as Astrovirus infection symptoms arrived in Gastroenterology area in a Spanish hospital.

The results were as follows:

		Ridascreen® Astrovirus ELISA Test		
		+	-	Total
Biotical Astrovirus Card	+	16	0	16
	-	1	11	12
	Total	17	11	28

<i>Biotical Astrovirus Card vs Ridascreen® Astrovirus ELISA Test</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	94.12%	71.31% - 99.85%
SPECIFICITY	100.00%	71.51% - 100.00%
PPV	100.00%	-
NPV	91.67%	62.16% - 98.66%
ACCURACY	96.43%	81.65% - 99.91%

The results showed a high sensitivity and specificity to detect Astrovirus using biotical Astrovirus Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Astrovirus Card; no cross reactivity against gastrointestinal organism occasionally present in faeces.

Adenovirus
Campylobacter coli
Campylobacter jejuni
Clostridium difficile
Cryptosporidium parvum
Enterovirus
Entamoeba histolytica
Escherichia coli O157:H7

Giardia lamblia
Helicobacter pylori
Listeria monocytogenes
Norovirus
Rotavirus
Salmonella enteritidis
Salmonella paratyphi
Salmonella typhi

Salmonella typhimurium
Shigella boydii
Shigella dysenteriae
Shigella flexneri
Shigella sonnei
Staphylococcus aureus
Yersinia enterocolitica

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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