ENGLISH - ESPAÑOL

Ref: RTB25NV
RTB25NVC
RTB25NVCN

biotical Norovirus Gen I + II Card

Rapid Test to detect
Norovirus Gen I + II in
card format

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biotical health

ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical Norovirus Gen I+II Card

ENGLISH

WHAT IS BIOTICAL NOROVIRUS GEN I+II CARD?

biotical Norovirus Gen I+II card is a rapid immunochromatography test for the qualitative simultaneous detection of *Norovirus genogroups I and II (Gen I and Gen II)* in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional is able to make a presumptive diagnosis of *Norovirus* infection.

ACUTE GASTROENTERITIS AND NOROVIRUS

Noroviruses are responsible for the majority of outbreaks of food poisoning gastroenteritis in the world. Gen I and Gen II are the genogroups that infect humans causing acute gastroenteritis symptoms: stomach pain, nausea, diarrhoea and vomiting. Usually, these symptoms appear 12 to 48 hours after being exposed to the virus and disappear 1 to 3 days later.

Norovirus affects people of all ages, and it could be severe, especially in children and seniors. It can be transmitted through the consumption of food or contaminated water, food handled by an infected person, through the contact with contaminated surfaces, or by person-to-person contact. Outbreaks often occur in closed or semi-closed communities, such as hospitals, prisons, residences and cruises.

Norovirus infection is highly contagious, so a quick diagnosis helps prevent the spread of infection. biotical Norovirus Gen I+II Card is an immunochromatographic qualitative test, with an easy use and with almost immediate results, for the detection of *Norovirus Genotypes I and II* directly from the patient's stool sample, providing results in 10 minutes.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25NV	biotical Norovirus Gen I+II Card	25 tests without positive or negative controls.
RTB25NVC	biotical Norovirus Gen I+II Card (+PC)	25 tests with positive control, without negative control.
RTB25NVCN	biotical Norovirus Gen I+II Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL NOROVIRUS CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In the **RED** T₁ line (Test *Norovirus Gen I*): with mouse monoclonal antibodies against *Norovirus Gen I*.
- In the **GREEN** C line (Control): with rabbit polyclonal antibodies against a specific protein.
- In the **RED** T₂ line (Test *Norovirus Gen II*): with mouse monoclonal antibodies against *Norovirus Gen II*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents of the T1 line (mouse monoclonal antibodies against *Norovirus Gen I*) conjugated with latex of **RED** polystyrene; with another preparation of reagents of the T2 line (mouse monoclonal antibodies against *Norovirus Gen II*) conjugated to **RED** polystyrene latex; and with another preparation for the control solution (specific binding protein) conjugated to **GREEN** polystyrene latex, forming three conjugate coloured complexes.

TEST LINES (T1 / T2):

- **If the sample is positive:** The antigens of the patient's stool sample will react with one (or both) of the red coloured conjugate complexes (anti-*Norovirus Gen I* monoclonal antibodies - latex red microspheres for the line T1 and anti-*Norovirus Gen II* monoclonal antibodies - red latex microspheres for the line T2) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti- *Norovirus Gen I* present in the T1 line and / or the antibodies anti-*Norovirus Gen II* present in the T2 line will capture the complex giving rise in a red line.
- **If the sample is negative:** If the sample is negative and therefore does not contain antigens (or these are present in a concentration below the detection limit), there will be no reaction with the red coloured complexes. The anti-*Norovirus Gen I* and anti-*Norovirus Gen II* antibodies present in the T1 and T2 lines will not capture the coloured antigen-complexes and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25NVC & RTB25NVCN references).
- Negative Control (only RTB25NVCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube **1** and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **2**, to collect faecal sample and add it to the stool collection tube (125mg aprox). For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **3**. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

1 Collection tube


2 Stool sample

3 Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

4 Cut the end of the cap.

5 Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 4.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 5. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

Norovirus Gen I **NEGATIVE**
Norovirus Gen II **NEGATIVE**

Norovirus Gen I **POSITIVE**
Norovirus Gen II **NEGATIVE**

Norovirus Gen I **NEGATIVE**
Norovirus Gen II **POSITIVE**

Norovirus Gen I **POSITIVE**
Norovirus Gen II **POSITIVE**

INVALID

INVALID

INVALID

INVALID

	biotical Test	Interpretation of the results
1	T1 = NO LINE C = GREEN T2 = NO LINE	NEGATIVE FOR NOROVIRUS GEN I + II Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>Norovirus Gen I</i> or <i>II</i> .
2	T1 = RED C = GREEN T2 = NO LINE	POSITIVE FOR NOROVIRUS GEN I In addition to the GREEN control line, a RED line also appears in the site marked with the letter T1 (test line) in the results window. <i>Norovirus Gen I</i> presence is detected in the sample. No presence of <i>Norovirus Gen II</i> .
3	T1 = NO LINE C = GREEN T2 = RED	POSITIVE FOR NOROVIRUS GEN II In addition to the GREEN control line, a RED line also appears in the site marked with the letter T2 (test line) in the results window. <i>Norovirus Gen II</i> presence is detected in the sample. No presence of <i>Norovirus Gen I</i> .

4	T1 = RED C = GREEN T2 = RED	<p style="text-align: center;">POSITIVE FOR NOROVIRUS GEN I + II</p> In addition to the GREEN control line, two RED lines also appears in the site marked with the letter T1 and T2 in the results window. <i>Norovirus Gen I and II</i> are detected in the sample.
5-8	Any other result	<p style="text-align: center;">INVALID</p> Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the tests lines (T1/T2) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25NVC** & **RTB25NVCN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

biotical Norovirus Gen I Positive control and biotical Norovirus Gen II Positive control are an external quality control for Norovirus genogroup I (Gen I) and genogroup II (Gen II) tests lines in biotical products. Use those controls swabs to check that the extraction reagents and the test are working properly.

biotical Norovirus Gen I Positive control: inactivated *Norovirus Gen I* antigens dried in swab containing a preservative. biotical Norovirus Gen II Positive control: inactivated *Norovirus Gen II* antigens dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biological device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

biological Test	Interpretation of the results
<p>T1 = RED LINE C = GREEN LINE T2 = NO LINE</p>	<p>RESULT: Norovirus Gen I POSITIVE CONTROL</p> <p>A GREEN line (C) and a RED line (T1) appears.</p>
<p>T1 = NO LINE C = GREEN LINE T2 = RED LINE</p>	<p>RESULT: Norovirus Gen II POSITIVE CONTROL.</p> <p>A GREEN line (C) and a RED line (T2) appears.</p>
	<p>INVALID RESULTS FOR POSITIVE CONTROL.</p> <p>Total absence of the control line, regardless the appearance or not of any of the test lines, or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE


Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube 1.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

biotal Test	Interpretation of the results
 <p>T1 = NO LINE C = GREEN LINE T2 = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>

RESULT: INVALID FOR NEGATIVE CONTROL

Total absence of the control line, and/or the card shows the test line.

Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new negative control swab.

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Norovirus* in the sample. To confirm the presence of infection by *Norovirus* it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques (biochemical methods and PCR).
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory techniques (biochemical methods and PCR). A negative result does not indicate absence of infection by *Norovirus*. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using an immunochromatographic test (biotical Norovirus Gen I+II Card, biotical) and these results were compared with a commercially available immunochromatographic test (Simple Norovirus, Operon) and confirmed by PCR.

The results were as follows:

		IC test: Simple Norovirus		
		+	-	Total
IC test: biotical Norovirus Gen I+II card (Norovirus Gen I)	+	2	0	2
	-	0	48	48
	Total	2	48	50

<i>Biotical Norovirus Gen I + II Card vs Simple Norovirus (Norovirus Gen I)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99.9%	15.81% - 100.00%
SPECIFICITY	>99.9%	92.60% - 100.00%
PPV	>99.9%	-
NPV	>99.9%	-
ACCURACY	>99.9%	92.89% - 100.00%

		PCR		
		+	-	Total
IC test: biotical Norovirus Gen I+II card (Norovirus Gen I)	+	2	0	2
	-	0	48	48
	Total	2	48	50

<i>Biotical Norovirus Gen I + II Card vs PCR (Norovirus Gen I)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99.9%	15.81% - 100.00%
SPECIFICITY	>99.9%	92.60% - 100.00%
PPV	>99.9%	-
NPV	>99.9%	-
ACCURACY	>99.9%	92.89% - 100.00%

		IC test: Simple Norovirus		
		+	-	Total
IC test: biotal Norovirus Gen I+II card (Norovirus Gen II)	+	10	0	10
	-	0	48	48
	Total	10	48	58

Biotal Norovirus Gen I + II Card vs Simple Norovirus (Norovirus Gen II)	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99.9%	69.15% - 100.00%
SPECIFICITY	>99.9%	92.60% - 100.00%
PPV	>99.9%	-
NPV	>99.9%	-
ACCURACY	>99.9%	93.84% - 100.00%

		PCR		
		+	-	Total
IC test: biotal Norovirus Gen I+II card (Norovirus Gen II)	+	8	0	8
	-	2	48	50
	Total	10	48	58

Biotal Norovirus Gen I + II Card vs PCR (Norovirus Gen II)	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	80.0%	44.39% to 97.48%
SPECIFICITY	>99.9%	92.60% to 100.00%
PPV	>99.9%	-
NPV	96.0%	87.42% to 98.81%
ACCURACY	96.6%	88.09% to 99.58%

The results showed a high sensitivity and specificity to detect *genogroup I and II Noroviruses* using biotal Norovirus Gen I+II Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Norovirus GI+GII Card, no cross reactivity against gastrointestinal organism occasionally present in faeces.

<i>Adenovirus</i>	<i>Hepatitis A</i>	<i>Salmonella paratyphi</i>
<i>Astrovirus</i>	<i>Giardia lamblia</i>	<i>Salmonella typhimurium</i>
<i>Campylobacter coli</i>	<i>Helicobacter pylori</i>	<i>Salmonella typhi</i>
<i>Campylobacter jejuni</i>	<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>
<i>Clostridium difficile</i>	<i>Norovirus GI (T2 line)</i>	<i>Shigella dysenteriae</i>
<i>Cryptosporidium parvum</i>	<i>Norovirus GII (T1 line)</i>	<i>Shigella flexneri</i>
<i>Enterovirus</i>	<i>Rotavirus</i>	<i>Shigella sonnei</i>
<i>Escherichia coli O111</i>	<i>RSV</i>	<i>Staphylococcus aureus</i>

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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