ENGLISH - ESPAÑOL

Ref: RTB25EV
RTB25EVC
RTB25EVCN

biotical Enterovirus Card

Rapid Test to detect
Enterovirus in card
format

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25 UN.

IU RTB25EV-RT25EVC-RTB25EVCN Rev. 01 / 2024-02

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical Enterovirus Card

ENGLISH

WHAT IS BIOTICAL ENTEROVIRUS CARD?

biotical *Enterovirus* card is a coloured chromatographic immunoassay test for the qualitative detection of *Enterovirus* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of *Enterovirus* infection.

INTRODUCTION

Enteroviruses belong to the Picornaviridae family and are the smallest known non-enveloped viruses capable of infecting both humans and animals. They are single-stranded RNA viruses.

It can be transmitted by the faecal-oral route and by mouth-to-mouth after contact with an infected person, or by direct contact with skin or ophthalmic lesions. Another possible route of infection, which usually causes epidemic outbreaks, is through contact of water, food or soil with contaminated faeces. The average incubation period for enteroviral infection can be from 3-10 to 30 days.

Once past the gastrointestinal barrier, the virus spreads through the bloodstream to every organ in the body, especially the heart, skin, and central nervous system. After infection, an infected person is capable of excreting large amounts of virus together with their faeces for up to 16 weeks.


The traditional division classifies this taxonomic group within the subgenus Poliovirus, Coxsackievirus (group A, B), Echovirus and a group of *Enteroviruses* named according to their serotype number (66-71, and the recently identified 73-75, 77,78). In 2003, the International Committee on Taxonomy of Viruses created a new taxonomic classification. They divided *Enteroviruses* into 5 species groups based on their molecular properties (HEV-A, HEV-B, HEV-C, HEV-D, and Poliovirus.).

TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25EV	biotical Enterovirus Card	25 test without positive or negative controls.
RTB25EVC	biotical Enterovirus Card (+PC)	25 test with positive control, without negative control.
RTB25EVCN	biotical Enterovirus Card (+PC +NC)	25 test with positive and negative control.

HOW DOES BIOTICAL ENTEROVIRUS CARD WORK?

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against *Enterovirus*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against *Enterovirus*) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to red polystyrene latex.

TEST LINES:

If the sample is positive: The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-*Enterovirus* monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-*Enterovirus* present in the T line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-*Enterovirus* antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25EVC & RTB25EVCN references).
- Negative Control (only RTB25EVCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.


Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.

Specimen preparation (see illustration):


- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

Collection tube

1 Stool sample

2 Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

3 Cut the end of the cap.


4 Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID RESULT

INVALID RESULT

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>Enterovirus</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>Enterovirus</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25EVC** & **RTB25EVCN**).

POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical *Enterovirus* Positive Control is an external quality control for biotical *Enterovirus* card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical *Enterovirus* Positive Control: inactivated *Enterovirus* antigens dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- 1 Add 15 drops of diluent and introduce the positive control swab into the testing tube.
- 2 Rotating 30 seconds and extract the liquid.
- 3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL - INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

Biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>

2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>
---	---	---

LIMITATIONS OF THIS TEST

- The test must be carried out within 2 hours after opening the sealed bag.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- After one month of infection, the number of viruses in faeces is decreasing, making the sample less reactive. Stool samples could be collected previously to the onset of symptoms or also at 24-48 hours.
- If the patient has been recently vaccinated (for example against Poliovirus), it could appear a positive result.
- biotical *Enterovirus* Card should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of *Enterovirus* in faecal samples; nevertheless, it can be due to several species. A positive result should be followed up with diagnostic procedures to confirm the results and the species. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation still persist, an *Enterovirus* determination should be carried out with another technique (for example microscopy, biochemical method).

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical Enterovirus Card, biotical) and another commercial *Enterovirus* assay (IDEIA *Enterovirus* assay, Dako and IMAGEN™ *Enterovirus*, Oxoid).

The results were as follows:

		IDEIA <i>Enterovirus</i> assay & IMAGEN™ <i>Enterovirus</i>		
		+	-	Total
Biotical Enterovirus Card	+	3	0	3
	-	0	32	32
	Total	3	32	35

<i>Biotical Enterovirus Card vs IDEIA Enterovirus assay, Dako y IMAGEN™ Enterovirus, Oxoid</i>	MIDDLE VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99.9%	29.2% - 100%
SPECIFICITY	>99.9%	89.1% - 100%
PPV	>99.9%	-
NPV	>99.9%	-
PRECISION	>99.9%	90.0% - 100%

The results showed a high sensitivity and specificity to detect *Enterovirus* using biotical Enterovirus Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Enterovirus Card, no cross reactivity against gastrointestinal organism occasionally present in faeces.

Adenovirus
Campylobacter coli
Campylobacter jejuni
Clostridium difficile
Cryptosporidium parvum
Entamoeba histolytica
Escherichia coli O157:H7
Giardia lamblia

Helicobacter pylori
Listeria monocytogenes
Norovirus
Rotavirus
Salmonella enteritidis
Salmonella paratyphi
Salmonella typhimurium
Salmonella typhi

Shigella boydii
Shigella dysenteriae
Shigella flexneri
Shigella sonnei
Staphylococcus aureus
Yersinia enterocolitica

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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REF

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REFERENCES / BIBLIOGRAFÍA /

- 1- RAJTAR, B. et al. "Enteroviruses in water environment-a potential threat to public health".Ann. Agric. Environ. Med., 2008, No 15, p. 199-203.
- 2- BAEK, K. et al. "Molecular and Epidemiological Characterization of Enteroviruses Isolated in Chungnam, Korea from 2005-2006." J. Microbiol. Biotechnol., 2009, 19(9), p. 1055-1064.
- 3- FONG, T. et al. "Enteric Viruses of Humans and Animals in Aquatic Environments: Health Risks, Detection, and Potential Water Quality Assessment Tools". Microbiology and Molecular Biology Reviews, June 2005, Vol. 69, No. 2: p. 357-371.
- 4- AFFFI, S. et al. "Isolation and Identification of Non-Polio Enteroviruses from Children in Different Egyptian Governorates", Australian Journal of Basic and Applied Sciences, 2009, Vol. 3, No. 4: pp. 3230-3238.

Sierra de Guadarrama, 1. Ed. Dublín
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es

biotical health S.L.U.