ENGLISH - ESPAÑOL

Ref: RTB25PG
RTB25PGC
RTB25PGCN

biotical

H. pylori Card

Rapid Test to detect
Helicobacter pylori in
card format

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biotical health
ready to test

IU RTB25PG-RTB25PGC-RTB25PGCN Rev. 05 / 2024-04

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical H. pylori Card

ENGLISH

WHAT IS BIOTICAL H. PYLORI CARD?

biotical H. pylori card is a rapid immunochromatography test for the qualitative detection of *Helicobacter pylori* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional can make a presumptive diagnosis to detect *Helicobacter pylori*.

WHAT IS HELICOBACTER PYLORI (H. PYLORI)?

Helicobacter pylori is a spiral-shaped bacterium that has bacterial flagella at one end, which gives it the mobility to colonize the stomach. It is one of the most prevalent human infections, affecting approximately 60% of the population, although in populations with lower socioeconomic status and worse health conditions, it reaches 80%.

The infection can develop in an inflammatory response in the stomach and chronic gastritis. Although only a small percentage of patients will develop symptoms despite suffering from gastritis, it seems to depend on both host factors and the strain of H. pylori which causes the infection.

In addition to causing chronic gastritis, this bacterium is related to gastric and duodenal ulcer, lymphoma and gastric cancer. H. pylori can be transmitted from one person to another by direct contact with saliva, vomit or stool. It can also be transmitted through contaminated food or water.

biotical H. pylori Card is a qualitative immunochromatographic test for the detection of H. pylori antigens in stool samples. The result is obtained in 10 minutes.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25PG	biotical H. pylori Card	25 Tests without positive or negative controls.
RTB25PGC	biotical H. pylori Card (+PC)	25 Tests with positive control, without negative control.
RTB25PGCN	biotical H. pylori Card (+PC +NC)	25 Tests with positive and negative control.

HOW DOES BIOTICAL H. PYLORI CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- **In RED line (C):** with rabbit polyclonal antibodies against a specific protein.
- **In RED line (T):** with mouse monoclonal antibodies against H. pylori.
- **Sample area (S):** the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against H. pylori) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to red polystyrene latex.

TEST LINES:

If the sample is positive: The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-H. pylori monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-H. pylori present in the T line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-H. pylori antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the red control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.

- Used tests should be managed as sanitary waste (sanitary waste container).
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25PGC & RTB25PGCN references).
- Negative Control (only RTB25PGCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION

Stool samples.


Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube **1** and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **2**, making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample (approx. 50mg) and add it to the collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **3**.

1 Collection tube


2 Stool sample


3 Shake

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

4 Cut the end of the cap.


5 Add 3 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 4.
- **Dispense 3 drops** from the testing tube in the circular window marked with the letter S 5. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID RESULT

INVALID RESULT

	biotical Test	Interpretation of the results
1	C = RED LINE T = NO LINE	NEGATIVE Only the C line appears across the results window in the site marked. No presence of <i>Helicobacter pylori</i> is detected in the sample.
2	C = RED LINE T = RED LINE	POSITIVE In addition to C line, T line also appears in the results window. <i>Helicobacter pylori</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control C line (regardless the appearance or not of T line) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A red line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25PGC** & **RTB25PGCN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotal H. pylori Positive Control is an external quality control for biotal H. pylori card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotal H. pylori Positive Control: inactivated H. pylori antigens dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by **rotating the positive control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes**. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = RED LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>C line and T line appears across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Presence of only one line or none of them.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the Biotical test Card device from its sealed bag just before using it.
- Dispense 4 drops from the testing tube in the corresponding circular window ③. Use a separate pipette and test for each control.
- Read the result at 10 minutes. Do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = RED LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Helicobacter pylori* in the sample. To confirm the presence of infection by *Helicobacter pylori* it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other invasive techniques (endoscopy).
- False positives may occur in the test with bloody samples or mucous stools, so the use of the test in this class of samples is not recommended.
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory or invasive techniques for confirmation. A negative result does not indicate absence of infection by *Helicobacter pylori*. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit):

Detection limit range is: 0.78ng/mL – 0.09ng/mL of H. pylori recombinant outer membrane protein.

Clinical sensitivity and specificity.

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical H. pylori Card, biotical) and qPCR technique (VIASURE *Helicobacter pylori* Real Time Detection kit, CerTest). The samples were taken from patients with the same as *Helicobacter pylori* infection symptoms arrived in Gastroenterology area in a Spanish hospital.

The results were as follows:

		qPCR test: Viasure Helicobacter pylori Real Time Detection kit		
		+	-	Total
Ic test: biotical H. pylori card	+	54	1	55
	-	1	60	61
	Total	55	61	116

<i>biotical H. pylori Card vs Viasure Helicobacter pylori Real Time Detection Kit</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	98,2%	90,3 - 100%
SPECIFICITY	98,4%	91,2 - 100%
PPV	98,2%	90,3 - 100%
NPV	98,4%	91,2 - 100%
ACCURACY	98,3%	93,9 - 99,8%

The results showed a high sensitivity and specificity to detect *Helicobacter pylori* using biotical H. pylori Card.

CROSS REACTIVITY

An evaluation was performed to determine possible biotical H. pylori Card cross-reactions; no cross-reaction was detected with gastrointestinal pathogens occasionally present in faeces.

Adenovirus
Astrovirus
 Human calprotectin
Campylobacter coli
Clostridium difficile GDH
Clostridium difficile Toxin A
Clostridium difficile Toxin B
Clostridium perfringens
Entamoeba histolytica
Escherichia coli O111
Escherichia coli O26
Escherichia coli O157
Giardia lamblia α -1 giardina

Giardia lamblia CPW1
 Bovine haemoglobin
 Pig haemoglobin
 Human haemoglobin
 Human lactoferrin
Legionella pneumophila
Listeria monocytogenes
Norovirus GI.1
Norovirus GI.4
Rotavirus
Salmonella enteritidis
Salmonella paratyphi
Salmonella typhi

Salmonella typhimurium
Staphylococcus aureus
Shigella boydii
Shigella dysenteriae
Shigella flexneri
Shigella sonnei
Streptococcus pneumococcal
Streptococcus pyogenes
 Human transferrin
Yersinia Enterocolitica O:3
Yersinia Enterocolitica O:9

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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