ENGLISH - ESPAÑOL

Ref: RTB25GD
RTB25GDC
RTB25GDCN

biotical

C. difficile GDH Card

Rapid Test to detect
Clostridium difficile
GDH in card format

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biotical health
ready to test

IU RTB25GD-RTB25GDC-RTB25GDCN Rev. 04 / 2024-03

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical C. difficile GDH Card

ENGLISH

WHAT IS BIOTICAL C. DIFFICILE GDH CARD USED FOR?

biotical C. difficile GDH Card is a rapid immunochromatography test for the qualitative detection of the enzyme glutamate dehydrogenase (GDH) *Clostridium difficile* in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a simple diagnosis and detect the possible presence of *Clostridium difficile* in infected people.

WHAT IS CLOSTRIDIUM DIFFICILE (C. DIFFICILE)?

Clostridium difficile (*C. difficile*) is an anaerobic Gram-positive bacterium. It is the leading pathogen linked to antibiotic-associated diarrhoea and/or pseudomembranous colitis in hospitalized patients, the leading cause of healthcare-associated diarrhoea in developed countries, and the leading source of nosocomial disease morbidity and mortality worldwide.

Most of the existing bacterial flora in the colon of a healthy adult is generally resistant to colonization by *C. difficile*. If the bacterial flora is altered (due to the excessive use of antibiotics, for example), resistance to colonization by these bacteria is lost and the microorganism is allowed to multiply in the colon. In addition, its ability to form spores allows this bacterium to persist in patients for long periods and facilitates its faecal-oral transmission.

This colonization causes the release of:

- **Toxins A and B:** toxigenic strains responsible for the clinical manifestations of *C. difficile* disease, which can range from non-serious watery diarrhoea, pseudomembranous colitis, toxic megacolon, to death (see biotical C. difficile Toxins A+B Card).
- **Clostridium difficile glutamate dehydrogenase (GDH):** enzyme produced in large quantities by toxigenic (toxins A and B) and non-toxigenic strains, which makes it an excellent marker to determine the presence of the microorganism.
- **Biotical Clostridium difficile GDH Card** will detect the presence of this enzyme.

Toxigenic culture (TC) is the reference technique used to determine *Clostridium difficile* infection. This method consists of culturing and isolating *C. difficile* from a stool sample, followed by identification of the toxin, which requires very laborious work to obtain a result. However, the biotical C. difficile GDH card is an easy-to-use qualitative immunochromatographic test with almost immediate results for the detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH) in stool samples.


REFERENCE TABLE

In the following table you can find all the product references related to the test purchased.

REFERENCE	NAME	CONTENT
RTB25GD	biotal C. difficile GDH Card	25 Test without positive or negative control.
RTB25GDC	biotal C. difficile GDH Card (+PC)	25 Test with positive control.
RTB25GDCN	biotal C. difficile GDH Card (+PC +NC)	25 Test with positive and negative control.

HOW DOES BIOTAL C. DIFFICILE GDH CARD WORK?

The cassette (card) has inside a test strip with a nitrocellulose membrane that has been previously fixed:

- In the **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In the **RED** line (T): with mouse monoclonal antibodies against antigens of GDH.
- In the **sample area (S)**: a reagent preparation for test line T (mouse monoclonal antibodies against GDH) conjugated with red polystyrene latex and another preparation for control line C (specific binding protein) conjugated with green polystyrene latex have been dispensed, forming two conjugated coloured complexes.

TEST LINE (T):

If the sample is positive: Antigens from the patient sample will react with the red coloured conjugate complex (anti-GDH monoclonal antibodies-red latex microspheres) and progress through the membrane by capillarity. The anti-GDH antibodies present on the membrane (test line) will capture the coloured complex of the test, giving rise to the red line.

If the sample is negative: if the sample is negative and therefore does not contain antigen (or the antigen is present at a concentration below the limit of detection), there is no reaction with the red coloured complex. Anti-GDH antibodies present on the membrane (test line) will not capture the red coloured antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or negative, to verify that the test has been carried out correctly (sufficient sample volume, appropriate flow rate, internal control of the reagents, etc.) the mixture will continue to move through the membrane towards the immobilized antibodies against the specific protein located on the line of control (C). These specific anti-protein antibodies present on the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE CONDITIONS

The product must be stored between 2 and 30°C in its original sealed container, to achieve optimum performance until the expiration date printed on the container. It should not be opened until the time of use. do not freeze.

GENERAL PRECAUTIONS

- For professional in vitro use only.
- Do not use after the expiration date.
- Samples must be considered potentially dangerous and must be handled in the same way as an infectious agent. A new test should be used for each sample to avoid contamination errors.
- Tests and samples must be disposed of in an appropriate medical waste container after use.
- The reagents contain preservatives. Avoid any contact with the skin or mucous membranes. Consult the safety data sheet, available on request.
- The components included in the kit are approved for use in the test as supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow good laboratory practices, wear protective clothing, wear disposal gloves, goggles, and mask. Do not eat, drink or smoke in the work area.
- Pale yellow or orange bands may appear in the result window of the test strips, at the positions where the coloured test and control bands should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or card (Card) of the test.
- Instructions for use.
- Tubes for sample dilution with diluent.
- Positive control (only in the reference RTB25GDC and RTB25GDCN).
- Negative control (only in the reference RTB25GDCN).

MATERIALS REQUIRED BUT NOT PROVIDED

- Container for sample collection.
- Disposable gloves.
- Stopwatch.

Kit contains materials for analysis of up to 25 samples or controls.

SAMPLE COLLECTION AND PREPARATION

Samples should be collected in a clean container. Samples must be kept cold (2-8°C), for a maximum of 7 days, until they are used. To preserve the samples for a long time, maximum 1 year, they must be kept frozen at -20°C. In this case, the sample must be completely thawed and brought to room temperature before it can be used in the test. Freeze and thaw cycles are not recommended, so only the amount of sample needed to perform the test should be thawed. Homogenize the sample vigorously before its preparation.

Sample preparation (see drawing):

- Open the sample dilution tube. ① and with the help of the stick take a sufficient amount of sample of the collected feces. To do this, the stick will be inserted once in 4 different areas of the sample ②, making sure that, at each insertion, the helix of the stick is covered with the sample to collect the appropriate amount of feces and then the sample will be introduced into the tube for sample dilution. Adding an additional amount of sample that exceeds the helix of the stick could cause erroneous results. For liquid samples, add approx. 125µL in the tube for sample dilution using a micropipette.
- Close the tube containing the sample and diluent. Shake it to facilitate the dispersion of the sample ③.

1 Sample dilution tube


2 Stool sample


3 shake vigorously

TEST PROCEDURE

Previously, the tests, stool samples and controls must be conditioned at room temperature (15-30 °C). do not open the containers until the moment of the test.

4 Cut the end of the cap.


5 Add 3 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 4.
- **Dispense 3 drops** from the testing tube in the circular window marked with the letter S 5. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID RESULT

INVALID RESULT

	biotal Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>C. difficile GDH</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>C. difficile GDH</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25GDC** & **RTB25GDCN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

biotal C. difficile GDH Positive control is external quality control for biotal C. difficile GDH Card. This control swab should be used to assess that the reagents and performance of the test are adequate..

biotal C. difficile GDH Positive control is made from swab-dried *C. difficile* Glutamate Dehydrogenase (GDH) inactivated antigens, contains a preservative. Not infective.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by **rotating the positive control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes**. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotal Test	Interpretation of the results
1	 <p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- 1 Add 15 drops of diluent and introduce the negative control swab into the testing tube.
- 2 Rotating 30 seconds and extract the liquid.
- 3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotal Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test must be performed a maximum of 2 hours after opening the bag that contains it.
- Indistinct brownish lines may appear. They are caused by an insufficient dilution of the sample and have no diagnostic value. If this occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other kinds of samples has not been tested, so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Clostridium difficile* in the sample, without determining whether they are toxigenic or non-toxigenic strains. To confirm *Clostridium difficile* infection, it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis that may include the toxicity of the strain found in the sample, it is necessary to confirm the test results using other diagnostic techniques (toxigenic culture).
- If a negative result is obtained in the test and the symptoms persist, it is advisable to repeat the test after a few days. A negative result does not indicate the absence of *Clostridium difficile* infection. The appearance of false negatives is possible due to not reaching the detection limit of the test.
- Bloody stool samples contain components that could cause non-specific reactions in the test. All bloody samples whose result is positive must be confirmed with other techniques.
- Mucous stool sample could cause non-specific reactions in the test. Any mucous stool sample with a positive result should be contrasted with other diagnostic techniques to confirm the result.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit):

The detection limit of biotical C. difficile GDH Card is 0.78 – 0.10 ng/mL of recombinant *Clostridium difficile* GDH protein.

Clinical sensitivity and specificity:

An evaluation was carried out, with stool samples. Comparing the results obtained with an immunochromatographic test (biotical C. difficile GDH Card, biotical) and an evaluation criterion. The evaluation criteria consisted of the parallel evaluation of two rapid tests (biotical C. difficile GDH Card) vs. another competitive rapid test: discrepant results were confirmed with the qPCR technique.

The results are shown below:

		qPCR		
		+	-	Total
IC test: biotical C. difficile GDH Card	+	51	0	51
	-	0	75	75
	Total	51	75	126

<i>Biotical Clostridium difficile</i> GDH Card vs evaluation criteria	Middle value	95% Confidence Interval
SENSITIVITY	100%	93,0 - 100%
SPECIFICITY	100%	95,2 - 100%
PPV	100%	93,0 - 100%
NPV	100%	95,2 - 100%
ACCURACY	100%	97,1 - 100%

The results showed that biotical C. difficile GDH Card has a high sensitivity and specificity to detect *Clostridium difficile* glutamate dehydrogenase (GDH).

CROSS REACTIVITY

An evaluation was carried out to determine possible cross-reactions of biotal Clostridium difficile GDH; no cross-reaction with organisms causing other gastrointestinal infections was detected.

<i>Adenovirus</i>	<i>Helicobacter pylori</i>	<i>Salmonella typhi</i>
<i>Astrovirus</i>	Bovine haemoglobin	<i>Salmonella typhimurium</i>
Human calprotectin	Pig haemoglobin	<i>Shigella boydii</i>
<i>Campylobacter coli</i>	Human haemoglobin	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	Bovine lactoferrin	<i>Shigella flexneri</i>
<i>Clostridium difficile Toxin A</i>	Human lactoferrin	<i>Shigella sonnei</i>
<i>Clostridium difficile Toxin B</i>	<i>Legionella pneumophila</i>	<i>Streptococcus pneumococcal</i>
<i>Clostridium perfringens</i>	<i>Listeria monocytogenes</i>	<i>Streptococcus pyogenes</i>
<i>Entamoeba histolytica</i>	<i>Norovirus GI</i>	Bovine transferrin
<i>Escherichia coli O111</i>	<i>Norovirus GII</i>	Human transferrin
<i>Escherichia coli O26</i>	<i>Rotavirus</i>	<i>Yersinia Enterocolitica O:3</i>
<i>Escherichia coli O157</i>	<i>Salmonella enteritidis</i>	<i>Yersinia Enterocolitica O:9</i>
<i>Clostridium Difficile GDH</i>	<i>Salmonella paratyphi A</i>	

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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REF

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REFERENCES / BIBLIOGRAFÍA /

1. WREN, M.W.D, et al. "Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009.
2. VAISHNAVI, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indican J. Med. Res. 131, April 2010, pp 487-499.

Sierra de Guadarrama, 1.
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es

biotical health S.L.U.