ENGLISH - ESPAÑOL

Ref: RTB25CD  
RTB25CDC  
RTB25CDCN

biotical

# C. difficile A+B Card

Rapid Test to detect  
Clostridium difficile  
A+B Toxins in card  
format.

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biotical health  
ready to test

# biotical C. difficile Toxins A+B Card

ENGLISH

## WHAT IS BIOTICAL C. DIFFICILE TOXINS A+B CARD?

---

biotical C. difficile Toxins A+B Card is a coloured chromatographic immunoassay test for the qualitative detection of *Clostridium difficile* Toxin A and B in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional can make a simple diagnosis and detect the possible presence of toxigenic strains of *Clostridium difficile* in infected people.

## WHAT IS CLOSTRIDIUM DIFFICILE (*C. difficile*)?

---

*Clostridium difficile* (*C. difficile*) is a Gram-positive anaerobic bacterium. It is the main pathogen linked to diarrhoea associated with antibiotics and / or pseudomembranous colitis in hospitalized patients, the main cause of diarrhoea associated with healthcare in developed countries, and the main origin of morbidity and mortality of nosocomial disease around the world.

Most of the bacterial flora in the colon of a healthy adult is generally resistant to colonization by *C. difficile*. If the bacterial flora is altered (due to the excessive use of antibiotics, for example), resistance to colonization by this bacterium is lost and the microorganism is allowed to multiply in the colon. In addition, its ability to form spores allows it to persist in patients for long periods and facilitates faecal-oral transmission. This colonization causes the release of:

- **Toxins A and B:** toxigenic strains responsible for the clinical manifestations of the clostridium disease, which can range from non-serious watery diarrhoea, pseudomembranous colitis, toxic megacolon, to death. biotical C. difficile Toxins A + B Card will detect the presence of these enzymes.
- **Clostridium difficile Glutamate Dehydrogenase (GDH):** enzyme produced in large quantities by toxigenic (toxins A and B) and non-toxigenic strains, which makes it an excellent marker to determine the presence of the microorganism. (See biotical C. difficile GDH Card).

Toxigenic culture (TC) is the reference technique used to determine *Clostridium difficile* infection. This method comprises the culture and isolation of *C. difficile* from stool sample, followed by identification of the toxin, which requires a very laborious work to obtain a result. However, biotical C. difficile Toxins A+B Card, is a qualitative immunochromatographic test of easy use and with almost immediate results, for the detection of toxigenic *Clostridium difficile* strains in stool samples.


## TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25CD	biotal C. difficile toxins A+B Card	25 tests without positive or negative control.
RTB25CDC	biotal C. difficile toxins A+B Card (+PC)	25 tests with positive control, without negative control.
RTB25CDCN	biotal C. difficile toxins A+B Card (+PC +NC)	25 tests with positive and negative control.

## HOW DOES BIOTAL C. DIFFICILE TOXINS A+B CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- **In the T1 RED test line (Toxin A):** with mouse monoclonal antibodies against *Clostridium difficile* TOXINS A.
- **In the C GREEN line (Control):** with rabbit polyclonal antibodies against a specific protein.
- **In the T2 RED test line (Toxin B):** with mouse monoclonal antibodies against *Clostridium difficile* TOXINS B.
- **Sample area (S):** the pad where the sample is placed has been dispensed with a preparation of reagents of the T1 line (mouse monoclonal antibodies against *Clostridium difficile* Toxin A) conjugated with latex of **RED** polystyrene; with another preparation of reagents of the T2 line (mouse monoclonal antibodies against *Clostridium difficile* Toxin B) conjugated to **RED** polystyrene latex; and with another preparation for the control solution (specific binding protein) conjugated to **GREEN** polystyrene latex, forming three conjugate colored complexes.

### TEST LINE (T1/T2):

- **If the sample is positive:** The antigens of the patient's stool sample will react with one (or both) of the red coloured conjugate complexes (anti-*Clostridium difficile* Toxin A monoclonal antibodies latex red microspheres for the line T1 and anti-*Clostridium difficile* Toxin B monoclonal antibodies - red latex microspheres for the line T2) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-*Clostridium difficile* Toxin A present in the T1 line and / or the antibodies anti-*Clostridium difficile* Toxin B present in the T2 line will capture the complex giving rise in a red line.
- **If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complexes. The anti-*Clostridium difficile* Toxin A and anti-*Clostridium difficile* Toxin B antibodies present in the T1 and T2 line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

### CONTROL LINE (C):

- **Regardless of whether the sample is positive or not,** to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

## STORAGE AND STABILITY

---

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

## GENERAL PRECAUTIONS

---

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- Used tests should be managed as sanitary waste (sanitary waste container).
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

## MATERIALS

---

### MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Tubes for the sample dilution with sample diluent.
- Positive Control (only RTB25CDC & RTB25CDCN references).
- Negative Control (only RTB25CDCN reference).

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

## SPECIMEN COLLECTION AND PREPARATION

---

Samples should be collected in a clean container. Samples must be kept cold (2-8°C), for a maximum of 1 day, until they are used. To preserve the samples for a long time, a maximum of 6 months, they must be kept frozen at -20°C. In this case, the sample must be completely thawed and brought to room temperature before it can be used in the test. Freeze and thaw cycles are not recommended, so only the amount of sample needed to perform the test should be thawed. Homogenize the sample vigorously before its preparation.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube ① and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample ②, to collect faecal sample and add it to the stool collection tube (125mg approx). For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion ③. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

TEST PROCEDURE

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
  - Remove the test card from its sealed bag just before using it.
  - Take the stool collection tube, **cut the end of the cap** ④
  - **Dispense 4 drops** from the testing tube in the circular window marked with the letter S ⑤. Avoid adding solid particles with the liquid.
  - **Read the results at 10 minutes.** The results will not be valid after 10 minutes.
- If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it does not work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)

Toxin A NEGATIVE  
Toxin B NEGATIVE

Toxin A POSITIVE  
Toxin B NEGATIVE


Toxin A NEGATIVE  
Toxin B POSITIVE

Toxin A POSITIVE  
Toxin B POSITIVE

INVALID

INVALID

INVALID

INVALID

	biological Test	Interpretation of the results
1	T1 = NO LINE C = GREEN T2 = NO LINE	<b>NEGATIVE FOR CLOSTRIDIUM DIFFICILE TOXINS A+B</b> Only the <b>GREEN</b> line appears across the results window in the site marked with the letter C (control line). No presence of CLOSTRIDIUM DIFFICILE TOXINS A+B is detected in the sample.
2	T1 = RED C = GREEN T2 = NO LINE	<b>POSITIVE FOR CLOSTRIDIUM DIFFICILE TOXIN A</b> In addition to the <b>GREEN</b> control line, a <b>RED</b> line also appears in the site marked with the letter T1 (test line) in the results window. CLOSTRIDIUM DIFFICILE TOXINS A presence is detected in the sample.
3	T1 = NO LINE C = GREEN T2 = RED	<b>POSITIVE FOR CLOSTRIDIUM DIFFICILE TOXIN B</b> In addition to the <b>GREEN</b> control line, a <b>RED</b> line also appears in the site marked with the letter T2 (test line) in the results window. CLOSTRIDIUM DIFFICILE TOXIN B presence is detected in the sample.

	biotal Test	Interpretation of the results
4	T1 = RED C = GREEN T2 = RED	<b>POSITIVE FOR CLOSTRIDIUM DIFFICILE TOXINS A + B</b> In addition to the GREEN control line, two RED lines also appears in the site marked with the letter T1 and T2 in the results window. CLOSTRIDIUM DIFFICILE TOXINS A + B presence is detected in the sample.
5-8	Any other result	<b>INVALID</b> Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the tests lines (T1 and T2) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25CDC** & **RTB25CDCN**).

## POSITIVE CONTROL: INTENDED USE & COMPOSITION

biotal C. difficile Toxin A and Toxin B Positive Control are an external quality control for the test line of Toxins A and B, respectively, of *Clostridium difficile* in biotal products. This control swab should be used to assess that the reagents and the test performance are adequate.

biotal C. difficile Toxin A Positive Control is manufactured with inactivated *Clostridium difficile* toxin A, dried on swab, contains a preservative.


biotal C. difficile Toxin B Positive Control is manufactured with inactivated *Clostridium difficile* toxin B, dried on swab, contains a preservative.

## POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by **rotating the positive control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biological device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes**. do not read the test result later than 10 minutes.

## POSITIVE CONTROL – INTERPRETATION OF RESULTS

biological Test	Interpretation of the results
 <p>T1 = <b>RED LINE</b> C = <b>GREEN LINE</b> T2 = <b>NO LINE</b></p>	<p><b>RESULT: Toxin A POSITIVE CONTROL</b> A <b>GREEN</b> line (C) and a <b>RED</b> line (T1) appears.</p>
 <p>T1 = <b>NO LINE</b> C = <b>GREEN LINE</b> T2 = <b>RED LINE</b></p>	<p><b>RESULT: Toxin B POSITIVE CONTROL.</b> A <b>GREEN</b> line (C) and a <b>RED</b> line (T2) appears.</p>
 <p>T1 = <b>RED LINE</b> C = <b>GREEN LINE</b> T2 = <b>RED LINE</b></p>	<p><b>RESULT: Toxin A &amp; B POSITIVE CONTROL</b> A <b>GREEN</b> line (C) and two <b>RED</b> lines (T1 &amp; T2) appears.</p>

biological Test	Interpretation of the results
	<p><b>INVALID RESULTS FOR POSITIVE CONTROL.</b></p> <p>Total absence of the control line, regardless the appearance or not of any of the test lines, or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

biological Negative Control is an external quality control for all immunochromatographic tests manufactured by biological. Use this control swab to check that the extraction reagents and the test are working properly.

biological Negative Control is manufactured with negative control buffer dried in swab containing a preservative.


## NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the negative control swab into the testing tube **1**.
- Mix the solution by **rotating the negative control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biological device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

biotical Test		Interpretation of the results
	T1 = NO LINE C = GREEN LINE T2 = NO LINE	<b>RESULT: CONTROL NEGATIVE</b> Only the control line (C) appears across the results window during the test performance.
		<b>RESULT: INVALID FOR NEGATIVE CONTROL</b> Total absence of the control line, and/or the card shows the test line.  Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new negative control swab.

## LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results
- A positive test result confirms the presence of Toxins A and B in the sample, without determining its origins. To confirm the presence of infection by *Clostridium difficile* it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques to identify the strain (toxigenic culture).
- False positives may occur in the test with bloody samples or mucous stools, so the use of the test in this class of samples is not recommended.
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. A negative result does not indicate absence of infection by *Clostridium difficile*.
- It is possible to appear false negatives by not reaching the detection limit of the test.
- Bloody stool samples contain components that could cause non-specific reactions in the test. All bloody samples whose result is positive must be confirmed with other techniques.
- Mucous stool sample could cause non-specific reactions in the test. Any mucous stool sample with a positive result should be contrasted with other diagnostic techniques to confirm the result.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity (detection limit):

- Detection Limit for toxin A: 8-1 ng/mL.
- Detection Limit for toxin B: 6.24-0.39 ng/mL.

### Clinical sensitivity and specificity (Toxin A+B):

An evaluation was carried out with stool samples, comparing the results obtained with the biotal C. difficile Toxins A+B Card immunochromatographic test against another commercial immunochromatographic test. Discordant results were confirmed by PCR techniques.

The results were as follows:

		Evaluation criteria (Toxin A)		
		+	-	Total
<b>biotal C. difficile Toxins A+B Card</b> (T1: Toxin A)	+	25	0	25
	-	2	99	101
	Total	27	99	126

<i>Biotal C. difficile Toxin A+B Card (Toxin A) vs Evaluation Criteria</i>	MIDDLE VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	92.59%	75.7% - 99.1%
SPECIFICITY	>99.9%	96.3% - 100%
VPP	>99.9%	-
VPN	98.0%	92.9% - 99.5%
ACCURACY	98.4%	94.4% - 99.8%

		Evaluation criteria (Toxin B)		
		+	-	Total
<b>biotal C. difficile Toxins A+B Card</b> (T2: Toxin B)	+	32	0	32
	-	2	92	94
	Total	34	92	126

<i>Biotal C. difficile Toxin A+B Card (Toxin B) vs Evaluation Criteria</i>	MIDDLE VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	94.1%	80.3% - 99.3%
SPECIFICITY	>99.9%	96.1% - 100%
VPP	>99.9%	-
VPN	97.9%	92.3% - 99.4%
ACCURACY	98.4%	94.4% - 99.8%

The results showed that biotal C. difficile Toxins A+B Card has a high sensitivity and specificity to detect *Clostridium difficile* Toxins A and B.

## CROSS REACTIVITY.

An evaluation was carried out to determine possible cross-reactions of biotical *Clostridium difficile* Toxins A+B; no cross-reaction with gastrointestinal pathogens, other organisms, substances and/or fecal markers occasionally present in feces was detected.

<i>Adenovirus</i>	<i>Escherichia coli</i> O26	<i>Norovirus</i> GI	<i>Shigella sonnei</i>
<i>Astrovirus</i>	<i>Escherichia coli</i> O157	<i>Norovirus</i> GII	<i>Streptococcus pneumoniae</i>
Human Calprotectin	<i>Giardia</i>	<i>Rotavirus</i>	<i>Streptococcus pyogenes</i>
<i>Campylobacter coli</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>	Human transferrin
<i>Campylobacter jejuni</i>	Bovine hemoglobin	<i>Salmonella paratyphi</i> A	<i>Yersinia Enterocolitica</i> O:3
<i>Clostridium difficile</i> GDH	Pork hemoglobin	<i>Salmonella typhi</i>	<i>Yersinia Enterocolitica</i> O:9
<i>Clostridium perfringens</i>	Human hemoglobin	<i>Salmonella typhimurium</i>	<i>C.difficile</i> Tox.A (Tox.B T2 line)
<i>Entamoeba dispar</i>	Human lactoferrin	<i>Shigella boydii</i>	<i>C.difficile</i> Tox.B (Tox.A T1 line)
<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Shigella dysenteriae</i>	
<i>Escherichia coli</i> O111	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>	

## SYMBOLS FOR IVD COMPONENTS & REAGENTS

In vitro diagnostic device

Consult instructions for use

Keep dry

Use by

Temperature limitation

Manufacturer

Contains sufficient for <n> test

DIL

Sample diluent

LOT

Batch code

REF

Catalogue number

\*\*\*

# biotical C. difficile Toxins A+B Card

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---

1. WREN, M.W.D, et al. "Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009.
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indican J. Med. Res. 131, April 2010, pp 487-499.

Sierra de Guadarrama, 1. Ed. Dublín  
28830 San Fernando de Henares. Madrid  
Tel.: + 34 91 677 43 08

[info@biotical.es](mailto:info@biotical.es)  
[www.biotical.es](http://www.biotical.es)