ENGLISH - ESPAÑOL

Ref: RTB25G  
RTB25GC  
RTB25GCN

# biotical Giardia Card

Rapid Test to detect  
Giardia in card format.

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biotical health  
ready to test

IU RTB25G-RTB25GC-RTB25GCN Rev. 04 / 2024-04

ISO 13485 CERTIFIED  
MEDICAL DEVICE  
QUALITY MANAGEMENT SYSTEM

# biotical Giardia Card

ENGLISH

## WHAT IS BIOTICAL GIARDIA CARD?

biotical Giardia Card is a coloured chromatographic immunoassay test for the qualitative detection of Giardia antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of giardiasis.

### GIARDIASIS

Giardia is one of the most frequent parasites that cause water-borne diseases in humans: cryptosporidiosis and giardiasis. Giardiasis is an intestinal disease caused by *Giardia intestinalis* (also called *Giardia lamblia* and *Giardia duodenalis*), a single-celled microscopic parasite that is found in the intestine of people and animals and is transmitted through the faeces of an infected person or animal.

This parasite is protected by an outer covering (cyst) that allows them to survive outside the body and in the environment for long periods of time. Its infective form is usually by faecal-oral transmission. Among the most common symptoms of both diseases are watery diarrhoea, cramps and stomach disorders. In immunocompromised people or in elderly people, complications can appear that lead to serious illness and even death.

Generally, the prophylactic measure to remove *Giardia* cysts from drinking water is disinfection by chlorination or by using sand filters, where the cysts would be retained.

The reference method for the diagnosis of giardiasis diseases is the identification of cysts in an examination with optical microscopy. With this complex technique, non-immediate results are obtained. Considering that giardiasis has treatment, it is important to do an immediate diagnosis.

biotical Giardia card, is a qualitative immunochromatographic test for the detection of Giardia antigens in stool samples. The result is obtained in 10 minutes.


## TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25G	biotical Giardia Card	25 tests without positive or negative controls.
RTB25GC	biotical Giardia Card (+PC)	25 tests with positive control, without negative control.
RTB25GCN	biotical Giardia Card (+PC +NC)	25 tests with positive and negative control.

## HOW DOES BIOTICAL GIARDIA CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against Giardia.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against Giardia) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex.

### TEST LINE (T):

**If the sample is positive:** The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-Giardia monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-Giardia present in the T line will capture the complex giving rise in a red line.

**If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-Giardia antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

### CONTROL LINE (C):

**Regardless of whether the sample is positive or not**, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

## STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

## GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in this test. Do not use any other commercial kit component.

- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

## MATERIALS

### MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25GC & RTB25GCN references).
- Negative Control (only RTB25GCN reference).

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.


The kit contains materials for the analysis of a maximum of 25 samples or controls.

## SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

### Specimen preparation (see illustration):

- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. Use only fresh samples or fresh frozen samples without preservatives or fixatives.


Collection tube

**1** Stool sample

**2** Shake vigorously

## TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

3 Cut the end of the cap.


4 Add 3 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 3 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID RESULT

INVALID RESULT

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	<b>NEGATIVE</b> Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of Giardia is detected in the sample.
2	C = GREEN LINE T = RED LINE	<b>POSITIVE</b> In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. Giardia presence is detected in the sample.
3-4	Any other result	<b>INVALID</b> Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25GC** & **RTB25GCN**).


## POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical Giardia Positive Control is an external quality control for biotical Giardia card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical Giardia Positive Control: inactivated *Giardia* antigens dried in swab containing a preservative.

## POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

## POSITIVE CONTROL - INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p><b>RESULT: CONTROL POSITIVE</b></p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p><b>INVALID RESULT</b></p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

## NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the negative control swab into the testing tube **1**.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the Biotical test Card device from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube in the corresponding circular window **3**. Use a separate pipette and test for each control.
- **Read the result at 10 minutes**. Do not read the test result later than 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = <b>GREEN LINE</b> T = <b>NO LINE</b></p>	<p><b>RESULT: CONTROL NEGATIVE</b></p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p><b>RESULT: INVALID FOR NEGATIVE CONTROL</b></p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

## LIMITATIONS

- The test must be carried out within 2 hours after opening the sealed bag.
- Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- After one week of infection, the number of parasites in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
- This test should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of Giardia ( $\alpha$ 1-giardin and/or CWP1) in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods or microscopy) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation still persist, a Giardia determination should be carried out with another technique (for example microscopy).
- Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
- Mucous stool samples may cause non-specific reactions in the test- Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

## PERFORMANCE CHARACTERISTICS

### Clinical sensitivity and specificity.

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical Giardia Card, biotical) vs the Evaluation Criteria.

Evaluation criteria: three rapid test were evaluated (biotical Giardia Card, Certest and other two rapid test from the competitors); the discrepant results were confirmed by qPCR technique (Viasure Cryptosporidium, Giardia & E. histolytica Real Time Detection Kits (CerTest)).

The results were as follows:

		Viasure Cryptosporidium, Giardia & E. histolytica Real Time Detection Kits (qPCR)		
		+	-	Total
Biotical Giardia Card	+	44	2	46
	-	1	79	80
	Total	45	81	126

<i>Biotical Giardia Card vs Viasure Cryptosporidium, Giardia &amp; E. histolytica Real Time Detection Kits (Giardia)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	97.8%	88.2 – 99.9%
SPECIFICITY	97.5%	91.4 – 99.7%
PPV	95.7%	85.2 – 99.5%
NPV	98.8%	93.2 – 100.0%
ACCURACY	97.8%	93.4 – 99.6%

The results showed a high sensitivity and specificity to detect Giardia using biotical Giardia Card.

## CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Giardia; no cross reactivity against gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

<i>Adenovirus</i>	<i>Entamoeba histolytica</i>	<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>
<i>Astrovirus</i>	<i>Escherichia coli O:157</i>	<i>Norovirus GI</i>	<i>Shigella dysenteriae</i>
<i>Calprotectin</i>	<i>Helicobacter pylori</i>	<i>Norovirus GII</i>	<i>Shigella flexneri</i>
<i>Campylobacter jejuni</i>	Bovine hemoglobin	<i>Resp. Syncytial Virus</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile antigen GDH</i>	Porcine hemoglobin	<i>Rotavirus</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile Toxin A</i>	Influenza A	<i>Salmonella enteritidis</i>	<i>Streptococcus pneumococcal</i>
<i>Clostridium difficile Toxin B</i>	Influenza B	<i>Salmonella paratyphi A</i>	Human transferrin
SARS-Cov-2	Human lactoferrin	<i>Salmonella typhi</i>	<i>Yersinia Enterocolitica O:3</i>
<i>Cryptosporidium parvum</i>	<i>Legionella pneumophila</i>	<i>Salmonella typhimurium</i>	<i>Yersinia Enterocolitica O:9</i>

## SYMBOLS FOR IVD COMPONENTS & REAGENTS

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