ENGLISH - ESPAÑOL

Ref: RTB25K
RTB25KC
RTB25KCN

biotical Crypto Card

Rapid Test to detect
Cryptosporidium in card
format

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biotical health
ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical Crypto Card

ENGLISH

WHAT IS BIOTICAL CRYPTO CARD?

biotical Crypto Card is a coloured chromatographic immunoassay test for the qualitative detection of *Cryptosporidium* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of cryptosporidiosis.

CRYPTOSPOROSIS

Giardia and *Cryptosporidium parvum* are the most frequent parasites that cause water-borne diseases in humans: Cryptosporidiosis and Giardiasis.

Cryptosporidium parvum is one of the most frequent causes of diarrheal processes in emerging countries. This parasite is in the majority of enteric pathogenic infections, being transmitted faecal-orally. Some of the most common symptoms of cryptosporidiosis are watery diarrhoea, severe stomach pains, weight loss, nausea, and sometimes fever.

In general, the prophylactic measure to remove *Giardia* cysts from drinking water is chlorine disinfection or using sand filters, where the cysts would be retained. This does not happen with *Cryptosporidium* cysts, because they are resistant to chlorine and their small size makes it difficult to purify using filters. The curative pharmacological treatment useful for giardiasis is not effective for cryptosporidiosis, which currently has no treatment.

biotical Crypto Card, is a qualitative immunochromatographic test for the detection of *Cryptosporidium* antigens in stool samples. The result is obtained in 10 minutes.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25K	biotical Crypto Card	25 test without positive or negative controls.
RTB25KC	biotical Crypto Card (+PC)	25 test with positive control, without negative control.
RTB25KCN	biotical Crypto Card (+PC +NC)	25 test with positive and negative control.

HOW DOES BIOTICAL CRYPTO CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In the **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In the **RED** line (T): with mouse monoclonal antibodies against *Cryptosporidium*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against *Cryptosporidium*) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex.

TEST LINES:

If the sample is positive: The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-*Cryptosporidium* monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-*Cryptosporidium* present in the T line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-*Cryptosporidium* antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25KC & RTB25KCN references).
- Negative Control (only RTB25KCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

NOTE: If you have purchased the kit without positive or negative controls, it is possible to purchase it independently.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.

Specimen preparation (see illustration):

- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

Collection tube


1 Stool sample

2 Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

3 Cut the end of the cap.


4 Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>Cryptosporidium</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>Cryptosporidium</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references ended with “C” or “CN”).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical Crypto Positive Control is an external quality control for biotical Crypto Card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical Crypto Positive Control: inactivated *Cryptosporidium parvum* antigens dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by **rotating the positive control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10-15 minutes.**

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

NEGATIVE CONTROL: INTENDED USE & COMPOSITION

biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- 1 Add 15 drops of diluent and introduce the negative control swab into the testing tube.
- 2 Rotating 30 seconds and extract the liquid.
- 3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube 1.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10-15 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test must be performed within 2 hours after opening the sealed pouch.
- Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
- biotical Crypto Card will only indicate the presence of *Cryptosporidium* in the specimen (qualitative detection) and should be used for the detection of *Cryptosporidium* antigens in human faeces specimens only. Neither the quantitative value nor the rate of increase in *Cryptosporidium* antigens concentration can be determined by this test.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- Samples should be collected as soon as possible once the symptoms appear. After one week of infection the number of parasites decrease, making the sample less reactive.
- Some stool samples can decrease the intensity of the control line.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods are recommended. A negative result does not at any time preclude the possibility of *Cryptosporidium* and/or *Giardia* infections it may be due to a very low antigen concentration.
- This test provides a presumptive diagnosis of *Cryptosporidium* infections. All results must be interpreted together with other clinical information and laboratory findings available to the physician.

EXPECTED VALUES

It is estimated that diseases caused by the presence of parasites cause 33% of global deaths, with intestinal parasitic infections being the most prevalent. Mainly caused by the continued consumption of unsafe water and unsafe food, intestinal parasites are responsible for more than 800 million episodes of diarrhoea and 4.5 million deaths directly or indirectly.

The enteric protozoa *Cryptosporidium parvum* is generally recognized as important causal agents of morbidity and mortality associated with episodes or diarrheal diseases in developed countries.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity.

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical Crypto Card, biotical) vs the Evaluation Criteria [Microscopy technics and qPCR confirmation method (Viasure Crypto qPCR)].

The results were as follows:

		Evaluation Criteria (Microscopy + qPCR)		
		+	-	Total
Biotical Crypto Card	+	25	0	25
	-	0	229	229
	Total	25	229	254

<i>Biotical Crypto Card vs Evaluation Criteria (Crypto)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99%	86,28% - 100%
SPECIFICITY	>99%	98,40% - 100%
PPV	>99%	-
NPV	>99%	-
ACCURACY	>99%	98,56% - 100%

The results showed a high sensitivity and specificity to detect *Cryptosporidium* using biotical Crypto Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Crypto Card, no cross reactivity against gastrointestinal organism occasionally present in faeces.

Campylobacter jejuni

Campylobacter coli

Clostridium difficile

Escherichia coli O157:H7

Entamoeba histolytica

Giardia lamblia

Helicobacter pylori

Listeria monocytogenes

Salmonella enteritidis

Salmonella paratyphi

Salmonella typhi

Salmonella typhimurium

Shigella boydii

Shigella dysenteriae

Shigella flexneri

Shigella sonnei

Staphylococcus aureus

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

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REF

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REFERENCES / BIBLIOGRAFÍA /

1. PENG M. M., XIAO L., FREEMAN A. R., "Genetic Polymorphism Among *Cryptosporidium parvum* Isolates. Evidence of Two Distinct Human Transmission Cycles", *Emerging Infectious Diseases*, Vol 3 No 4, Oct-Dec. 1997, pp 567-573.
2. MARSHALL, M.M., et al., "Waterborne Protozoan Pathogens", *Clinical Microbiology Review*, Jan. 1997, pp 67-85.
3. GUERRANT, R.L., "Cryptosporidiosis: An Emerging, Highly Infectious Threat", *Emerging Infectious Diseases*, Vol 3 No 1, January-March 1997, pp 51-57.
4. FAYER, R. and UNGAR, L.P., "*Cryptosporidium* spp. and Cryptosporidiosis", *Microbiological Reviews*, Dec. 1986, Vol. 50, No. 4, pp. 458- 483.

Sierra de Guadarrama, 1. Ed. Dublín.
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es

biotical health S.L.U.