ENGLISH - ESPAÑOL

Ref: RTB25EH
RTB25EHC
RTB25EHCN

biotical Entamoeba Card

Rapid Test to detect
Entamoeba in card
format

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biotical health
ready to test

IU RTB25EH-RTB25EHC-RTB25EHCN Rev. 02 / 2024-01

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical Entamoeba Card

ENGLISH

WHAT IS BIOTICAL ENTAMOEBEA CARD?

biotical Entamoeba Card is a rapid immunochromatography test for the qualitative detection of *Entamoeba histolytica* and *Entamoeba dispar* antigens in stool samples. Thanks to its easy handling and high sensitivity, the healthcare professional can make a presumptive diagnosis of amebiasis.

AMEBIASIS

Amebiasis is an infection of the gastrointestinal tract due to *Entamoeba histolytica*, a protozoan capable of invading the intestinal mucosa and can spread to other organs. *Entamoeba dispar* (amoeba morphologically similar to *E. histolytica*), is the species without invasive potential and is the cause of most asymptomatic infections.

Amebiasis uses to occur in areas where sanitary conditions are inadequate. On the entire planet, each year about 40 to 50 million people develop amebiasis, and between 40,000 and 70,000 of these people die as a result of the disease, making amebiasis the fourth leading cause of death due to protozoan infection.

Entamoeba exists in two forms: trophozoite (active parasite) or cyst (latent parasite). Many people have no symptoms or have few symptoms; however, they excrete cysts in stool and, therefore, they can spread the infection.

The disease can manifest itself as acute, chronic or as an asymptomatic infection.

Possible symptoms of amebiasis include flatulence, painful abdominal cramps, intermittent diarrhoea or constipation and, in most severe cases, dysentery (loose stools with mucus and blood), formation of amebomas (large lumps in the colon) and abscesses in the liver. The infection occurs in contaminated water, fresh vegetables and fruits or other undercooked or washed food. It is also possible that faeces cysts are transported to food through flies or cockroaches.

The diagnosis is usually analysing the stool and, if necessary, performing a colonoscopy, ultrasound or blood test. biotical Entamoeba card is a qualitative immunochromatographic test for the detection of and Entamoeba histolytic and dispar antigens in stool samples. The result is obtained in 10 minutes.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25EH	biotical Entamoeba Card	25 tests without positive or negative controls.
RTB25EHC	biotical Entamoeba Card (+PC)	25 tests with positive control, without negative control.
RTB25EHCN	biotical Entamoeba Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL ENTAMOEBEA CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In the **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In the **RED** line (T): with mouse monoclonal antibodies against *Entamoeba*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against *Entamoeba*) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex.

TEST LINES:

If the sample is positive: The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-*Entamoeba* monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-*Entamoeba* present in the T line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-*Entamoeba* antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.

- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25EHC & RTB25EHCN references).
- Negative Control (only RTB25EHCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.

The kit contains materials for the analysis of a maximum of 25 samples or controls.


SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.


Specimen preparation (see illustration):

- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample (aprox. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.


- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 3.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 4. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>Entamoeba histolytica</i> or <i>Entamoeba dispar</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>Entamoeba histolytica</i> or <i>dispar</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25EHC** & **RTB25EHCN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical Entamoeba Positive Control is an external quality control for biotical Entamoeba card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical Entamoeba Positive Control: inactivated *Entamoeba* antigens dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL: INTERPRETATION OF RESULTS

	biotal Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the negative control swab into the testing tube **1**.
- Mix the solution by **rotating the negative control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes**. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test must be carried out within 2 hours after opening the sealed bag.
- Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- This test should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- After one week of infection, the number of parasites in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
- Positive results determine the presence of *Entamoeba histolytica* and/or *E. dispar* in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods or microscopy). A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, an *Entamoeba histolytica* and/or *E. dispar* determination should be carried out with another technique (for example, microscopy).
- Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
- Mucous stool samples may cause non-specific reactions in the test- Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation was performed, using faecal samples, comparing the results obtained by an immunochromatographic test (biotical Entamoeba, biotical) and other commercial qPCR tests (VIASURE Entamoeba histolytica Real Time Detection Kit and VIASURE Entamoeba dispar Real Time Detection Kit, Certest).

The results were as follows:

		qPCR test: VIASURE Entamoeba histolytica Real Time detection Kit		
		+	-	Total
IC test: Biotical Entamoeba Card	+	4	1	5
	-	1	107	108
	Total	5	108	113

<i>Biotical Entamoeba Card vs VIASURE Entamoeba Histolytica Real Time detection Kit</i>	MEAN VALUE E. histolytica	95% CONFIDENCE INTERVAL
SENSITIVITY	80.0%	28.4 - 99.5%
SPECIFICITY	99.1%	94.9 - 100%
PPV	80.0%	35.1 - 96.7%
NPV	99.1%	94.9 - 99.8%
ACCURACY	98.2%	93.8 - 99.8%

		qPCR test: VIASURE Entamoeba Dispar Real Time detection Kit		
		+	-	Total
IC test: Biotical Entamoeba Card	+	15	1	19
	-	6	107	113
	Total	21	108	129

<i>Biotical Entamoeba Card vs VIASURE Entamoeba Dispar Real Time detection Kit</i>	MEAN VALUE E. dispar	95% CONFIDENCE INTERVAL
SENSITIVITY	71.4%	47.8 - 88.7%
SPECIFICITY	99.1%	95.0 - 100%
PPV	93.8%	67.7 - 99.1%
NPV	94.8%	90.2 - 97.3%
ACCURACY	94.6%	89.2 - 97.8%

		qPCR test: VIASURE Entamoeba Histolytica Real Time Detection Kit + VIASURE Entamoeba Dispar Real Time detection Kit		
		+	-	Total
IC test: Biotical Entamoeba Card	+	19	2	21
	-	7	214	221
	Total	26	216	242

<i>Biotical Entamoeba Card vs VIASURE Entamoeba Histolytica Real Time detection Kit + VIASURE Entamoeba Histolytica Real Time detection Kit</i>	MEAN VALUE E. histolytica + E. dispar	95% CONFIDENCE INTERVAL
SENSITIVITY	73.1%	52.2 – 88.4%
SPECIFICITY	99.1%	96.7 – 99.9%
PPV	90.5%	70.1 – 97.5%
NPV	96.8%	94.2 – 98.3%
ACCURACY	96.3%	93.1 – 98.3%

The results showed a high sensitivity and specificity to detect *Entamoeba histolytica* & *Entamoeba dispar* using biotical Entamoeba Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Entamoeba Card; no cross reactivity against gastrointestinal organism occasionally present in faeces.

Campylobacter jejuni

Campylobacter coli

Clostridium difficile

Escherichia coli O157:H7

Cryptosporidium parvum

Giardia lamblia

Helicobacter pylori

Listeria monocytogenes

Salmonella enteritidis

Salmonella paratyphi

Salmonella typhi

Salmonella typhimurium

Shigella boydii

Shigella dysenteriae

Shigella flexneri

Shigella sonnei

Staphylococcus aureus

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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REF

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REFERENCES / BIBLIOGRAFÍA /

BLESSMANN J., et al. "Longitudinal Study of Intestinal Entamoeba histolytica Infections in Asymptomatic Adult Carriers". Journal of clinical microbiology, Oct. 2003, p. 4745–4750 Vol. 41, No. 10.

GONIN P., et al. "Detection and Differentiation of Entamoeba histolytica and Entamoeba dispar Isolates in Clinical Samples by PCR and Enzyme-Linked Immunosorbent Assay". Journal of clinical microbiology, Jan. 2003, p. 237–241 Vol. 41, No. 1

Sierra de Guadarrama, 1. Ed. Dublín
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es

biotical health S.L.U.