ENGLISH - ESPAÑOL

Ref: RTB25SC
RTB25SCN

biotical Strep A Card

Rapid Test to detect
Strep A in card format

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biotical health
ready to test

IU RTB25SC-RTB25SCN Rev. 03 / 2024-03

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical STREP A Card

ENGLISH

WHAT IS BIOTICAL STREP A CARD?

biotical Strep A Card is a rapid immunochromatography test for the qualitative detection of Group A *Streptococcus* antigens (GAS) in throat swab samples and culture of isolated colonies suspected of Group A *Streptococcus*. Thanks to its easy handling and high sensitivity, the healthcare professional is able to make a presumptive diagnosis of respiratory infection caused by GAS.

INFECTIONS BY STREPTOCOCCUS GROUP A (GAS)

Group A *Streptococcus* (Strep) is a bacterium that is commonly found in the throat and causes diseases such as tonsillitis and pharyngitis. Symptoms of this infection include fever, stomach pain and inflammation and redness in the tonsils, and may be moderate or severe. The treatment in case of infection by Group A *Streptococcus* is the taking of antibiotics, so an early detection helps to treat the infection before and get a faster cure. Similarly, confusion with infections of viral origin is avoided for which the taking of antibiotics is not prescribed facilitating the rational use of them.

biotical Strep A card is an immunochromatographic qualitative test, easy to use and with almost immediate results, for the detection of Group A *Streptococcus* directly in nasopharyngeal swab samples, providing results in several minutes.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25SC	biotical Strep A Card (+PC)	25 tests with positive control, without negative control.
RTB25SCN	biotical Strep A Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL STREP A CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In the **BLUE** line (C): with rabbit polyclonal antibodies against a specific protein.
- In the **RED** line (T): with mouse monoclonal antibodies against Group A *Streptococcus*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a test solution (mouse monoclonal antibodies against Group A *Streptococcus*) conjugated to red polystyrene latex and a control solution (specific binding protein) conjugated to blue polystyrene latex, forming two conjugate coloured complexes.

TEST LINE (T):

If the sample is positive: The antigens of the patient's sample will react with the red coloured conjugate complex (anti-Group A *Streptococcus* monoclonal antibodies - red latex microspheres) and advance by capillarity through the membrane. The anti-Group A *Streptococcus* antibodies present in the membrane of the test line (T) will capture the coloured test complex resulting in the red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-Group A *Streptococcus* antibodies present in the membrane of the test line (T) will not capture the red coloured antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the blue control line will always appear.

STORAGE AND STABILITY

The product should be stored at 2-30°C in its original sealed container for optimum performance until the expiry date printed on the container. Do not open until ready for use. Do not freeze. Once opened, reagents A and B should be stored according to the conditions indicated on the label. Under these storage conditions, the reagents retain their activity throughout their shelf life (as indicated on the respective label). All kit components should be conditioned at room temperature (15-30°C) before testing.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Read the instructions for use carefully before using the kit.
- Do not use after expiration date.
- Do not use the kit if the sticker covering the outer box is damaged or broken.

- Do not use the product if the box or foil envelope is open or damaged on receipt.
- Do not use Reagent A and Reagent B if they are open or damaged upon receipt.
- Do not change the lids between Reagent A and Reagent B.
- Do not use the test (or positive / negative control if present in the final reference) if the desiccant included in each foil pouch is missing or damaged.
- All kit components, except Reagent A and Reagent B, are for single use only. Do not reuse.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing. These containers should be disposed of in accordance with local or national laws and regulations.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing (PPE), use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area. Wash hands after testing is completed.
- Clean up spills with a suitable disinfectant.
- The reagents contain preservatives (12-14% sodium nitrite). Contact with skin or mucous membranes should be avoided. According to Regulation (EC) N° 1907/2006 (REACH), biotical Strep A Card contains substances and/or mixtures that meet the hazard classification criteria set out in Regulation (EC) N° 1272/2008 (CLP) or in concentrations exceeding the values specified in the applicable reporting regulations. Safety data sheets are not included with this product but are available on request from biotical health.
- A certificate of analysis is not included with this product but can be requested from biotical health if required.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Reagent A (2M Sodium Nitrite).
- Reagent B (0.15M Acetic Acid).
- Sampling swabs.
- Disposable pipettes.
- Testing tubes
- Positive Control.
- Negative Control (only RTB25SCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Blood-agar.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Throat swab method:

- Remove the swab from its packing.
- Collect specimen with a sterile swab from the tonsils and/or the back of the throat. Take care to avoid the teeth, gums, tongue or cheek surfaces.
- Insert the swab into the tonsils and/or the back of the throat.
- Process the swab as soon as possible after collecting the specimen.

Suspected Group A *Streptococcus* colonies:

The Blood Agar is the method of choice for the isolation of Group A *Streptococcus* (anaerobic atmosphere, 72 hours /37°C). After 72 hours of incubation in blood the typical Group A *Streptococcus* colonies will be translucent white showing beta haemolysis

- Examine Blood plates after incubation period. Select *Streptococcus* typical colonies. Use a sterile swab to pick up 1 or 3 colonies suspected (showing characteristic beta haemolysis).
- Follow the instructions in the test procedure section to test the swab.

Immediately continue with the test procedure.


Each sample vial (containing the sample) and cassette used for the analysis of the sample must be perfectly recognisable, including the identity of the patient.

TEST PROCEDURE


Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- 1 Add 5 drops of Reagent A.
- 2 Add 5 drops of Reagent B.


- 3 Mix the solution with the swab at least 30 seconds.

- Add 5 drops of Reagent A (light pink) ① and 5 drops Reagent B (pale yellow) ② in a testing tube (the solution should turn light yellow (colourless)) immediately put the swab into the tube.
- Mix the solution by rotating the swab forcefully against the side of the tube at least 30 seconds ③. Best results are obtained when the specimen is vigorously extracted in the solution. Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- Dispense 3 drops from the testing tube, into the circular window marked with the letter S ④. Use a separate pipette and test for each sample or control.
- Read the results at 10 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotal Test	Interpretation of the results
1	C = BLUE LINE T = NO LINE	NEGATIVE FOR STREP A Only the BLUE line appears across the results window in the site marked with the letter C (control line). No presence of <i>Group A Streptococcus</i> (Strep A) is detected in the sample.
2	C = BLUE LINE T = RED LINE	POSITIVE FOR STREP A In addition to the BLUE control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>Group A Streptococcus</i> (Strep A) presence is detected in the sample.

		INVALID
3-4	Any other result	Total absence of the control line (BLUE) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A blue line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections ([POSITIVE CONTROL](#) & [NEGATIVE CONTROL](#)) only apply if the corresponding positive and negative controls are available (sold separately or with kits references [RTB25SC](#) & [RTB25SCN](#)).

POSITIVE CONTROL: INTENDED USE & COMPOSITION

The biotal STREP A Positive Control is an external quality control for biotal STREP A card. Use this control swab to check that the sample reagents and the test card are working properly.

biotal STREP A Positive Control is manufactured with inactivated *Group A Streptococcus* antigens extract dried in swab containing a preservative.

POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

Results window

BLUE C line (Control)
RED T line (Test Strep A)

Sample window (S)

- 1 Add 5 drops of Reagent A.
- 2 Add 5 drops of Reagent B.
- 3 Mix the solution with the swab at least 30 seconds.
- 4 Add 4 drops in the circular window (S).

- Add 5 drops of Reagent A (light pink) ① and 5 drops Reagent B (pale yellow) ② in a testing tube (the solution should turn light yellow (colourless)) immediately put the positive control swab into the tube.
- Mix the solution by rotating the swab forcefully against the side of the tube at least 30 seconds ③. Best results are obtained when the specimen is vigorously extracted in the solution. Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- Dispense 4 drops from the testing tube, into the circular window marked with the letter S ④. Use a separate pipette and test for each sample or control.
- Read the results at 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotal Test	Interpretation of the results
1	 <p>C = BLUE LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A BLUE line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (BLUE), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add **5 drops of Reagent A (light pink)** ① and **5 drops Reagent B (pale yellow)** ② in a testing tube (the solution should turn light yellow (colourless)) immediately put the negative control swab into the tube.
- Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds** ③. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ④. Use a separate pipette and test for each sample or control.
- **Read the results at 10 minutes.**

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = BLUE LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test must be carried out within 2 hours after opening the sealed bag.
- This test does not differentiate between carriers and acute infection. Pharyngitis may be caused by organisms other than Group A *Streptococcus*.
- The test should be used only with throat swabs or colonies taken directly from a plate. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper throat swab specimens must be obtained.
- The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- Positive results determine the presence of Group A *Streptococcus* infection. However, positive results require comparison with other diagnostic methods to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the throat samples is lower than the detection limit value. If the symptoms or situation still persist, a Group A *Streptococcus* determination should be carried out, on a sample from an enrichment culture.
- Excess blood or mucus in the swab sample may interfere with the performance of the test and give a false positive result. Avoid touching the tongue, cheeks, teeth and any bleeding areas of the mouth with the swab when collecting samples.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity.

An evaluation, with throat specimens from patients with pharyngitis symptoms, was performed comparing the results obtained by an immunochromatographic test (biotical Strep A Card, biotical) and another rapid test commercially available (OSOM® Strep A Test, Genzyme Diagnostics).

The results were as follows:

		OSOM® Strep A Test (immunochromatographic test)		
		+	-	Total
Biotical Strep A Card	+	4	0	4
	-	0	20	20
	Total	4	20	24

	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99.9%	39.8% - 100%
SPECIFICITY	>99.9%	83.2% - 100%
PPV	>99.9%	-
NPV	>99.9%	-
ACCURACY	>99.9%	85.8% - 100 %

The results showed a high sensitivity and specificity to detect Group A *Streptococcus* using biotical Strep A Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Strep A; no cross reactivity against organisms that cause other respiratory infections.

Influenza type A

Influenza type B

Adenovirus

Streptococcus Grupo D: Enterococcus

SYMBOLS FOR IVD COMPONENTS & REAGENTS

In vitro diagnostic device

Consult instructions for use

Keep dry

Use by

Temperature limitation

Manufacturer

Contains sufficient for <n> test

DIL

Sample diluent

LOT

Batch code

REF

Catalogue number

HAZARD IDENTIFICATION.

REAGENT A: ATTENTION

Danger indication:

- Harmful if swallowed.

In case of ingestion, contact the Toxicological Information Center and go to your emergency doctor.

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