ENGLISH - ESPAÑOL

Ref: RTB25RVC
RTB25RVCN

biotical RSV Card

Rapid Test to detect
RSV in card format

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biotical health
ready to test

IU RTB25RVC-RTB25RVCN Rev. 03 / 2024-02

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical RSV Card

ENGLISH

WHAT IS BIOTICAL RSV CARD?

biotical RSV card is a rapid immunochromatography test for the qualitative detection of Respiratory Syncytial Virus (RSV) antigens in nasal swab sample, or through nasopharyngeal washing or nasopharyngeal aspirate. Thanks to its easy handling and high sensitivity, the healthcare professional is able to make a presumptive diagnosis of infections caused by RSV.

INFECTIONS BY RESPIRATORY SYNCYTIAL VIRUS (RSV)

The Respiratory Syncytial Virus (RSV) is one of the main respiratory diseases in children younger than 1 year, affecting the lungs, respiratory tracks and also presenting fever and cough. In the case of most healthy children, RSV causes like a cold that usually lasts between 5 and 7 days, but in some cases, it can be followed by bronchiolitis or pneumonia. They can be contracted by adults and can be serious in the elderly if it affects the lower respiratory tract.

Because it is a virus, antibiotics are not useful to fight it, so it is very important for the healthcare professional to make a correct identification.

biotical RSV card is a qualitative immunochromatographic test, easy to use and with almost immediate results, for the detection of RSV directly in nasopharyngeal swab samples, providing results in 10 minutes.

TABLE OF REFERENCES


In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25RVC	biotical RSV Card (+PC)	25 test with positive control, without negative control.
RTB25RVCN	biotical RSV Card (+PC +NC)	25 test with positive and negative control.

HOW DOES BIOTICAL RSV CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against RSV.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a test solution (mouse monoclonal antibodies against RSV conjugated to red polystyrene latex and a control solution (specific binding protein) conjugated to green polystyrene latex, forming two conjugate coloured complexes.

TEST LINE:

If the sample is positive: the antigens of the patient's sample will react with the red coloured conjugate complex (anti- RSV monoclonal antibodies - red latex microspheres) and advance by capillarity through the membrane. The anti- RSV antibodies present in the T line will capture the coloured test complex resulting in the red line.

If the sample is negative: if the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti- RSV antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Sample diluent (Reagent B).
- Swabs for sample taking.
- Disposable pipettes.
- Testing tubes
- Positive Control.
- Negative Control (only RTB25RVCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Vortex.

The kit contains materials for the analysis of a maximum of 25 samples or controls.


SPECIMEN COLLECTION AND PREPARATION.

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Nasal swab method:


- Remove the swab from its packing.
- Collect specimen with a sterile swab from one nostril.
- Insert the swab into the nostril to the nasopharynx, rotating against the nasal wall (to ensure swab contains cells as well as mucus).
- Repeat procedure using another nostril.
- Process the swab as soon as possible after collecting the specimen.

Nasopharyngeal aspirate method (suction device, sterile suction catheter):

For Adults:

- Place the irrigator up to the nose.
- Let the sterile saline water run into the nose (2.5 mL). It will run out the opposite side.
- Tilt and twist the irrigator side to side and up and down directing the water flow into all portions of the nasal cavity.
- Collect the wash in a clean specimen container, tilt the head forward and allow the water with mucus to run out of the nostril into the specimen container. Repeat the mucus collection for the other nostril and collect it into the same container.

For children:

- Use an aspiration bulb or bulb syringe to instill the saline water into one nostril leaning the children head.
- Aspirate the mix of mucus-saline water into the bulb and transfer it into a clean container.
- Repeat for the other nostril and transfer the fluid into the same specimen container.

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

Procedure A using nasal swab samples:

Results window

GREEN C Line (Control)


RED T line (Test RSV)

Sample window (S)

- 1 Add 15 drops of sample diluent
- 2 Put the swab into the tube, rotating 30 seconds and extract the liquid.
- 3 Add 4 drops in the circular window (S).

- Add **15 drops ①** of **sample diluent** and immediately put the swab into the tube.
- Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds ②**. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

Procedure B using nasopharyngeal wash or aspirate samples:


- Add **6 drops of the nasopharyngeal wash or aspirate ①** samples with a pipette and **9 drops of sample diluent ②** in a testing tube. **Mix with vortex at least 30 seconds ③** to homogenize. Best results are obtained when the specimen is vigorously extracted in the solution.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S. Use a separate pipette and test for each sample or control ④.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent B until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	<p style="text-align: center;">NEGATIVE FOR RSV</p> <p>Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of RSV is detected in the sample.</p>
2	C = GREEN LINE T = RED LINE	<p style="text-align: center;">POSITIVE FOR RSV</p> <p>In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. RSV presence is detected in the sample.</p>
3-4	Any other result	<p style="text-align: center;">INVALID</p> <p>Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.</p>

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25RVC** & **RTB25RVCN**).


POSITIVE CONTROL: INTENDED USE & COMPOSITION

The biotical RSV Positive Control is an external quality control for biotical RSV card. Use this control swab to check that the sample reagents and the test card are working properly.

biotical RSV Positive Control is manufactured with inactivated Respiratory Syncytial Virus antigens dried in swab containing a preservative.

POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by **rotating the positive control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes**. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE.</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT.</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE.</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL.</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- The test has been validated with samples of nasal swabs, washing and sucking nasopharyngeal. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of RSV in the sample. To confirm the presence of infection by RSV it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques (cellular culture and PCR).
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory techniques (cellular culture and PCR). A negative result does not indicate absence of infection by RSV. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity.

An evaluation, with nasal samples, was performed comparing the results obtained by an immunochromatographic test (biotical RSV, biotical) and another commercially available immunochromatographic test (BinaxNOW® RSV, Alere).

The results were as follows:

		BinaxNoW® RSV (immunochromatographic test)		
		+	-	Total
Biotical RSV Card	+	18	0	18
	-	1	10	11
	Total	19	10	29

<i>Biotical RSV Card vs BinaxNoW® RSV</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	94,7%	73,97 – 99,87
SPECIFICITY	>99%	69,15 – 100%
PPV	>99%	–
NPV	>90,9%	59,75 – 98,54%
ACCURACY	>96,5%	82,24% – 99,91%

The results showed a high sensitivity and specificity to detect RSV using biotical RSV card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical RSV; no cross reactivity against organisms that cause other respiratory infections.

Influenza type A

Influenza type B

Adenovirus

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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