ENGLISH - ESPAÑOL

Ref: RTB25APC  
RTB25APCN

biotical  
**Adeno  
Respiratory  
Card**

Rapid Test to detect  
Adeno Respiratory in  
card format

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**biotical health**  
ready to test

IU RTB25APC-RTB25APCN Rev. 03 / 2024-04

ISO 13485 CERTIFIED  
MEDICAL DEVICE  
QUALITY MANAGEMENT SYSTEM

# biotical Adeno Respiratory Card

ENGLISH

## WHAT IS BIOTICAL ADENO RESPIRATORY CARD?

biotical Adenovirus Respiratory card is a rapid immunochromatography test for the qualitative detection of *Adenovirus* antigens in nasal swab, washed or nasopharyngeal aspirate samples. Thanks to its easy handling and high sensitivity, the healthcare professional is able to make a presumptive diagnosis of infection caused by *Adenovirus*.

## INFECTIONS BY ADENOVIRUS

*Adenovirus* are a group of viruses that cause respiratory diseases such as colds, bronchiolitis or pneumonia. In children, Adenovirus usually cause infections in respiratory tract.

Adenovirus infections can occur at any age and they are spread when a person comes into contact with another individual or an object that is contaminated with the virus. There is a higher prevalence in children between six months and two years of age attending day-care centres.

Most *Adenovirus* infections are mild and have the typical symptoms of a cold: runny nose, sore throat and headache, fever, severe cough, swollen lymph nodes, feeling sick and conjunctivitis.

biotical Adeno Respiratory card is an immunochromatographic qualitative test, easy to use and with almost immediate results, for the detection of *Adenovirus* directly in nasopharyngeal swab samples, providing results in 10 minutes.

## TABLE OF REFERENCES


In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25APC	biotical Adeno Respiratory Card (+PC)	25 test with positive control, without negative control.
RTB25APCN	biotical Adeno. Respiratori Card (+PC +NC)	25 test with positive and negative control.

## HOW DOES BIOTICAL ADENOVIRUS RESP. CARD WORK?

---

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In the **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In the **RED** line (T): with mouse monoclonal antibodies against Respiratory Adenovirus.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a test solution (mouse monoclonal antibodies anti-Adenovirus) conjugated to red polystyrene latex and a control solution (specific binding protein) conjugated to green polystyrene latex, forming two conjugate coloured complexes.

### TEST LINE:

**If the sample is positive:** The antigens of the patient's sample will react with the red coloured conjugate complex (anti-Adenovirus monoclonal antibodies - red latex microspheres) and advance by capillarity through the membrane. The anti- Adenovirus antibodies present in the in the membrane of the test line (T) will capture the coloured test complex resulting in the red line.

**If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti- Adenovirus antibodies present in the in the membrane of the test line (T) will not capture the red coloured antigen-complex (not formed) and the red line will not appear.

### CONTROL LINE (C):

**Regardless of whether the sample is positive or not,** Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

## STORAGE AND STABILITY

---

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

## GENERAL PRECAUTIONS

---

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.

- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

## MATERIALS

---

### MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Sample diluent (Reagent B).
- Sampling swabs.
- Disposable pipettes.
- Testing tubes
- Positive Control.
- Negative Control (only RTB25APCN reference).

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Vortex.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

## SPECIMEN COLLECTION AND PREPARATION


---

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

### Specimen preparation (see illustration):

#### Nasal swab method:


- Remove the swab from its packing.
- Collect specimen with a sterile swab from one nostril.
- Insert the swab into the nostril to the nasopharynx, rotating against the nasal wall (to ensure swab contains cells as well as mucus).
- Repeat procedure using another nostril.
- Process the swab as soon as possible after collecting the specimen.

Nasopharyngeal aspirate method (suction device, sterile suction catheter):

**For Adults:**

- Place the irrigator up to the nose.
- Let the sterile saline water run into the nose (2.5 mL). It will run out the opposite side.
- Tilt and twist the irrigator side to side and up and down directing the water flow into all portions of the nasal cavity.
- Collect the wash in a clean specimen container, tilt the head forward and allow the water with mucus to run out of the nostril into the specimen container. Repeat the mucus collection for the other nostril and collect it into the same container.



**For children:**

- Use an aspiration bulb or bulb syringe to instill the saline water into one nostril leaning the children head.
- Aspirate the mix of mucus-saline water into the bulb and transfer it into a clean container.
- Repeat for the other nostril and transfer the fluid into the same specimen container.

## TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

Procedure A using nasal swab samples:


1 Add 15 drops of sample diluent

2 Put the swab into the tube, rotating 30 seconds and extract the liquid.

3 Add 4 drops in the circular window (S).

- Add **15 drops ①** of **sample diluent** and immediately put the swab into the tube.
- Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds ②**. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

**Procedure B using nasopharyngeal wash or aspirate samples:**


- Add **6 drops of the nasopharyngeal wash or aspirate ①** samples with a pipette and **9 drops of sample diluent ②** in a testing tube. **Mix with vortex at least 30 seconds ③** to homogenize. Best results are obtained when the specimen is vigorously extracted in the solution.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S. Use a separate pipette and test for each sample or control ④.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent B until seeing the liquid running through the reaction zone.

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID

INVALID

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	<b>NEGATIVE for Adeno Respiratory</b> Only the <b>GREEN</b> line appears across the results window in the site marked with the letter C (control line). No presence of Adeno Respiratory is detected in the sample.
2	C = GREEN LINE T = RED LINE	<b>POSITIVE for Adeno Respiratory</b> In addition to the <b>GREEN</b> control line, a <b>RED</b> line also appears in the site marked with the letter T (test line) in the results window. Adeno Respiratory presence is detected in the sample.
3-4	Any other result	<b>INVALID</b> Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

### QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references [RTB25APC](#) & [RTB25APCN](#)).


## POSITIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Adeno Respiratory Positive Control is an external quality control for biotal Adenovirus Respiratory card. Use this control swab to check that the sample reagents and the test card are working properly.

biotal Adeno Respiratory Positive Control is manufactured with inactivated Adenovirus Respiratory antigens dried in swab containing a preservative.

## POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by **rotating the positive control swab** forcefully against the side of the testing tube at least **30 seconds** **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S **3**.
- **Read the result at 10 minutes**. do not read the test result later than 10 minutes.

## POSITIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = RED LINE</p>	<p><b>RESULT: CONTROL POSITIVE</b></p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p><b>INVALID RESULT</b></p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

## NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p><b>RESULT: CONTROL NEGATIVE</b></p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p><b>RESULT: INVALID FOR NEGATIVE CONTROL.</b></p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

## LIMITATIONS OF THIS TEST

- The test must be carried out up to 2 hours after opening the bag that contains it.
- The test has been validated with nasal swabs and nasopharyngeal aspirate and washes samples.
- The use with samples from other sources (for example: urine, sputum or saliva) has not been tested and therefore its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of the pathogen in the sample. To confirm Adenovirus infection, it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient in order to issue a valid clinical diagnosis.
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days and confirm the absence of Adenovirus using PCR or cell culture assays. It is possible the appearance of false negatives if the sample does not reach the detection limit of the test.

## PERFORMANCE CHARACTERISTICS

### Clinical sensitivity and specificity

An evaluation, with nasal samples, was performed comparing the results obtained by an immunochromatographic test (biotical Adenovirus Respiratory Card, biotical), another commercial available immunochromatographic test (Adenovirus Respi, CorisBioConcept) and a commercial available Immunofluorescence test (PathoDx® Adenovirus, Remel).

The results were as follows:

		Adenovirus Respi (immunochromatographic test)		
		+	-	Total
Biotical Adeno Respiratory Card	+	20	0	20
	-	0	5	5
	Total	20	5	25

		Pathodx® Adenovirus (Immunofluorescence test)		
		+	-	Total
Biotical Adeno Respiratory Card	+	20	0	20
	-	0	5	5
	Total	20	5	25

<i>Biotical Adeno Respiratory Card vs Adenovirus Respi &amp; Pathodx® Adenovirus</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99%	83,2 - 100%
SPECIFICITY	>99%	47,8 - 100%
PPV	>99%	-
NPV	>99%	-
ACCURACY	>99%	86,3 - 100%

The results showed a high sensitivity and specificity to detect Adenovirus Respiratory using biotical Adeno Respiratory Card.

## CROSS REACTIVITY

---

An evaluation was performed to determine the cross reactivity of biotical Adeno Respiratory Card; no cross reactivity against organisms that cause other respiratory infections.

*Influenza type A*

*Influenza type B*

*Respiratory Syncytial Virus*

## SYMBOLS FOR IVD COMPONENTS & REAGENTS

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## REFERENCES / BIBLIOGRAFÍA /

---

MARCELA ECHAVARRIA, JOSE L. SANCHEZ, et al. "Rapid Detection of Adenovirus in Throat Swab Specimens by PCR during Respiratory Disease Outbreaks among Military Recruits", Journal of Clinical Microbiology, Feb. 2003, Vol. 41, No. 2, p. 810–812

MARILYN J. AUGUST AND ANN L. WARFORD; "Evaluation of a Commercial Monoclonal Antibody for Detection of Adenovirus Antigen" Journal of Clinical Microbiology, Nov. 1987, Vol. 25, No. 11; p. 2233-2235

Sierra de Guadarrama, 1. Ed. Dublín  
28830 San Fernando de Henares. Madrid  
Tel.: + 34 91 677 43 08

[info@biotical.es](mailto:info@biotical.es)  
[www.biotical.es](http://www.biotical.es)

biotical health S.L.U.