ENGLISH - ESPAÑOL

Ref: RTB25VRC
RTB25VRCN

biotical RSV + Adeno Resp. Card

Rapid Test to detect
RSV + Adeno Resp. in
card format

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biotical health

ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical

ENGLISH

RSV + Adeno Resp. Card

WHAT IS BIOTICAL RSV + ADENO RESP. CARD?

biotical RSV + Adeno Respiratory card is a rapid immunochromatography test for the simultaneous qualitative detection of Respiratory Syncytial Virus (RSV) and Respiratory Adenovirus antigens in nasopharyngeal swab samples. Thanks to its easy handling and high sensitivity, the healthcare professionals are able to make a presumptive diagnosis of infection caused by RSV and / or Respiratory Adenovirus.

INFECTIONS BY ADENOVIRUS AND SYNCYTIAL RESPIRATORY VIRUS (RSV)

The Respiratory Syncytial Virus (RSV) is one of the main respiratory diseases in children younger than 1 year, affecting the lungs, respiratory tracks and also presenting fever and cough. In the case of most healthy children, RSV causes like a cold that usually lasts between 5 and 7 days, but in some cases, it can be followed by bronchiolitis or pneumonia. They can be contracted by adults and can be serious in the elderly if it affects the lower respiratory tract.

Adenoviruses are a group of viruses that cause respiratory illnesses such as colds, bronchiolitis, or pneumonia. In children, Adenoviruses generally cause infections in the respiratory tract. We can find other clinical pictures such as gastroenteritis or conjunctivitis, depending on the Adenovirus serotype that causes the infection.

Adenovirus infections can occur at any age and are spread when a person comes into contact with another individual or an object that is contaminated with the virus. There is a higher prevalence in children between six months and two years old who attend nurseries.

Because both infections are viral, antibiotics are not useful to combat it, so it is very important that the health professional make a correct identification.

biotical RSV + Adeno Respiratory card, is a qualitative immunochromatographic test, easy to use and with almost immediate results, for the detection of RSV and / or Adenovirus directly in nasopharyngeal swab samples, providing results in 10 minutes.

TABLE OF REFERENCES


In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25VRC	biotical RSV + Adeno Resp. Card (+PC)	25 test with positive control, without negative control.
RTB25VRCN	biotical RSV + Adeno Resp. Card (+PC +NC)	25 test with positive and negative control.

HOW DOES BIOTICAL RSV+ADENOVIRUS RESPIRATORY CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In the T₁ **RED** line (RSV): with mouse monoclonal antibodies against RSV.
- In the C **GREEN** line (Control): with rabbit polyclonal antibodies against a specific protein.
- In the T₂ **RED** line (Adeno Resp.): with mouse monoclonal antibodies against Adenovirus.
- In the **sample area (S)**: the pad where the sample is placed, it has been dispensed: a preparation of reagents of the T1 line (mouse monoclonal antibodies against RSV) conjugated with latex of red polystyrene; with another preparation of reagents of the T2 line (mouse monoclonal antibodies against Adenovirus) conjugated with latex of red polystyrene; and with another preparation for the control line (specific binding protein) conjugated with green polystyrene latex; forming three conjugate coloured complexes.

TEST LINES (T₁ and T₂):

If the sample is positive: the antigens of the patient's sample will react with one (or both) of the red coloured conjugate complexes (anti-RSV monoclonal antibodies-latex red microspheres for the line T1 and anti-Adenovirus monoclonal antibodies-red latex microspheres for the line T2) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-RSV present in the T1 line and / or the antibodies anti-Adenovirus present in the T2 line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or these are present in a concentration below the detection limit), there will be no reaction with the red coloured complexes. The anti-RSV and anti-Adenovirus antibodies present in the T1 and T2 lines will not capture the coloured antigen-complexes and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Reagent B (sample diluent).
- Swabs for sample taking.
- Disposable pipettes.
- Testing tubes.
- Instructions for use.
- Positives Controls (RSV & Adeno Resp.).
- Negative Control (only RTB25VRCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Vortex.

The kit contains materials for the analysis of a maximum of 25 samples or controls.


SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Nasal swab method:


- Remove the swab from its packing.
- Collect specimen with a sterile swab from one nostril.
- Insert the swab into the nostril to the nasopharynx, rotating against the nasal wall (to ensure swab contains cells as well as mucus).
- Repeat procedure using another nostril.
- Process the swab as soon as possible after collecting the specimen.

Nasopharyngeal aspirate method (suction device, sterile suction catheter):

For Adults:

- Place the irrigator up to the nose.
- Let the sterile saline water run into the nose (2.5 mL). It will run out the opposite side.
- Tilt and twist the irrigator side to side and up and down directing the water flow into all portions of the nasal cavity.
- Collect the wash in a clean specimen container, tilt the head forward and allow the water with mucus to run out of the nostril into the specimen container. Repeat the mucus collection for the other nostril and collect it into the same container.

For children:

- Use an aspiration bulb or bulb syringe to instill the saline water into one nostril leaning the children head.
- Aspirate the mix of mucus-saline water into the bulb and transfer it into a clean container.
- Repeat for the other nostril and transfer the fluid into the same specimen container.

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

Procedure A using nasal swab samples:

Results window

RED line (Test RSV)
GREEN Line (Control)
RED line (Test Adeno)

Sample window (S)


1 Add 15 drops of sample diluent

2 Put the swab into the tube, rotating 30 seconds and extract the liquid.

3 Add 4 drops in the circular window (S).

- Add **15 drops ①** of **sample diluent** and immediately put the swab into the tube.
- Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds ②**. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 10 minutes**. The results will not be valid after 10 minutes.

Procedure B using nasopharyngeal wash or aspirate samples:

- Add **6 drops of the nasopharyngeal wash or aspirate ①** samples with a pipette and **9 drops of sample diluent ②** in a testing tube. **Mix with vortex at least 30 seconds ③** to homogenize. Best results are obtained when the specimen is vigorously extracted in the solution.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S. Use a separate pipette and test for each sample or control ④.
- **Read the results at 10 minutes**. The results will not be valid after 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS


(Please refer to the illustration below)


RSV NEGATIVE
ADENOVIRUS NEGATIVE

RSV POSITIVE
ADENOVIRUS NEGATIVE

RSV NEGATIVE
ADENOVIRUS POSITIVE

RSV POSITIVE
ADENOVIRUS POSITIVE

INVALID

INVALID

INVALID

INVALID

	biotical Test	Interpretation of the results
1	T1 = NO LINE C = GREEN LINE T2 = NO LINE	NEGATIVE FOR RSV & ADENO RESPIRATORY Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of RSV or Adeno Respiratory are detected in the sample.
2	T1 = RED LINE C = GREEN LINE T2 = NO LINE	POSITIVE FOR RSV – NEGATIVE FOR ADENO RESPIRATORY In addition to the GREEN control line, a RED line also appears in the site marked with the letter T1 (RSV detection line) in the results window. RSV presence is detected in the sample. No presence of Adeno Respiratory is detected.
3	T1 = NO LINE C = GREEN LINE T2 = RED LINE	POSITIVE FOR ADENO RESPIRATORY – NEGATIVE FOR RSV In addition to the GREEN control line, a RED line also appears in the site marked with the letter T2 (Adeno Respiratory detection line) in the results window. Adeno Respiratory presence is detected in the sample. No presence of RSV is detected.

biotical Test		Interpretation of the results
4	T1 = RED LINE C = GREEN LINE T2 = RED LINE	<p>POSITIVE FOR RSV & ADENO RESPIRATORY</p> <p>In addition to the GREEN control line, a RED line also appears in the site marked with the letter T1 (RSV detection line) and a RED line appears in the site marked with the letter T2 (Adeno Respiratory detection line) in the results window. RSV and Adeno Respiratory presence are detected in the sample.</p>
5-8	Any other result	<p>INVALID</p> <p>Total absence of the control line (GREEN) regardless the appearance or not of any test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.</p>

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T1 – T2) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits reference **RTB25VRCN**).

POSITIVE CONTROLS: INTENDED USE & COMPOSITION


The biotical RSV Positive Control and biotical Adeno Respiratory Positive Control are two external quality controls for biotical RSV + Adeno Resp. Card. Use these controls swabs to check that the sample reagents and the test card are working properly.

biotical RSV Positive Control: inactivated Respiratory Syncytial Virus antigens dried in swab containing a preservative.

biotical Adeno Respiratory Positive Control: inactivated Respiratory Adenovirus antigens dried in swab containing a preservative.


POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube ①.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

biotical Test	Interpretation of the results
 <p>T1 = RED LINE C = GREEN LINE T2 = NO LINE</p>	<p>RESULT: RSV POSITIVE CONTROL.</p> <p>A GREEN line (C) and a RED line (T1) appear across the results window during the test performance due to presence of the positive control swab.</p>
 <p>T1 = NO LINE C = GREEN LINE T2 = RED LINE</p>	<p>RESULT: Adeno Respiratory POSITIVE CONTROL.</p> <p>A GREEN line (C) and a RED line (T2) appear across the results window during the test performance due to presence of the positive control swab.</p>
	<p>INVALID RESULTS FOR POSITIVE CONTROL.</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of any of the test lines (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.


NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. Do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

biotical Test	Interpretation of the results
 <p>T1 = NO LINE C = GREEN LINE T2 = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE.</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
	<p>RESULT: INVALID FOR NEGATIVE CONTROL.</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- The test has been validated with samples of nasal swabs, washing and sucking nasopharyngeal. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of RSV/Adenovirus in the sample. To confirm the presence of infection by RSV/Adenovirus it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques (cellular culture and PCR).
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory techniques (cellular culture and PCR). A negative result does not indicate absence of infection by RSV/Adenovirus. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity (RSV):

An evaluation, with nasal samples, was performed comparing the results obtained by an immunochromatographic test (biotical RSV + Adeno Resp., biotical) and another commercially available immunochromatographic test (BinaxNOW® RSV, Alere).

The results for RSV were as follows:

		BinaxNOW® RSV, Alere (immunochromatographic)		
		+	-	Total
Biotical RSV + Adeno Resp. Card (T1: RSV)	+	18	0	18
	-	1	10	11
	Total	19	10	29

<i>Biotical RSV + Adeno Resp. Card vs BinaxNOW® RSV, Alere (RSV)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	94,5%	74,0% - 99,9%
SPECIFICITY	>99,9%	74,0% - 99,9%
PPV	>99,9%	-
NPV	90,9%	59,8% - 98,5%
ACCURACY	96,6%	82,2% - 99,9%

Clinical sensitivity and specificity (Adenovirus):

An evaluation, with nasal samples, was performed comparing the results obtained by an immunochromatographic test (biotical RSV + Adeno Resp., biotical), another commercially available immunochromatographic test (Adenovirus Respi, CorisBioConcept) and a commercially available Immunofluorescence test (PathoDx®Adenovirus, Remel).

The results for Adeno Respiratory were as follows:

		Adenovirus Respi, CorisBioConcept (immunochromatographic)		
		+	-	Total
Biotical RSV + Adeno Resp. Card (T2: ADENO RESPIRATORY)	+	20	0	20
	-	0	5	5
	Total	20	5	25

<i>Biotical RSV + Adeno Resp. Card vs Adenovirus Respi, CorisBioConcept (ADENOVIRUS RESPIRATORY)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99,9%	83.2% - 100.0%
SPECIFICITY	>99,9%	47.8% - 100.0%
PPV	>99,9%	-
NPV	>99,9%	-
ACCURACY	>99,9%	86.3% - 100.0%

		PathoDx®Adenovirus, Remel (Immunofluorescence)		
		+	-	Total
Biotical RSV + Adeo Resp. Card (T2: ADENO RESPIRATORY)	+	20	0	20
	-	0	5	5
	Total	20	5	25

<i>Biotical RSV + Adeno Resp. Card vs PathoDx®Adenovirus, Remel (ADENOVIRUS RESPIRATORY)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99,9%	83,2% - 100,0%
SPECIFICITY	>99,9%	47,8% - 100,0%
PPV	>99,9%	-
NPV	>99,9%	-
ACCURACY	>99,9%	86,9% - 100,0%

The results showed a high sensitivity and specificity to detect RSV and Respiratory Adenovirus using biotical RSV + Adeno Resp. Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical RSV + Adeno Resp Card; no cross reactivity against organisms that cause other respiratory infections:

For RSV:

Influenza type A

Influenza type B

Adenovirus

For Adeno Resp.:

Influenza type A

Influenza type B

RSV

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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MARCELA ECHAVARRIA, JOSE L. SANCHEZ, et al. "Rapid Detection of Adenovirus in Throat Swab Specimens by PCR during Respiratory Disease Outbreaks among Military Recruits", Journal of Clinical Microbiology, Feb. 2003, Vol. 41, No. 2, p. 810–812.

MARILYN J. AUGUST AND ANN L. WARFORD; "Evaluation of a Commercial Monoclonal Antibody for Detection of Adenovirus Antigen" Journal of Clinical Microbiology, Nov. 1987, Vol. 25, No. 11; p. 2233-2235.

ANN R. FALSEY and EDWARD E. WALSH. "Respiratory Syncytial Infection in Adults", Clinical Microbiology Reviews, July 2000, Vol. 13, No.3, p. 371-384.

DIANE C. HALSTEAD, SANDRA TODD, and GALE FRITCH; "Evaluation of Five Methods for Respiratory Syncytial Virus Detection" Journal of Clinical Microbiology, May 1990, Vol. 28, No. 5; p. 1021-1025.

Sierra de Guadarrama, 1. Ed. Dublín
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es

biotical health S.L.U.