ENGLISH - ESPAÑOL

Ref: RTB25ABC  
RTB25ABCN

biotical

# Influenza A+B Card

Rapid Test to detect  
Influenza A+B in card  
format

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biotical health  
ready to test

ISO 13485 CERTIFIED  
MEDICAL DEVICE  
QUALITY MANAGEMENT SYSTEM

# biotical Influenza A+B Card

ENGLISH

## WHAT IS BIOTICAL INFLUENZA A+B CARD?

biotical Influenza A+B card is a rapid immunochromatography test for the qualitative simultaneous detection of Influenza A and Influenza B antigens from nasal swab, nasopharyngeal wash or aspirate specimens. Thanks to its easy handling and high sensitivity, the healthcare professional can make a presumptive diagnosis of infection caused by Influenza A and / or Influenza B.

## WHAT IS THE INFLUENZA VIRUS?

Influenza is a highly contagious respiratory disease caused by Influenza viruses (or Influenzaviruses) that generally infect the upper respiratory tract, throat, and bronchi. Sometimes, although less commonly, it also affects the lungs. Although the infection is usually overcome in a week, there are certain populations that are at a higher risk for developing complications (children, the elderly, pregnant women, immunosuppressed).

A quick and accurate diagnosis means reducing the risk of transmission and outbreak, preventing inappropriate use of antibiotics and ensuring early prescription of antiviral medication for high-risk patients. biotical Influenza A+B card is an immunochromatographic qualitative test, easy to use and with almost immediate results, for the detection of Influenza virus A subtypes (including subtypes A / H1N1, A / H3N2, A / H5N1, etc.) e Influenza virus B directly from the patient's nasopharyngeal sample.


## TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25ABC	biotical Influenza A+B Card (+PC)	25 test with positive control, without negative control.
RTB25ABCN	biotical Influenza A+B Card (+PC +NC)	25 test with positive and negative control.

## HOW DOES BIOTICAL INFLUENZA A+B CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In **RED** T<sub>1</sub> line (Influenza A): with mouse monoclonal antibodies against Influenza A.
- In **GREEN** C line (Control): with rabbit polyclonal antibodies against a specific protein.
- In **RED** T<sub>2</sub> line (Influenza B): with mouse monoclonal antibodies against Influenza B.
- **Sample area (S)**: the pad where the sample is placed, it has been dispensed: a preparation of reagents of the T1 line (mouse monoclonal antibodies against Influenza A) conjugated with latex of red polystyrene; with another preparation of reagents of the T2 line (mouse monoclonal antibodies against Influenza B) conjugated with latex of red polystyrene; and with another preparation for the control line (specific binding protein) conjugated with green polystyrene latex; forming three conjugate coloured complexes.

### TEST LINES:

**If the sample is positive:** the antigens of the patient's sample will react with one (or both) of the red coloured conjugate complexes (anti-Influenza A monoclonal antibodies-latex red microspheres for the line T1 and anti-Influenza B monoclonal antibodies-red latex microspheres for the line T2) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-Influenza A present in the T1 line and / or the antibodies anti-Influenza B present in the T2 line will capture the complex giving rise in a red line.

**If the sample is negative:** if the sample is negative and therefore does not contain antigens (or these are present in a concentration below the detection limit), there will be no reaction with the red coloured complexes. The anti-Influenza A and anti-Influenza B antibodies present in the T1 and T2 lines will not capture the coloured antigen-complexes and the red lines will not appear.

### CONTROL LINE (C):

**Regardless of whether the sample is positive or not,** to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

## STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

## GENERAL PRECAUTIONS

---

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

## MATERIALS

---

### MATERIALS PROVIDED

- Cassette or Test Card.
- Reagent B (sample diluent).
- Swabs for sample taking.
- Disposable pipettes.
- Testing tubes.
- Instructions for use.
- Positive Control (Influenza A+B in one swab).
- Negative Control (only RTB25ABCN reference).

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Vortex.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

## SPECIMEN COLLECTION AND PREPARATION


---

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

### Specimen preparation (see illustration):

#### Nasal swab method:

- Remove the swab from its packing.
- Collect specimen with a sterile swab from one nostril.
- Insert the swab into the nostril to the nasopharynx, rotating against the nasal wall (to ensure swab contains cells as well as mucus).
- Repeat procedure using other nostril.
- Process the swab as soon as possible after collecting the specimen.

#### Nasopharyngeal aspirate method (suction device, sterile suction catheter):

##### For Adults:

- Place the irrigator up to the nose.
- Let the sterile saline water run into the nose (2.5 mL). It will run out the opposite side.
- Tilt and twist the irrigator side to side and up and down directing the water flow into all portions of the nasal cavity.
- Collect the wash in a clean specimen container, tilt the head forward and allow the water with mucus to run out of the nostril into the specimen container. Repeat the mucus collection for the other nostril and collect it into the same container.


##### For children:

- Use an aspiration bulb or bulb syringe to instil the saline water into one nostril leaning the children head.
- Aspirate the mix of mucus-saline water into the bulb and transfer it into a clean container.
- Repeat for the other nostril and transfer the fluid into the same specimen container.

## TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

### Procedure A using nasal swab samples:

- Add **15 drops** ① of sample diluent and immediately put the swab into the tube.
- Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds** ②. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

### Procedure B using nasopharyngeal wash or aspirate samples:


- Add 6 drops of nasopharyngeal wash or aspirate **1** with a pipette and 9 drops of sample diluent **2** in a test tube. Mix with vortex for at least 30 seconds **3** to homogenize. Better results are obtained if the sample is extracted vigorously.
- Remove the test card from its package just before use.
- Dispense 4 drops from the test tube into the circular window marked with the letter S. Use a separate pipette and card for each sample or control **4**.
- Read the results after 10 minutes. The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	T1 = NO LINE C = GREEN LINE T2 = NO LINE	<b>NEGATIVE FOR INFLUENZA A &amp; INFLUENZA B</b> Only the <b>GREEN</b> line appears across the results window in the site marked with the letter C (control line). No presence of influenza A or influenza B are detected in the sample.

biotical Test		Interpretation of the results
2	T1 = RED LINE C = GREEN LINE T2 = NO LINE	<b>POSITIVE FOR INFLUENZA A – NEGATIVE FOR INFLUENZA B</b> In addition to the GREEN control line, a RED line also appears in the site marked with the letter T1 (influenza A detection line) in the results window. Influenza A presence is detected in the sample. No presence of influenza B is detected.
3	T1 = NO LINE C = GREEN LINE T2 = RED LINE	<b>POSITIVE FOR INFLUENZA B – NEGATIVE FOR INFLUENZA A</b> In addition to the GREEN control line, a RED line also appears in the site marked with the letter T2 (influenza B detection line) in the results window. Influenza B presence is detected in the sample. No presence of influenza A is detected.
4	T1 = RED LINE C = GREEN LINE T2 = RED LINE	<b>POSITIVE FOR INFLUENZA A &amp; INFLUENZA B</b> In addition to the GREEN control line, a RED line also appears in the site marked with the letter T1 (influenza A detection line) and a RED line appears in the site marked with the letter T2 (INFLUENZA B detection line) in the results window. Influenza A and influenza B presence are detected in the sample.
5-8	Any other result	<b>INVALID</b> Total absence of the control line (GREEN) regardless the appearance or not of any test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T1 – T2) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL & NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25ABC & RTB25ABCN**).

## POSITIVE CONTROLS: INTENDED USE & COMPOSITION

biotical Influenza A+B Positive Control is an external quality control for biotical Influenza A+B card. Use this control swab to check that the sample reagents and the test card are working properly.

biotical Influenza A+B Positive Control: inactivated Influenza A+B antigens dried in swab containing a preservative.

## POSITIVE CONTROL: TESTING PROCEDURE


Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- 1 Add 15 drops of diluent and introduce the positive control swab into the testing tube.
- 2 Rotating 30 seconds and extract the liquid.
- 3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. Do not read the test result later than 10 minutes.

## POSITIVE CONTROL – INTERPRETATION OF RESULTS

biotical Test	Interpretation of the results
 <p>T1 = RED LINE C = GREEN LINE T2 = RED LINE</p>	<p><b>RESULT: INFLUENZA A+ B POSITIVE CONTROL</b></p> <p>A GREEN line (C) and two RED lines (T1 – T2) appear across the results window during the test performance due to presence of the positive control swab.</p>
	<p><b>INVALID RESULTS FOR POSITIVE CONTROL</b></p> <p>Total absence of the control line (GREEN), regardless the appearance or not of any of the test lines (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

## NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1 Add 15 drops of diluent and introduce the negative control swab into the testing tube.


2 Rotating 30 seconds and extract the liquid.

3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube 1.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. Do not read the test result later than 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

biotical Test	Interpretation of the results
 <p>T1 = NO LINE C = GREEN LINE T2 = NO LINE</p>	<p><b>RESULT: CONTROL NEGATIVE</b></p> <p>Only the control line (C) appears across the results window during the test performance.</p>

**RESULT: INVALID FOR NEGATIVE CONTROL**

Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.

Review the procedure and repeat the control with a new test and a new negative control swab.

## LIMITATIONS OF THIS TEST

- The test must be carried out within 2 hours after opening the sealed bag.
- biotical Influenza A+B card should be used only with nasal swab, nasopharyngeal wash and aspirate samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasal specimens must be obtained.
- The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- Positive results determine the presence of type A and/or type B respiratory infection. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the nasal samples is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture or PCR.

## PERFORMANCE CHARACTERISTICS

### Clinical sensitivity and specificity:

An evaluation was performed comparing the results obtained by an immunochromatographic test (biotical Influenza A+B, biotical) and another commercially available immunochromatographic test (Binax-NOW® Influenza A&B, Alere).

The results for RSV were as follows:

		Binax-NOW® Influenza A&B, Alere (immunochromatographic)		
		+	-	Total
Biotical Influenza A + B Card	+	5	0	5
	-	0	6	6
	Total	5	6	11

<i>Biotical Influenza A+B Card vs Binax-NOW® Influenza A&amp;B, Alere</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99,9%	47,8% - 100,0%
SPECIFICITY	>99,9%	54,1% - 100,0%
PPV	>99,9%	-
NPV	90,9%	-
ACCURACY	96,6%	71,5% - 100,0%

The results showed a high sensitivity and specificity to detect Influenza type A and B using biotical Influenza A+B Card.

## CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Influenza A+B; no cross reactivity against organisms that cause other respiratory infections:

For Influenza type A:

*Syncytial Respiratory Virus*

*Influenza type B*

For Influenza type B:

*Syncytial Respiratory Virus*

*Influenza type A*

## SYMBOLS FOR IVD COMPONENTS & REAGENTS.

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---

Daniel E. Noyola, et al. Comparison of a New Neuraminidase Detection Assay with an Enzyme Immunoassay, Immunofluorescence, and Culture for Rapid Detection of Influenza A and B Viruses in Nasal Wash Specimens. *Journal of Clinical Microbiology*, Mar. 2000, p. 1161–1165 Vol. 38, No. 3

Adriana Weinberg and Miranda L. Walker. Evaluation of Three Immunoassay Kits for Rapid Detection of Influenza virus A and B. *Clinical And Diagnostic Laboratory Immunology*, Mar. 2005, p. 367–370 Vol. 12, No. 3

Andreea C. Cazacu, Gail J. Demmler et al. Comparison of a New Lateral-Flow Chromatographic Membrane Immunoassay to Viral Culture for Rapid Detection and Differentiation of Influenza A and B Viruses in Respiratory Specimens. *Journal of Clinical Microbiology*, Aug. 2004, p. 3661–3664 Vol. 42, No. 8.

Kenneth E. Irmen and James J. Kelleher. Use of Monoclonal Antibodies for Rapid Diagnosis of Respiratory Viruses in a Community Hospital. *Clinical and Diagnostic Laboratory Immunology*, May 2000, p. 396–403 Vol. 7, No. 3.

Sierra de Guadarrama, 1. Ed. Dublín  
28830 San Fernando de Henares. Madrid  
Tel.: + 34 91 677 43 08

[info@biotical.es](mailto:info@biotical.es)  
[www.biotical.es](http://www.biotical.es)

biotical health S.L.U.