ENGLISH - ESPAÑOL

Ref: RTB25CoV2C
RTB25CoV2CN

biotical SARS CoV-2 Ag Card

Rapid Test to detect
SARS CoV-2 Ag in card
format

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biotical health

ready to test

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical SARS CoV-2 Ag Card

ENGLISH

WHAT IS BIOTICAL SARS COV-2 AG CARD?

The biotical SARS-CoV-2 Ag Card is a rapid immunochromatography test for the qualitative detection of the structural nucleoprotein antigen (N protein) of the SARS-CoV-2 coronavirus, from nasopharyngeal swab samples and nasal swabs samples from patients with suspected infection. by COVID-19.

Thanks to its easy handling and high sensitivity, the health professional or the layman in the field will be able to make a presumptive diagnosis of infection caused by the SARS-CoV-2 coronavirus.

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
TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25COV2C	biotical SARS CoV-2 Ag Card (+PC)	25 test with positive control, without negative control.
RTB25COV2CN	biotical SARS CoV-2 Ag Card (+PC +NC)	25 test with positive and negative control.

HOW DOES BIOTICAL SARS COV-2 Ag CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against SARS-CoV-2 N antigen.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of test line reagents (antibodies versus SARS-CoV-2 Ag) conjugated with red polystyrene latex and other preparation for the control line (specific binding protein) conjugated with green polystyrene latex has been dispensed, forming two coloured conjugate complex.

TEST LINE:

If the sample is positive: the antigens of the SARS-CoV-2 virus, from the patient's nasopharyngeal or nasal sample will react with the red conjugate complex (anti-SARS-CoV-2-red latex microspheres) and capillarity advance through the membrane. Anti-SARS-CoV-2 antibodies present in the test line membrane (T) will capture the coloured complex of the test resulting in the red line.

If the sample is negative: if the sample is negative and therefore does not contain viral antigens (or the antigen is present at a concentration below the detection limit), there is no reaction with the red coloured complex. Anti-SARS-CoV-2 antibodies present in the test line membrane (T) will not capture the red (unformed) coloured antigen-complex and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test was successful (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the antibodies immobilized against the specific protein, located on the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.

- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Monodose sample vials (with sample diluent).
- Sampling swabs.
- Lab rack.
- Positive Control.
- Negative Control (only RTB25COV2CN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Scissors.

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION.


Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8°C) for 8 hours prior to testing.

Select the type of sample to be used. This kit admits the use of nasopharyngeal and nasal samples. Although it has been shown that the nasopharyngeal sample is the one that contains a greater viral load, its extraction properly is much more complex and is only recommended for the scope of professional diagnosis, so it should always be a health professional who extracts it.

Specimen preparation (see illustration):


Nasal samples:

- Sound the nose before extraction to clean the nostrils.
- Carefully take out the sampling swab of its container.
- Insert the swab in vertical position by the nasal hole from 2 to 3 cm. Rotate the swab rubbing with the cotton part by the internal walls of the nasal cavity for at least 10 seconds to collect cells and snot.
- Remove the swab carefully and repeat the procedure with the other nasal orifice and the same swab.
- Remove the swab carefully. Do not let the cotton tip contact any surface you can contaminate. Perform the test with the extracted sample as soon as possible.

Nasopharyngeal samples (Recommended professional use):

- Carefully take out the sampling swab of its container.
- Introduce it by a nasal hole towards the nasopharynx, through the floor of the nasal cavity.
- Once contacted with the nasopharynx (the swab should be entered 3-5cm until touching the end of the nasal cavity), rotate the swab gently to capture both cells and snot.
- Repeat, if necessary, the procedure with the other nasal orifice.
- Process the swab as soon as possible, after the collection of the sample.



A correct sampling and conservation is essential for obtaining a correct result.

TEST PROCEDURE


It is necessary to temper the components of the test (sample diluent and cassette with the reactive strip) at room temperature (15 to 30°C). It is recommended not to open the containers until the moment of performing the test.

**Familiarize yourself with the components of the test:
sample vial and cassette, before performing the test.**

CASSETTE OF TEST

VIAL FOR SAMPLE

1. Open the **pre-dosed vial with sample diluent**, and immediately introduce the swab with the corresponding sample, carefully not to contaminate any external surface with the sample.
2. Rotate the swab at the bottom of the sample diluent to allow it to be mixed, for **at least 30 seconds**. Better results will be obtained if the mixture is made vigorously, but extreme precautions to avoid splashing.
3. Once mixed, extract the swab trying to **squeeze the highest volume possible of liquid**, pressing it against the tube walls or cut it with scissors leaving the cottony end inside the vial.

Discard the swab and close the sample vial with blue and white plugs.

4. Remove the **white stopper** to release access to the dropper and remove the SARS-COV-2 Ag test from its aluminium envelope and dispense in the sample window (marked with the letter S), **3 drops** of the sample, pressing softly and continuously the side of the vial of samples.

Close the sample vial with white stopper.


5. Wait **10-15 minutes to read the results**. The results will not be considered valid overcome this time.

NOTE: Sample diluent presents an inactivation capacity of the SARS Virus COV-2 of 99,68% and 99,98% after 1 and 10 minutes of contact, respectively. Even so, extreme precautions and treat the swab used and the diluent and vial for samples as if they contain biohazard species.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE FOR SARS CoV-2 Ag Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of SARS COV-2 Ag is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE FOR SARS CoV-2 Ag In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. SARS COV-2 Ag presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kit reference **RTB25CoV2CN**).


POSITIVE CONTROL: INTENDED USE & COMPOSITION

The biotical SARS COV-2 Ag Positive Control is an external quality control for biotical SARS COV-2 Ag card. Use this control swab to check that the sample reagents and the test card are working properly.

biotical SARS COV-2 Ag Positive Control is manufactured with inactivated SARS CoV-2 nucleoprotein antigens dried in swab containing a preservative.

POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Open the pre-dosed vial with sample diluent and **introduce the control swab**. Rotate the swab at the bottom of the sample diluent **for at least 30 seconds**.
- **Extract the swab**, trying to squeeze the highest volume possible of liquid, pressing it against the tube walls. Discard the swab and close the sample vial with blue and white plugs.
- Remove the white stopper to access the dropper and remove the test card from its aluminium envelope. Dispense in the sample window (marked with the letter S) **3 drops of the control** pressing softly and continuously the side of the vial. Close the white plug.
- **Wait 10 minutes to read the results**. The results will not be considered valid overcome this time.

POSITIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Open the pre-dosed vial with sample diluent and **introduce the control swab**. Rotate the swab at the bottom of the sample diluent **for at least 30 seconds**.
- **Extract the swab**, trying to squeeze the highest volume possible of liquid, pressing it against the tube walls. Discard the swab and close the sample vial with blue and white plugs.
- Remove the white stopper to access the dropper and remove the test card from its aluminium envelope. Dispense in the sample window (marked with the letter S) **3 drops of the control** pressing softly and continuously the side of the vial. Close the white plug.
- **Wait 10-15 minutes to read the results**. The results will not be considered valid overcome this time

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test should be performed up to 2 hours after the containing bag has been opened.
- Only nasal and nasopharyngeal samples freshly collected or preserved in transport media VTM, UTM and saline buffer have been validated. The use of other samples (oropharyngeal swabs, saliva, sputum,

urine, etc.) or other storage buffers has not been determined. The validity of the results depends directly on a correct collection and conservation of the sample. The use of storage buffers or their inappropriate use can reduce the sensitivity of the test due to excessive dilutions of the original sample or incompatibilities with the selected sample buffer

- biotical SARS-CoV-2 Ag card will only indicate the presence of SARS-CoV-2 in the sample (qualitative detection) and should be used only for the detection of SARS-CoV-2 antigens in nasopharyngeal samples. This test cannot determine either the quantitative value or rate of increase in SARS-CoV-2 antigen concentration.
- If the test result is negative and clinical symptoms persist, additional tests with other clinical methods are recommended. A negative result does not at any time prevent the possibility of SARS-CoV-2 infection as it may be due to a very low antigen concentration.
- This test provides a presumptive diagnosis of SARS-CoV-2 infections. All results should be interpreted by a physician along with other clinical information and other available laboratory results.
- False negative and positive results may occur in symptomatic and asymptomatic patients, within the sensitivity and specificity values provided in this manual.
- As a confirmatory diagnostic test, it is recommended to perform a Real Time-PCR.
- Any positive or negative result of the biotical SARS CoV-2 Ag rapid test does not identify or rule out co-infections with other pathogens.
- For professional in vitro use only. Exceptionally, due to the increase in cases of COVID infection, the Spanish Agency for Medicines and Medical Devices (AEMPS) has issued express authorization for the commercialization and use for its use as a self-diagnosis of this test, through sale in pharmacies, with a certain duration (authorization 1001-21-AE).
- The supplied sterile swabs are for the sole and exclusive use of obtaining nasal or nasopharyngeal samples. They must not be reused. It is convenient not to touch the head of the swab with the hands or other surfaces to avoid contamination.

WORLD HEALTH ORGANIZATION (WHO) RECOMMENDATIONS

The World Health Organization (WHO) has published on September 11, 2020 a guide aimed at improving the use and understanding of rapid tests for the detection of SARS CoV-2 antigen that can help organizations and countries interested in integrating this type of test into their diagnosis and surveillance system against SARS CoV-2, and especially its use in the identification of new outbreaks.

These recommendations include:

- To ensure the performance of the test, it must be carried out by qualified personnel, strictly following the instructions for use. The test should preferably be performed between days 5 and 7 from the onset of symptoms, this being the period with the highest viral load.
- Whenever possible, positive results must be confirmed by PCR techniques. A positive result has a high chance of being confirmed by these techniques, so it is advisable to isolate the patient while waiting for results.
- A negative result should not completely rule out SARS CoV-2 infection, so it is advisable to carry out confirmatory tests, especially in individuals who show symptoms. In special detection situations, such as the screening of new outbreaks and if the capacity for confirmation by PCR is insufficient, confirmatory tests should be prioritized over samples with a negative result.
- Antigen detection tests are especially important in cases: early detection of infections and identification of new outbreaks, monitoring of contagion trends, community transmission, analysis of contacts close to infections and risk personnel who may asymptotically carry the disease (It has been shown that the viral load of asymptomatic patients can be equivalent to that of symptomatic patients).
- SARS CoV-2 antigen detection tests work satisfactorily with samples with high viral load ($Ct \leq 25$ or $>10^6$ genome copies per mL). One way of evaluating the sensitivity for this class of tests is to check if the value of this is higher than 90% for samples with a Ct cycle threshold equal or lower than 28 or 30, since the virus concentration is expected to be abundant in this kind of samples.

You can consult the full text of the report published by the World Health Organization at the following address: <https://www.who.int/publications/i/item/antigen-detection-in-the-diagnosis-of-sars-cov-2infection-using-rapid-immunoassays>.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity: Test detection limit

The detection limit of the SARS CoV-2 Ag biotical rapid test has been determined at 1.25x10² UFP/mL nCoV-2019 D614G(S) or 2ng/mL of recombinant viral protein.

Due to the peculiarities of the technique in the processes related to the manufacture of biological raw materials, it is impossible to establish an exact guaranteed detection limit, therefore a standard detection limit of 2ng/mL is determined and an acceptance range between 4 ng/mL and 0.5 ng/mL.

Clinical sensitivity and specificity for nasopharyngeal samples

Multi-center validation with nasopharyngeal samples was performed and comparing the results obtained by this test (immunochromatographic) with a PCR test (reference technique). Samples employed had a value of CT<28 (nasopharyngeal samples with high viral load, being the expected in the first phases in which symptoms of the disease appear, days 5-7 from the beginning of the disease).

The results were as follows:

		Reference Technique (qPCR Test)		
		+	-	Total
Biotical SARS CoV-2 Ag Card (nasopharyngeal samples)	+	93	1	94
	-	7	455	462
	Total	100	456	556

<i>Biotical SARS COV-2 Ag Card (nasopharyngeal samples) vs Reference Technique (qPCR)</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	93,0%	86,1 - 97,1%
SPECIFICITY	99,8%	98,8 - 100%
PPV	98,9%	94,2 - 100%
NPV	98,5%	96,9 - 99,4%
ACCURACY	98,6%	97,2 - 99,4%

Clinical sensitivity and specificity for nasal samples

A validation with nasal samples was carried out and comparing the results obtained by this test (immunochromatographic) with a PCR test (reference technique). In this case, the samples used has not been discriminated on according to its CT value (Samples have been included with CT greater than 30, with very low viral load, for which the use of this Test class has not been recommended).

The results were as follows:

		Reference Technique (qPCR Test)		
		+	-	Total
Biotical SARS CoV-2 Ag Card (nasal samples)	+	129	7	136
	-	27	827	854
	Total	156	834	990

Biotical SARS COV-2 Ag Card (nasal samples) vs Reference Technique (qPCR)	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	82,7%	75,8 - 88,3%
SPECIFICITY	99,2%	98,3 - 99,7%
PPV	94,9%	89,7 - 97,9%
NPV	96,8%	95,4 - 97,9%
ACCURACY	96,6%	95,2 - 97,6%

If we only select the samples with CT<28, for which the use of this test class is recommended by the World Health Organization, the sensitivity obtained increases to 96,9% (IC: 91,2 – 99,4%) and the specificity to 99,2% (IC: 98,3 – 99,7%).

The results obtained both with nasal sample and nasopharyngeal sample, indicate a high sensitivity and specificity of the test, with very high general precision values.

HOOK EFFECT

Biotical SARS COV-2 AG rapid test showed no inhibitory effect at higher concentrations 100,000xLoD.

INTERFERING SUBSTANCES

No interference has been found in the results of the Biotical Rapid Test SARS COV Ag, with the following substances and concentrations.

INTERFERING SUBSTANCE	CONCENTRATION (mg/ML)	RESULT
Metronidazole	3,0	Negative
Ampicilin	3,0	Negative
Oseltamivir	3,0·10 ⁻²	Negative
Amantadine	0,3	Negative
Ribavirin	3,0	Negative
Codeine (Toseina)	0,2	Negative
Benzocaine (Angileptol)	3,0·10 ⁻²	Negative

INTERFERING SUBSTANCE	CONCENTRATION (mg/ML)	RESULT
Cloperastine (Flutox)	0,3	Negative
Carbocisteine (Iniston mucolítico)	3,0·10 ⁻²	Negative
Loratadine	0,3	Negative
Dexchlorpheniramine (Polaramine)	0,3	Negative
Loperamide hydrochloride (Fortasec)	0,15	Negative
Heparin (Hibor)	350,0 (IU/mL)	Negative
Almagato (Almax)	3,0	Negative

INTERFERING SUBSTANCE	CONCENTRATION (mg/ML)	RESULT
Fosfamyacin (Monurol)	$3,0 \cdot 10^{-3}$	Negative
Acetylcysteine (Fluimucil)	$3,0 \cdot 10^{-2}$	Negative
Dexketoprofen trometamol (Enantyum)	3,0	Negative
Hemoglobina humana	5,0 µg/mL	Negative
Calprotectina humana	5,0 µg/mL	Negative
Transferrina humana	0,5 µg/mL	Negative
Levofloacin	$3,0 \cdot 10^{-2}$	Negative
Ciprofoxacin	0,3	Negative
Rifampicin (Rifaldin)	0,3	Negative
Phenoxymethylpenicillin potasium	3,0	Negative
Ambroxol hydrochloride (Mucosan)	0,3	Negative
Ebastine (Ebastel)	3,0	Negative
Acetyl salicylic (Adiro)	0,3	Negative
Ibuprofen (espidifren)	0,3	Negative
Paracetamol (Dolocatil)	5,0	Negative
Metamizole (Nolotil)	5,0	Negative
Prednisone	0,3	Negative
Omeprazole	$2,0 \cdot 10^{-2}$	Negative
Macrogol 3350 (Movicol)	3,0	Negative
Lysine Carbocysteinate (Pectox)	$3,0 \cdot 10^{-2}$	Negative
Hydroxyzine dihydrochloride	0,3	Negative
Lorazepam	$3,0 \cdot 10^{-3}$	Negative
Amoxicillin	3,0	Negative
Mercaptopurine	0,3	Negative
Biotine	100,0 µg/mL	Negative
Human lactoferrin	5,0 µg/mL	Negative
Mucina	5,0	Negative
Human blood	50,0	Negative

CROSSED REACTIONS

An evaluation was carried out to determine possible cross-reactions of Biotical Sars-COV-2 AG; Cross reaction was not detected with organisms, pathogens or substances that could cause infections.

<i>Adenovirus</i>	<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Shigella dysenteriae</i>
<i>Adenovirus</i>	<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Shigella dysenteriae</i>
<i>Astrovirus</i>	<i>Escherichia coli</i> O157	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
Human calprotectin	<i>Giardia</i>	<i>Norovirus</i> GI	<i>Shigella sonnei</i>
<i>Campylobacter jejuni</i>	<i>Helicobacter pylori</i>	<i>Norovirus</i> GII	<i>Streptococcus pneumococcal</i>
<i>Clostridium difficile</i> antigen GDH	Bovine hemoglobin	<i>Rotavirus</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile</i> toxin A	Pork hemoglobin	<i>Salmonella enteritidis</i>	Human transferrin
<i>Clostridium difficile</i> toxin B	Human hemoglobin	<i>Salmonella paratyphi A</i>	Respiratory Syncytial Virus
Coronavirus strain 229E	<i>Influenza</i> type A	<i>Salmonella typhi</i>	<i>Yersinia enterocolitica</i> O:3
Coronavirus strain NL63	<i>Influenza</i> type B	<i>Salmonella typhimurium</i>	<i>Yersinia enterocolitica</i> O:9
Coronavirus strain OC43	Human lactoferrin	<i>Shigella boydii</i>	

Five positive samples for HKU1 human Coronavirus were tested in an external evaluation carried out with the Biotical SARS-COV-2 AG Card and the results obtained were negative. This suggests that there is no crossed reaction with HKU1. In addition, the percentage of similarity of the HKU1 nucleocapside protein sequence with the protein sequence of the SARS-COV-2 nucleocapside is approximately 35%, which is considered a very low homology. This suggests that there is no crossed reaction with HKU1.

Three positive samples for Rhinovirus were tested in an external evaluation carried out with the Certest Sars-COV-2 test and the results obtained were negative. This suggests that there is no crossed reaction with rhinovirus.

Biotical Sars-COV-2 AG shows some crossed reaction with SARS and very low with MERS.

NEW VARIANTS DISCOVERED

Recently new strains of the SARS CoV-2 virus, as the variants appeared in the United Kingdom, Brazil or Japan have been discovered. These variants have mutations in the gene sequences encoding for the SPIKE protein (S), so they have variations in these proteins that generally affect their capacity for infection.

As a result of these variations, the diagnostic tests that focus their detection in this gene can see their sensitivity capacity altered compared to these strains. Biotical SARS CoV-2 Ag Card specifically detects antigenic regions of the nucleocapside protein (NP) so that their sensitivity should not be altered for the detection of these and other new variants of the virus with mutations in protein S.

In addition, analysis of sequence in silico has been carried out to verify that the mutated regions of CDA one of the most widespread variants is outside the Biotical Test SARS CoV-2 Ag detection region.

VARIANT	ORIGIN	MUTATIONS IN PROTEIN N (POSITIONS)	AFFECTS THE TARGET REGION OF THE TEST? (247-364)	CONCLUSION
ALPHA	United Kingdom (VUI 202012/01)	3 (D3L) 235 (S235F)	NO	There should not be incompatibilities of detection.
BETA	South Africa (501Y.V2)	205 (T205I)	NO	There should not be incompatibilities of detection.
GAMMA	Brazil (GR/501Y.V3 - B.1.1.28)	80 (P80R) 203 (R203K) 204 (G204R)	NO	There should not be incompatibilities of detection.
DELTA	India (VUI-21APR-01 - B.1.617.2)	63 (D63G) 203 (R203M) 377 (D377Y)	NO	There should not be incompatibilities of detection.
EPSILON	USA GH/452R.V1 (B.1.429)	205 (T205I)	NO	There should not be incompatibilities of detection.
OMICRON	South Africa GR/484A.V2 (B.1.1.529)	13 (P12L), 203 (R203K) 204 (G204R), and deletions in 31, 32 and 33 (Deletions E, R and S)	NO	There should not be incompatibilities of detection.

In all cases, the result was that none of the analysed variants should alter the correct functioning of the test

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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REF

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REFERENCES / BIBLIOGRAFÍA /

Palestino, G.; García-Silva, I.; González-Ortega, O. and Rosales-Mendoza, S. Can nanotechnology help in the fight against COVID-19? Expert Review of Anti-infective Therapy. <https://doi.org/10.1080/14787210.2020.1776115>. Jun. 2020.

Moeis MR.; Rahayu AP.; Ihsani N. and Pertiwi W. Coronavirus-SARS-CoV-2: Biology and Problems in rRT-PCR Detection. Borneo Journal of Pharmacy. Vol. 3; special Issue 1:136-145. Jun. 2020. <https://doi.org/10.33084/bjop.v3iSpecial-1.1429>

Arnedo-Pena, A.; Sabater-Vidal, S.; Meseguer-Ferrer, N.; Pac-Sa, M. R.; Mañes-Flor, P.; Gascó-Laborda, J. C.; Larrea, R. M.; Tirado-Balaguer, M. D.; Romeu-García, M. A.; Gil-Fortuño, M.; Safont-Adsuara, L.; Blasco, A.; Gomila-Sar, B.; Moreno-Muñoz, M. R. and Bellido-Blasco, J. COVID-19 secondary attack rate and risk factors in household contacts in Castellon (Spain): Preliminary report. Revista Enfermedades Emergentes. 2020;19(2):64-70.

Younes, N.; Al-Sadeq, D. W.; AL-Jighefee, H.; Younes, S.; Al-Jamal, O.; Daas, H. I.; Yassine, H. M. and Nasrallah, G. K. Challenges in Laboratory Diagnosis of the Novel Coronavirus SARS-CoV-2. Viruses (Multidisciplinary Digital Publishing Institute). 12 (6): 582. <https://doi.org/10.3390/v12060582>. May. 2020.

Rincón, A.; Moreso, F.; López-Herradón, A.; Fernández-Robres, M.A.; Cidraque, I.; Nin, J.; Méndez, O.; López, M.; Pájaro, C.; Satorra, A.; Stuard, S. and Ramos, R. The keys to control a coronavirus disease 2019 outbreak in a haemodialysis unit. Clinical Kidney Journal. 2020: 1–8.

Fragkou, P. C.; Papaevangelou, V.; Antoniadou, A.; Kavvatha, D.; Ploussi, A.; Pantazis N.; Sirmipilantze, T.; Psarrakis, C.; Pournaras, S. A.; Tsiodras, S. and Kelekis, A. Preliminary Data of a Quantitative Point of Care Test for SARS-CoV-2 Antibodies From Greece. In vivo. 34: xxx-xxx (2020).

Ezhilan, M.; Suresh, I. and Nesakumar, N. SARS-CoV, MERS-CoV and SARS-CoV-2: A Diagnostic Challenge. Measurement. 168 (2021): 108335. <https://doi.org/10.1016/j.measurement.2020.108335>. Aug. 2020.8.

Antigen detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays. Interim guidance. 11 September 2020. World Health Organization.

Public Health England. Investigation of novel SARS-CoV-2 variant: Technical briefing 1. <https://www.gov.uk/government/publication/investigation-of-novel-sars-cov-2-variant-of-concern-20201201> (Accessed 12/01/2021).

Tegally H et al. Emergence and rapid spread of a new severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2) lineage with multiple spike mutations in South Africa. medRxiv 2020. F01: 10.1101/2020.12.21.20248640.

Sierra de Guadarrama, 1. Ed. Dublín
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es

biotical health S.L.U.