biotical

# S. pneumoniae Card

Rapid Test to detect  
S. pneumoniae in card  
format

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ENGLISH

## WHAT IS BIOTICAL S. PNEUMONIAE CARD?

biotical S. pneumoniae card is a coloured chromatographic immunoassay test for the qualitative detection of *Streptococcus* of the subspecies *pneumoniae* in urine samples. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of pneumonia in infected persons.

## WHAT IS PNEUMOCOCCAL?

Pneumococcus (*Streptococcus pneumoniae*) is a Gram-positive bacterium that can cause infections that can lead to life-threatening diseases in children and the elderly. Pneumococcus is responsible for many cases of meningitis, pneumonia, bacteraemia, arthritis, otitis media, sinusitis and conjunctivitis. These infections are often spread through sneezing, coughing or contact with dirty hands. They can be serious and pose bigger health risks in very young children or in immunosuppressed patients. Based on differences in the polysaccharide capsules of the pneumococcal cell wall, *Streptococcus pneumoniae* is classified into over 90 serotypes, which present different antigenic properties and induce different inflammatory responses.

The biotical S. *pneumoniae* card allows early diagnosis of *Streptococcus pneumoniae* infection through detection of a specific soluble antigen present in the urine of patients with *Streptococcus pneumoniae* infection. The result is obtained quickly both in the early phase and in other stages of the disease, conveniently using a urine sample.


## TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

| REFERENCE | NAME                                  | CONTENT                                     |
|-----------|---------------------------------------|---|
| RTB25SNCN | biotical S. pneumoniae Card (+PC +NC) | 25 test with positive and negative control. |

## HOW DOES BIOTICAL *S. PNEUMONIAE* CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In **BLUE** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against *S. pneumoniae*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a test solution (mouse monoclonal antibodies against *S. pneumoniae*) conjugated to red polystyrene latex and a control solution (specific binding protein) conjugated to blue polystyrene latex, forming two conjugate coloured complexes.

### TEST LINE:

**If the sample is positive:** the antigens of the patient's sample will react with the red coloured conjugate complex (anti-*S. pneumoniae* monoclonal antibodies - red latex microspheres) and advance by capillarity through the membrane. The anti-*Campylobacter* antibodies present in the T-line will capture the coloured test complex resulting in the red line.

**If the sample is negative:** if the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-*S. pneumoniae* antibodies present in the T-line will not capture the red coloured antigen-complex (not formed) and the red line will not appear.

### CONTROL LINE (C):

**Regardless of whether the sample is positive or not**, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the membrane of the control line (C). These specific antiprotein antibodies present in the membrane will capture the control conjugate complex and the blue control line will always appear.

## STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

## GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

## MATERIALS

### MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Sample diluent (Reagent).
- Disposable pipettes.
- Testing tubes.
- Positive Control.
- Negative Control.

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.


The kit contains materials for the analysis of a maximum of 25 samples or controls.

## SPECIMEN COLLECTION AND PREPARATION

Urine specimens should be collected in standard containers. The samples can be stored at room temperature (15-30°C) if assayed within 24 hours of collection. Alternatively, specimens may be stored at 2-8°C for up to 14 days or -20°C for longer periods before testing. When necessary, urine specimens should be shipped in leak-proof containers at 2-8°C or frozen. Allow all specimens to equilibrate to room temperature before testing.

## TEST PROCEDURE


Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add **3 drops ①** of **urine sample** in a testing tube. Use a separate testing tube and pipette for each sample or control.
- Add **1 drop ②** of **sample diluent (reagent)** and mix in order to homogenize the sample.
- Remove the test card from its sealed bag just before using it.
- **Dispense 3 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 15 minutes.** The results will not be valid after 10 minutes.

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)

|     | biotical Test                               | Interpretation of the results   |
|-----|---|---|
| 1   | C = <b>BLUE LINE</b><br>T = <b>NO LINE</b>  | <b>NEGATIVE FOR S. PNEUMONIAE</b><br>Only the <b>BLUE</b> line appears across the results window in the site marked with the letter C (control line).<br>No presence of <i>Streptococcus pneumoniae</i> is detected in the sample.  |
| 2   | C = <b>BLUE LINE</b><br>T = <b>RED LINE</b> | <b>POSITIVE FOR S. PNEUMONIAE</b><br>In addition to the <b>BLUE</b> control line, a <b>RED</b> line also appears in the site marked with the letter T (test line) in the results window.<br><i>Streptococcus pneumoniae</i> presence is detected in the sample.   |
| 3-4 | Any other result                            | <b>INVALID</b><br>Total absence of the control line ( <b>BLUE</b> ) regardless the appearance or not of the test line ( <b>RED</b> ) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor. |

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A blue line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits reference **RTB25SNCN**).


## POSITIVE CONTROL: INTENDED USE & COMPOSITION

The biotal *S. pneumoniae* Positive Control is an external quality control for biotal *S. pneumoniae* card. Use this control swab to check that the sample reagents and the test card are working properly.

biotal *S. pneumoniae* Positive Control is manufactured with inactivated *Streptococcus pneumoniae* antigens dried in swab containing a preservative.

## POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add **13 drops of Reagent (sample diluent) slowly** ① in a testing tube.
- Immediately put the positive control swab into the tube ②. Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds**. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 15 minutes.** The results will not be valid after 10 minutes.

## POSITIVE CONTROL – INTERPRETATION OF RESULTS

|   | biotal Test   | Interpretation of the results   |
|---|---|---|
| 1 |  <p>C = BLUE LINE<br/>T = RED LINE</p> | <p><b>RESULT: POSITIVE CONTROL</b></p> <p>A BLUE line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>   |
| 2 |                                        | <p><b>INVALID RESULT</b></p> <p>Total absence of the control line (BLUE), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p> |


## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotal Negative Control is an external quality control for all immunochromatographic tests manufactured by biotal. Use this control swab to check that the extraction reagents and the test are working properly.

biotal Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

## NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 13 drops of Reagent (sample diluent) slowly ① in a testing tube.
- Immediately put the negative control swab into the tube ②. Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds**. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 15 minutes**. The results will not be valid after 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

|   | biotical Test   | Interpretation of the results   |
|---|---|---|
| 1 |  <p>C = BLUE LINE<br/>T = NO LINE</p> | <p><b>RESULT: CONTROL NEGATIVE</b></p> <p>Only the control line (C) appears across the results window during the test performance.</p> <p>NOTE: Depend on used test card, the control line is blue, red or blue. Consult instructions for use of the product.</p>   |
| 2 |                                      | <p><b>RESULT: INVALID FOR NEGATIVE CONTROL</b></p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p> |

## LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- The test has been validated with human urine samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Streptococcus pneumoniae* in the sample. To confirm the presence of infection by *Streptococcus pneumoniae* it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other diagnostic techniques (culture and serological).
- False positives may occur in samples pretreated with boric acid, so the use of the test in this class of samples is not recommended.
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. Negative test results may not be conclusive and may require the use of other laboratory techniques (culture and serological). A negative result does not indicate absence of infection by *Streptococcus pneumoniae*. It is possible to appear false negatives by not reaching the detection limit of the test.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity (detection limit)

Detection limit value of biotical *S. pneumoniae* Card is: 0.25 ng/mL of CWPS.

### Clinical sensitivity and specificity

An evaluation was performed, with urine samples, comparing the results obtained by the immunochromatographic test (biotical *S. pneumoniae* Card, biotical) and another commercial rapid test (BinaxNOW® *Streptococcus pneumoniae* Antigen Card, Alere).

The results were as follows:

|                                    |       | Binax now® <i>Streptococcus pneumoniae</i> Antigen card, Alere (immunochromatographic) |    |       |
|------------------------------------|-------|--|----|-------|
|                                    |       | +  | -  | Total |
| Biotical <i>S. pneumoniae</i> Card | +     | 33   | 1  | 34    |
|                                    | -     | 4  | 73 | 77    |
|                                    | Total | 37   | 74 | 111   |

| <i>Biotical S. pneumoniae</i> Card vs<br><i>Binax now® S. pneumoniae</i> Antigen card, Alere | MEAN VALUE | 95% CONFIDENCE INTERVAL |
|--|------------|-------------------------|
| SENSITIVITY  | 89,2%      | 74,6% - 97,0%           |
| SPECIFICITY  | 98,7%      | 92,7% - 100,0%          |
| PPV  | 97,1%      | 82,4% - 99,6%           |
| NPV  | 94,8%      | 87,9% - 97,9%           |
| ACCURACY   | 95,5%      | 89,8% - 98,5%           |

The results showed a high sensitivity and specificity to detect *Streptococcus pneumoniae* using biotical *S. pneumoniae* Card.

## CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical *S. pneumoniae* Card; no cross reactivity against pathogens present in urine:

*Legionella pneumophila*

## SYMBOLS FOR IVD COMPONENTS & REAGENTS

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