# biotical Legionella Card

Rapid Test to detect  
Legionella in card  
format

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25 UN.

IU RTB25LPCN Rev. 03 / 2024-02

ISO 13485 CERTIFIED  
MEDICAL DEVICE  
QUALITY MANAGEMENT SYSTEM

# biotical Legionella Card

ENGLISH

## WHAT IS BIOTICAL LEGIONELLA CARD?

biotical Legionella Card is a rapid immunochromatography test for the qualitative detection of Legionella pneumophila serogroup 1 in urine samples. Thanks to its easy handling and high sensitivity, the healthcare professional can make a presumptive diagnosis of Legionellosis in infected humans.

## WHAT IS LEGIONELLOSIS?

Legionellosis, or legionnaire’s disease, is a severe form of pneumonia caused by Legionella pneumophila that is usually associated with 10-15% mortality in healthy adults and 40-80% in immunosuppressed patients. The symptoms are like those of a flu, in addition to a dry cough and rapid progress towards pneumonia. If left untreated, legionnaire’s disease usually worsens in the first week, so early diagnosis is key.

The incubation period usually lasts between 2-10 days. The most frequent form of Legionella transmission is the inhalation of contaminated aerosols. Aerosol sources that have been linked to Legionella transmission include cooling towers for air conditioning, hot and cold-water systems, humidifiers, and whirlpool facilities.

biotical Legionella card allows early diagnosis of Legionella pneumophila serogroup 1 infections by detecting a specific soluble antigen present in the urine of patients affected by legionellosis. The test is rapid: the results are obtained in 15 minutes, and a urine sample is used, which is an advantage for the sampling, transport, and subsequent detection of early and advanced stages of the disease.

## TABLE OF REFERENCES


In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25LPCN	biotical Legionella Card (+PC +NC)	25 tests with positive and negative control.

## HOW DOES BIOTICAL LEGIONELLA CARD WORK?

---

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In **BLUE** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against *Legionella pneumophila*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against *Legionella*) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to blue polystyrene latex.

### TEST LINES:

**If the sample is positive:** The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-*Legionella* monoclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti- *Legionella* present in the T line will capture the complex giving rise in a red line.

**If the sample is negative:** If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red-coloured complex. The anti-*Legionella* antibodies present in the T line will not capture the red-coloured antigen-complex (not formed) and the red lines will not appear.

### CONTROL LINE (C):

**Regardless of whether the sample is positive or not**, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the blue control line will always appear.

## STORAGE AND STABILITY

---

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

## GENERAL PRECAUTIONS

---

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- Used tests should be managed as sanitary waste (sanitary waste container).

- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

## MATERIALS

### MATERIALS PROVIDED

- Cassette or Test Card.
- Plastic pipettes
- Testing tubes
- Reagent (sample and controls diluent)
- Instructions for use.
- Positive Control.
- Negative Control.

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.


The kit contains materials for the analysis of a maximum of 25 samples or controls.

## SPECIMEN COLLECTION AND PREPARATION

Urine specimens should be collected in standard containers. The samples can be stored at room temperature (15-30 °C) if assayed within 24 hours of collection. Alternatively, specimens may be stored at 2-8 °C for up to 14 days or at -10 °C to -20 °C for longer periods before testing. When necessary, urine specimens should be shipped in leak-proof containers at 2-8 °C or frozen. Allow all specimens to equilibrate to room temperature before testing.

## TEST PROCEDURE


Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Using a separate testing tube or vial for each sample. Add **3 drops of urine sample** ①.
- Add **1 drop of Reagent** into the testing tube or vial and mix ②. Homogenize the sample.
- Remove biotical card test from its sealed bag just before using it.
- Use a **separate pipette and device for each sample** or control. Dispense **3 drops** from the testing tube in the circular window marked with the letter S ③.
- **Read the results at 15 minutes.** Do not read the test result later than 15 minutes.

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = <b>BLUE LINE</b> T = <b>NO LINE</b>	<b>NEGATIVE</b> Only the <b>BLUE</b> line appears across the results window in the site marked with the letter C (control line). No presence of <i>Legionella pneumophila</i> is detected in the sample.
2	C = <b>BLUE LINE</b> T = <b>RED LINE</b>	<b>POSITIVE</b> In addition to the <b>BLUE</b> control line, a <b>RED</b> line also appears in the site marked with the letter T (test line) in the results window. <i>Legionella pneumophila</i> presence is detected in the sample.
3-4	Any other result	<b>INVALID</b> Total absence of the control line ( <b>BLUE</b> ) regardless the appearance or not of the test line ( <b>RED</b> ) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A blue line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available.


## POSITIVE CONTROLS: INTENDED USE & COMPOSITION

biotical Legionella Positive Control is an external quality control for biotical Legionella card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical Legionella Positive Control: inactivated *Legionella* antigens dried in swab containing a preservative.

## POSITIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add **13 drops of Reagent (sample diluent) slowly** ① in a testing tube.
- Immediately put the positive control swab into the tube ②. Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds**. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 15 minutes**. The results will not be valid after 10 minutes.

## POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = BLUE LINE T = RED LINE</p>	<p><b>RESULT: CONTROL POSITIVE</b></p> <p>A BLUE line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p><b>INVALID RESULT</b></p> <p>Total absence of the control line (BLUE), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

## NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 13 drops of Reagent (sample diluent) slowly ① in a testing tube.
- Immediately put the negative control swab into the tube ②. Mix the solution by **rotating the swab forcefully against the side of the tube at least 30 seconds**. Best results are obtained when the specimen is vigorously extracted in the solution. **Extract as much liquid as possible** from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
- Remove the test card from its sealed bag just before using it.
- **Dispense 4 drops** from the testing tube, into the circular window marked with the letter S ③. Use a separate pipette and test for each sample or control.
- **Read the results at 15 minutes**. The results will not be valid after 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = <b>BLUE LINE</b> T = <b>NO LINE</b></p>	<p><b>RESULT: CONTROL NEGATIVE</b></p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p><b>RESULT: INVALID FOR NEGATIVE CONTROL</b></p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

## LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- biotical Legionella Card has been validated using human urine samples only. Other samples (e.g. plasma, serum or other body fluids) that may contain Legionella antigen have not been evaluated. The test cannot be used on environmental samples (i.e. potable water).
- Positive results determine the presence of *L. pneumophila* serogroup 1 in urine samples; nevertheless, a positive result should be followed up with additional laboratory techniques (culture or serology) to confirm the results. Therefore, culture results, serology and antigen detection methods should be used in conjunction with clinical findings to make an accurate diagnosis.
- This test will not detect infections caused by other *L. pneumophila* serogroups and by other *Legionella* species. A negative antigen result does not exclude infection with *L. pneumophila* serogroup 1. Culture is recommended for suspected pneumonia to detect causative agents other than *L. pneumophila* serogroup 1 and to recover *L. pneumophila* serogroup 1 when antigen is not detected in urine.
- Excretion of *Legionella* antigen in urine may vary depending on the individual patient. Antigen excretion may begin as early as 3 days after onset of symptoms and persist for up to 1 year afterwards. A

positive biotical Legionella card result can occur due to current or past infection and therefore is not definitive for infection without other supporting evidence.

- It has been proven that although the test is compatible with boric acid, its presence increases the reactivity and can lead to false positives, so it is not recommended as a preservative.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity (detection limit):

Detection limit range is: 12.5ng/mL (pool of several serovars of Legionella pneumophila).

### Clinical sensitivity and specificity.

An evaluation was performed, with urine samples, comparing the results obtained by an immunochromatographic test (biotical Legionella Card, biotical) and other commercial rapid test (Binax NOW® Legionella Urinary Antigen Card, Alere).

The results were as follows:

		Ic test: Binax noW® legionella Urinary Antigen card		
		+	-	Total
biotical Legionella Card	+	31	1	32
	-	0	117	117
	Total	31	118	149

<i>biotical Legionella vs Binax NOW® Legionella Urinary Antigen Card</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99.9%	88.8 - 100%
SPECIFICITY	99.2%	95.4 - 100%
PPV	96.9%	81.5 - 99.5%
NPV	>99.9%	-
ACCURACY	99.3%	96.3 - 100%

The results showed a high sensitivity and specificity to detect *L. pneumophila* serogrupo 1. using biotical Legionella Card.

## CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Legionella Card; no cross reactivity against other pathogens occasionally present in urine:

*Streptococcus Pneumoniae*

## SYMBOLS FOR IVD COMPONENTS & REAGENTS

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