ENGLISH - ESPAÑOL

Ref: RTB25E
RTB25EC
RTB25ECN

biotical

E. Coli O157:H7 Card

Rapid Test to detect
E. Coli O157:H7 in card
format

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biotical health
ready to test

IU RTB25E-RTB25EC-RTB25ECN Rev. 02 / 2024-02

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical

E. coli O157:H7 Card

ENGLISH

WHAT IS BIOTICAL E. COLI O157:H7 CARD?

biotical E. coli O157:H7 Card is a rapid immunochromatography test for the qualitative detection of *Escherichia coli O157:H7* (*E. coli O157:H7*) in stool samples and suspicious colonies of *E. coli O157:H7* in stool culture. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of infection caused by *E. coli O157:H7* and could be used for the identification of suspected isolates of *E. coli O157:H7* from a selective medium (stool culture).

WHAT IS E. COLI O157: H7?

Escherichia coli (*E. coli*) is a bacterium present in the distal intestine of homeothermic organisms. Most strains of *E. coli* are harmless, but some, such as O157:H7, can cause serious food poisoning.

E. coli O157:H7 belongs to the Enterohemorrhagic group *Escherichia coli*, EHEC (also called Shiga Toxin producing *Escherichia coli*, STEC) for its ability to produce the “Shiga” toxin. This toxin can cause serious poisoning through contaminated products such as raw or undercooked meat, unpasteurized milk, water and fresh products contaminated by faecal matter, and by contagion between people through the faecal-oral route.

Some of the symptoms of the “Shiga” toxin could be: diarrhoea (which may vary from mild and liquid to severe and bloody), abdominal cramps, nausea, vomiting and fever. The incubation period varies between three and eight days, with an average of three to four days. Most patients recover within ten days, but in a small percentage of cases (especially young children and the elderly) the infection can lead to a life-threatening disease, such as Haemolytic Uremic Syndrome (SHU). HUS is characterized by acute renal failure, haemolytic anaemia and thrombocytopenia (platelet deficiency).

To confirm the diagnosis and identify specific toxins, such as those produced by *E. coli O157:H7*, blood, faeces, sometimes urine or other infected material are usually taken and sent to the laboratory for bacterial culture. This procedure may take several days. However, biotical E. coli O157:H7 Card, is a qualitative immunochromatographic test for the detection of *Escherichia coli O157:H7* in stool samples and suspicious colonies of *E. coli O157:H7* in stool culture, fast and easy to handle with which results are obtained in a very short time.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25E	biotical E. Coli O157:H7 Card	25 Tests without positive or negative controls.
RTB25EC	biotical E. Coli O157:H7 Card (+PC)	25 Tests with positive control, without negative control.
RTB25ECN	biotical E. Coli O157:H7 Card (+PC +NC)	25 Tests with positive and negative control.

HOW DOES BIOTICAL E. COLI O157:H7 CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against E. coli O157:H7.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a test solution (antibodies against E. coli O157:H7) conjugated to red polystyrene latex and a control solution (specific binding protein) conjugated to green polystyrene latex, forming two conjugate-coloured complexes.

TEST LINES (T):

If the sample is positive: The antigens of the patient's sample will react with the red-coloured conjugate complex (anti E. coli O157:H7 antibodies-red latex microspheres) and advance by capillarity through the membrane. The anti E. coli O157:H7 antibodies present in the membrane of the test line (T), will capture the coloured test complex resulting in the red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti E. coli O157:H7 antibodies present in the membrane of the test line (T), will not capture the red coloured antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25EC & RTB25ECN references).
- Negative Control (only RTB25ECN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Agar MacConkey-sorbitol (SMACA)

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION


Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- Take out the cap of the collection tube **1** and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **2**, to collect faecal sample (approx. 125mg) and add it to the collection tube. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **3**.

1 Collection tube

2 Stool sample

3 Shake vigorously

Suspected E. coli o157:H7 colonies in stool culture:

The Sorbitol MacConkey Agar (SMACA) is the method of choice for the isolation of *E. coli* O157:H7 (aerobic atmosphere, 24 hours/37 °C). After 24 hours of incubation in Sorbitol MacConkey the typical *E. coli* O157:H7 colonies will be small and colourless.

- Examine Sorbitol MacConkey plates after 1 day incubation. Select *E. coli* O157:H7 typical colonies. Take out the cap of the collection tube 1. Use the stick or an inoculating needle to pick up 3 or 4 suspected *E. coli* O157:H7 colonies 2 and add them to the collection tube.
- Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion 3.

E.Coli O157:H7 colonies in SMACA

1 Colection tube


2 E.Coli O157:H7 colonies in plate


3 shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

4 Cut the end of the cap.


5 Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 4.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 5. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID RESULT

INVALID RESULT

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>E. coli O157:H7</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>E. coli O157:H7</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25EC** & **RTB25ECN**).

POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotal E. coli O157:H7 Positive Control is an external quality control for biotal E. COLI O157:H7 card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotal E. coli O157:H7 Positive Control: inactivated *E. coli* O157:H7 antigens dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

Results window

GREEN C line (Control)

RED T line (Test E. coli O157:H7)

Sample window (S)

- 1 Add 15 drops of diluent and introduce the positive control swab into the testing tube.
- 2 Rotating 30 seconds and extract the liquid.
- 3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube 1.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

NEGATIVE CONTROL: INTENDED USE & COMPOSITION


The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

Results window

GREEN C line (Control)

RED T line (Test E. coli O157:H7)


Sample window (S)

1 Add 15 drops of diluent and introduce the negative control swab into the testing tube.

2 Rotating 30 seconds and extract the liquid.

3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube ①.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test must be performed within 2 hours after opening the sealed pouch.
- Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
- biotical E. coli O157:H7 Card will only indicate the presence of *E. coli* in the specimen (qualitative detection) and should be used for the detection of *E. coli* antigens in human faeces specimens only. Neither the quantitative value nor the rate of increase in *E. coli* antigens concentration can be determined by this test.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- Samples should be collected as soon as possible once the symptoms appear. After one week of infection the number of bacteria decrease, making the sample less reactive.
- Some stool samples can decrease the intensity of the control line.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods are recommended. A negative result does not at any time preclude the possibility of *E. coli* infection as it may be due to a very low antigen concentration.
- This test provides a presumptive diagnosis of *E. coli* infections. All results must be interpreted together with other clinical information and laboratory findings available to the physician.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value of biotical E. coli O157:H7 Card is 1.87×10^4 CFU/mL.

Clinical sensitivity and specificity

biotical E. coli O157:H7 Card was evaluated to determine sensibility in selective enrichment culture and samples, specificity with producers' organisms of Shiga toxins, non-Shiga toxins producers and other Enterobacteriaceae species in the Reference Laboratory for Escherichia coli (LREC), Universidad Santiago de Compostela, Lugo (SPAIN).

14 STEC strains (O157:H7 antigen), 4 Non STEC strains (O157), 9 STEC strains (non O157), 4 other Enterobacteriaceae spp.

The results were as follows:

		CULTURE			
		+ STEC (O157 antigen)	- Non STEC (O157)	- Non O157	Total
Biological E. coli O157:H7 Card	+	14	4	2	20
	-	0	0	11	11
	Total	14	4	13	31

	SENSITIVITY	SPECIFICITY	PPV	NPV
<i>Biological E. coli O157:H7 Card vs culture</i>	>99%	85%	70%	>99%

The results showed a high sensitivity and specificity to detect *Escherichia coli O157:H7* using biotical E. coli O157:H7 Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical E. coli O157:H7 Card; no cross reactivity against gastrointestinal organism occasionally present in faeces.

Campylobacter coli
Campylobacter jejuni
Citobacter freundii
Clostridium difficile
E. Coli O22:H8
E. Coli O91:H-
E. Coli O103:H2
E. Coli O111:H21
E. Coli O145:H-

E. Coli O171:H2
E. Coli O174:H8
Klebsiella pneumoniae
Helicobacter pylori
Listeria monocytogenes
Morganella morganii
Proteus mirabilis
Salmonella enteritidis
Salmonella paratyphi

Salmonella typhi
Salmonella typhimurium
Shigella boydii
Shigella dysenteriae
Shigella flexneri
Shigella sonnei
Staphylococcus aureus
Yersinia enterocolitica

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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