ENGLISH - ESPAÑOL

Ref: RTB25SA
RTB25SAC
RTB25SACN

biotical Salmonella Card

Rapid Test to detect
Salmonella in card
format

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biotical health
ready to test

IU RTB25SA-RTB25SAC-RTB25SACN Rev. 04 / 2024-04

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical Salmonella Card

ENGLISH

WHAT IS BIOTICAL SALMONELLA CARD?

biotical Salmonella Card is a rapid immunochromatography test for the qualitative detection of *Salmonella* antigens in stool samples and suspicious colonies of *Salmonella* in stool culture. Thanks to its easy handling and high sensitivity, the healthcare professional is able to make a presumptive diagnosis of *Salmonella* infection (salmonellosis) and it could be used to identify suspected *Salmonella* isolates from a selective medium (stool culture).

WHAT IS SALMONELLOSIS?

Salmonella infection, or salmonellosis, is a foodborne disease caused by a bacterium of the *Salmonella* gender. Most infections are caused due to undercooked contaminated foods (usually veal, poultry, eggs or milk), by cross contamination (other foods that have been in contact with contaminated food due to contaminated hands and contaminated kitchen utensils), or by person-to-person contact.

Salmonellosis usually causes the following symptoms: nausea, vomiting, abdominal pain, diarrhoea (sometimes bloody diarrhoea), fever and headache, which appear between 12-72 hours (but sometimes they can take up to 7 days) after infection. These symptoms, together with the dehydration caused by the illness, are aggravated in higher risk people (children and infants, immunosuppressed, elderly, etc.).

Considering that there are many diseases that can cause the same symptomatic picture as salmonellosis, it is important to make an accurate diagnosis.

Biotal Salmonella Card is a qualitative immunochromatographic test for the detection of *Salmonella* in stool samples and suspicious colonies of *Salmonella* in stool culture.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25SA	biotical Salmonella Card	25 tests without positive or negative controls.
RTB25SAC	biotical Salmonella Card (+PC)	25 tests with positive control, without negative control.
RTB25SACN	biotical Salmonella Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL SALMONELLA CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In **RED** line (T): with mouse monoclonal antibodies against *Salmonella*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse polyclonal antibodies against *Salmonella*) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex.

TEST LINES:

If the sample is positive: The antigens of the patient's stool sample will react with the red coloured conjugate complexes (anti-*Salmonella* polyclonal antibodies-latex red microspheres) and it will advance by capillary through the membrane. As the sample migrates, so do the conjugate complexes. The antibodies anti-*Salmonella* present in the T line will capture the complex giving rise in a red line.

If the sample is negative: If the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with the red coloured complex. The anti-*Salmonella* antibodies present in the T line will not capture the red coloured antigen-complex (not formed) and the red lines will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.

- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- Pale yellow or orange bands may appear in the test strip results window, at the positions where the coloured test and control lines should appear, before the test is performed. These bands are completely normal and correspond to the capture antibodies of the test. The existence of these bands does not alter the state of conservation and functionality of the kit.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25SAC & RTB25SACN references).
- Negative Control (only RTB25SACN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Hektoen Enteric Agar (HEA) / Salmonella Shigella Agar (SS).

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION

Stool samples.

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

- Take out the cap of the collection tube **1** and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **2**, making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample and add it to the collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **3**.

1 Collection tube

2 Stool sample


3 Shake vigorously

Suspected *Salmonella* colonies in stool culture:

Method 1: Hektoen Enteric Agar (HEA):

The Hektoen Enteric Agar (HEA) is the method of choice for the isolation of *Salmonella* (aerobic atmosphere, 24 hours/37 °C). After 24 hours of incubation in Hektoen the green *Salmonella* colonies.

- Examine Hektoen plates after 1 day incubation. Select *Salmonella* typical colonies. Take out the cap of the collection tube ①. Use the stick or an inoculating needle to pick up 3 or 4 suspected *Salmonella* colonies ② and add them to the collection tube.
- Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion ③.


Salmonella colonies in HEA

If the content of *Salmonella* in the faeces is not large enough, false negatives can occur because the antigen is not detected. It is recommended to carry out the enrichment in the Rappaport-Vassiliadis medium (aerobic atmosphere, 24h/37°C). After 24h of incubation, it is spread on Hektoen and incubated (aerobic atmosphere, 24h/37°C).

Method 2: Salmonella Shigella Agar (SS):

Salmonella Shigella (SS) Agar is another method for the isolation of *Salmonella* (aerobic atmosphere, 24-48 hours/35-37 °C). After the incubation period in SS the typical *Salmonella* colonies will be transparent. SS agar is a highly inhibitory medium that sometimes prevents the growth of certain *Salmonella* species, so it is recommended to use this medium accompanied by another less inhibitory medium such as HEA.

- Examine SS plates after incubation. Select *Salmonella* typical colonies. Take out the cap of the collection tube ①. Use the stick or an inoculating needle to pick up 3 or 4 suspected *Salmonella* colonies ② and add them to the collection tube.
- Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion ③.

Salmonella colonies in SS

① Collection tube


② *Salmonella* colonies in plate

③ Shake vigorously

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

4 Cut the end of the cap.


5 Add 3 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 4.
- **Dispense 3 drops** from the testing tube in the circular window marked with the letter S 5. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

NEGATIVE

POSITIVE

INVALID RESULT

INVALID RESULT

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>Salmonella</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>Salmonella</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL & NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25SAC** & **RTB25SACN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical Salmonella Positive Control is an external quality control for biotical SALMONELLA card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical Salmonella Positive Control: inactivated *Salmonella* antigens dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube **1**.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds **2**. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S **3**.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>


NEGATIVE CONTROL: INTENDED USE & COMPOSITION

biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- **Add 15 drops of diluent** and immediately introduce the negative control swab into the testing tube 1.
- Mix the solution by **rotating the negative control swab** forcefully against the side of the testing tube at least **30 seconds** 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 3.
- **Read the result at 10 minutes.** do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Salmonella* in the sample, without determining its cause or species. To confirm the presence of infection it is necessary for a specialist to evaluate the results together with the other clinical evidence of the patient. In order to issue a valid clinical diagnosis, it is necessary to confirm the results of the tests using other laboratory techniques (serological methods, biochemicals or PCR confirmation).
- False positives may occur in the test with bloody samples or mucous stools, so the use of the test in this class of samples is not recommended.
- If a negative test result is obtained and symptoms persist, it is recommended to repeat the test after a few days. It is possible to appear false negatives by not reaching the detection limit of the test.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit range for the different species are:

- *Salmonella enteritidis*: $5 \cdot 10^7 - 3.12 \cdot 10^6$ UFC/mL of inactivated culture.
- *Salmonella typhimurium*: $10 \cdot 10^6 - 6.25 \cdot 10^5$ UFC/mL of inactivated culture.
- *Salmonella typhi*: $10 \cdot 10^7 - 6.25 \cdot 10^6$ UFC/mL of inactivated culture.

Clinical sensitivity and specificity:

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical Salmonella Card, biotical) and qPCR technique (VIASURE Salmonella, Campylobacter & Y. enterocolitica Real Time PCR Detection kit).

The results were the following:

		qPCR test: Viasure Salmonella, Campylobacter & Y. enterocolitica PCR Real Time Detection kit		
		+	-	Total
Biotical Salmonella Card	+	47	2	49
	-	15	103	118
	Total	62	105	167

<i>biotical Salmonella Card vs VIASURE Salmonella, Campylobacter & Y. enterocolitica Real Time PCR Detection kit</i>	MEAN VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	75.8%	63.3 - 85.8%
SPECIFICITY	98.1%	93.5 - 99.8%
PPV	95.9%	86.0 - 99.5%
NPV	87.3%	79.9 - 92.7%
ACCURACY	89.8%	84.2 - 94.0%

The results showed a high sensitivity and specificity to detect *Salmonella* using biotical Salmonella Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Salmonella Card; no cross reactivity against gastrointestinal organism occasionally present in faeces.

<i>Adenovirus</i>	<i>Escherichia coli</i> O:026	<i>Norovirus</i> GII
<i>Astrovirus</i>	<i>Escherichia coli</i> O157	<i>Rotavirus</i>
<i>Calprotectin</i>	<i>Entamoeba dispar</i>	<i>Shigella boydii</i>
<i>Campylobacter coli</i>	<i>Entamoeba histolytica</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	<i>Giardia lamblia</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i> ag. GDH	<i>Helicobacter pylori</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i> Toxin A	Hemoglobin	<i>Streptococcus pneumoniae</i>
<i>Clostridium difficile</i> Toxin B	Lactoferrin	<i>Streptococcus pyogenes</i>
<i>Clostridium perfringens</i>	<i>Legionella</i>	Transferrin
<i>Cryptosporidium</i>	<i>Listeria monocytogenes</i>	<i>Yersinia enterocolitica</i> O:3
<i>Escherichia coli</i> O:111	<i>Norovirus</i> GI	<i>Yersinia enterocolitica</i> O:9

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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REF

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REFERENCES / BIBLIOGRAFÍA /

1. TABUBAKAR, I. "A systematic review of the clinical, public health and cost-effectiveness of rapid diagnostic tests for the detection and identification of bacterial intestinal pathogens in faeces and food". Health Technology Assessment 2007; Vol. 11: No. 36.
2. OCHOA REPÁRAZ, J. "Protective ability of subcellular extracts from Salmonella Enteritidis and from a rough isogenic mutant against salmonellosis in mice". Vaccine 23 (2005), 1495-1501.

Sierra de Guadarrama, 1. Ed. Dublín.
28830 San Fernando de Henares. Madrid
Tel.: + 34 91 677 43 08

info@biotical.es
www.biotical.es