ENGLISH - ESPAÑOL

Ref: RTB25ST
RTB25STC
RTB25STCN

biotical Salmonella typhi Card

Rapid Test to detect
Salmonella typhi in card
format

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biotical health
ready to test

IU RTB25ST-RTB25STC-RTB25STCN Rev. 01 / 2024-02

ISO 13485 CERTIFIED
MEDICAL DEVICE
QUALITY MANAGEMENT SYSTEM

biotical *Salmonella typhi* Card

ENGLISH

WHAT IS BIOTICAL SALMONELLA TYPHI CARD?

biotical *Salmonella typhi* Card is a coloured chromatographic immunoassay test for the qualitative detection of *Salmonella typhi* in stool samples and *Salmonella typhi* suspected colonies in stool culture. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of *Salmonella typhi* infection and it could be used to identify suspected *Salmonella typhi* isolates from selective media (stool culture).

ENTERIC FEVER

Enteric, typhoid, or paratyphoid fever is a major health problem in many endemic countries.

Symptoms of enteric fever include fever, weakness, lethargy, muscle pain, headache, loss of appetite, diarrhoea or constipation, pink spots on the cheeks and, in the most severe cases, enlarged liver and spleen, altered mental status, and meningitis (fever, stiff neck, and seizures). The difference between enteric fever caused by *Salmonella typhi* and enteric fever caused by *Salmonella paratyphi* is that the former usually appears between one and three weeks after exposure and varies in severity.

The high mortality and morbidity caused by the different salmonellas is also due to their resistance to multiple drugs, with reduced susceptibility to fluoroquinolones, which makes their rapid detection crucial. With the biotical *Salmonella typhi* Card, the health professional will be able to make a rapid presumptive diagnosis of enteric fever, and/or use the test to identify suspected *Salmonella typhi* isolates from a stool culture.


TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25ST	biotical <i>Salmonella typhi</i> Card	25 tests without positive or negative controls.
RTB25STC	biotical <i>Salmonella typhi</i> Card (+PC)	25 tests with positive control, without negative control.
RTB25STCN	biotical <i>Salmonella typhi</i> Card (+PC +NC)	25 tests with positive and negative control.

HOW DOES BIOTICAL SALMONELLA TYPHI CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously pre-coated with:

- In the **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In the **RED** line (T): with mouse monoclonal antibodies against *Salmonella typhi*.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against *Salmonella typhi*) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex.

TEST LINE (T):

If the sample is positive: the antigens of the patient's sample will react with the RED coloured conjugate complex (anti- *Salmonella typhi* monoclonal antibodies red latex microspheres for the line T) and advance by capillarity through the membrane. The anti-*Salmonella typhi* antibodies present in the T line will capture the coloured conjugate complex resulting in a RED line.

If the sample is negative: if the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with red coloured complex. The anti-*Salmonella typhi* antibodies present in the T line will capture the coloured antigen-complex (not formed) and the red line will not appear.

CONTROL LINE (C):

Regardless of whether the sample is positive or not, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the GREEN control line will always appear.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

MATERIALS

MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25STC & RTB25STCN references).
- Negative Control (only RTB25STCN reference).

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- Hektoen Enteric Agar (HEA) / Salmonella Shigella Agar (SS)

The kit contains materials for the analysis of a maximum of 25 samples or controls.

SPECIMEN COLLECTION AND PREPARATION


Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):


- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample ①, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion ②. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

Collection tube

① Stool sample

② Shake vigorously


Suspected *Salmonella typhi* colonies in stool culture:

1 Collection tube

2 *Salmonella typhi* colonies in plate


3 Shake vigorously

Method 1: Hektoen Enteric Agar (HEA):

The Hektoen Enteric Agar (HEA) is the method of choice for the isolation of *Salmonella typhi* (aerobic atmosphere, 24 hours/37 °C). After 24 hours of incubation in Hektoen the green *Salmonella typhi* colonies.

- Examine Hektoen plates after 1 day incubation. Select *Salmonella typhi* typical colonies. Take out the cap of the collection tube 1. Use the stick or an inoculating needle to pick up 3 or 4 suspected *Salmonella typhi* colonies 2 and add them to the collection tube.
- Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion 3.


Salmonella typhi colonies in HEA

If the content of *Salmonella* in the faeces is not large enough, false negatives can occur because the antigen is not detected. It is recommended to carry out the enrichment in the Rappaport-Vassiliadis medium (aerobic atmosphere, 24h/37°C). After 24h of incubation, it is spread on Hektoen and incubated (aerobic atmosphere, 24h/37°C).

Method 2: Salmonella Shigella Agar (SS):

Salmonella Shigella (SS) Agar is another method for the isolation of *Salmonella typhi* (aerobic atmosphere, 24-48 hours/35-37 °C). After the incubation period in SS the typical *Salmonella typhi* colonies will be transparent. SS agar is a highly inhibitory medium that sometimes prevents the growth of certain *Salmonella* species, so it is recommended to use this medium accompanied by another less inhibitory medium such as HEA.


- Examine SS plates after incubation. Select *Salmonella typhi* typical colonies. Take out the cap of the collection tube 1. Use the stick or an inoculating needle to pick up 3 or 4 suspected *Salmonella typhi* colonies 2 and add them to the collection tube.
- Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion 3.

Salmonella typhi colonies in SS

TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.


- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** 4.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S 5. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	NEGATIVE Only the GREEN line appears across the results window in the site marked with the letter C (control line). No presence of <i>Salmonella typhi</i> is detected in the sample.
2	C = GREEN LINE T = RED LINE	POSITIVE In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window. <i>Salmonella typhi</i> presence is detected in the sample.
3-4	Any other result	INVALID Total absence of the control line (GREEN) regardless the appearance or not of the test line (RED) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references **RTB25STC** & **RTB25STCN**).


POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotical *Salmonella typhi* Positive control is an external quality control for biotical *Salmonella typhi* Card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotical *Salmonella typhi* Positive control is manufactured with *Salmonella typhi* inactivated extract dried in swab containing a preservative.

POSITIVE CONTROLS: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube ①.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

POSITIVE CONTROL – INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = RED LINE</p>	<p>RESULT: CONTROL POSITIVE</p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p>INVALID RESULT</p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

NEGATIVE CONTROL: INTENDED USE & COMPOSITION

biological Negative Control is an external quality control for all immunochromatographic tests manufactured by biological. Use this control swab to check that the extraction reagents and the test are working properly.

biological Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

NEGATIVE CONTROL: TESTING PROCEDURE

Allow tests and controls to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1 Add 15 drops of diluent and introduce the negative control swab into the testing tube.

2 Rotating 30 seconds and extract the liquid.

3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube 1.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes. do not read the test result later than 10 minutes.

NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	<p>C = GREEN LINE T = NO LINE</p>	<p>RESULT: CONTROL NEGATIVE</p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p>RESULT: INVALID FOR NEGATIVE CONTROL</p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- Poorly defined brown lines may appear. These lines are caused by insufficient dilution of the sample and have no diagnostic value. If it occurs, the test must be repeated in a new cassette, increasing the dilution factor of the sample.
- The intensity of the test line can depend on whether the antigen concentration is close to the detection limit value of the test or not.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Salmonella typhi* in the sample, nevertheless it must be confirmed together with other laboratory techniques (biochemical and serological methods or by PCR). The infection must be confirmed by a specialist, after evaluating the clinical tests and all the laboratory results.
- False negative may occur in the test in samples with a low concentration of *Salmonella Typhi* (below the limit of detection value). Repeat test if symptoms persist.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Cut-off value of biotical *Salmonella typhi* Card is 1.25×10^7 UFC/mL.

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotical *Salmonella typhi* Card, biotical) and an ELISA test (*Salmonella typhi* ELISA kit, ICL).

An evaluation was carried out, with culture of *Salmonella typhi*, with an immunochromatographic test (biotical *Salmonella typhi* Card, biotical) and these results were compared with the results obtained by a commercialized immunochromatographic test (Singlepath® *Salmonella*, Merck).

The results were as follows:

		IC Test: Singlepath® <i>Salmonella</i>		
		+	-	Total
IC Test: Biotical <i>Salmonella typhi</i> Card	+	10	0	10
	-	0	25	25
	Total	10	25	35

<i>Biotical Salmonella typhi Card vs Singlepath® Salmonella, Merck</i>	MIDDLE VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99.9%	69.2% - 100%
SPECIFICITY	>99.9%	86.3% - 100%
PPV	>99.9%	-
NPV	>99.9%	-
ACCURACY	>99.9%	90.0% - 100%

The results showed a high sensitivity and specificity to detect *Salmonella typhi* using biotical *Salmonella typhi* Card.

CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical *Salmonella typhi* Card; no cross reactivity against gastrointestinal pathogens occasionally present in faeces.

Campylobacter coli

Campylobacter jejuni

Clostridium difficile

Escherichia coli O157:H7

Helicobacter pylori

Listeria monocytogenes

Shigella boydii

Shigella dysenteriae

Shigella flexneri

Shigella sonnei

Streptococcus aureus

Yersinia enterocolitica

SYMBOLS FOR IVD COMPONENTS & REAGENTS

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