ENGLISH - ESPAÑOL

Ref: RTB25Y9  
RTB25Y9C  
RTB25Y9CN

# biotical Yersinia enterocolitica O:9 Card

Rapid Test to detect  
Yersinia enterocolitica  
O:9 in card format

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ready to test

IU RTB25Y9-RTB25Y9C-RTB25Y9CN Rev. 02 / 2024-01

ISO 13485 CERTIFIED  
MEDICAL DEVICE  
QUALITY MANAGEMENT SYSTEM

# biotical *Yersinia enterocolitica* O:9 Card

ENGLISH

## WHAT IS BIOTICAL *YERSINIA ENTEROCOLITICA* O:9 CARD?

biotical *Yersinia enterocolitica* O:9 Card is a coloured chromatographic immunoassay test for the qualitative detection of *Yersinia enterocolitica* O:9 in stool samples and *Yersinia enterocolitica* O:9 suspected colonies in stool culture. Thanks to its easy handling and high sensitivity, the healthcare professional will be able to make a presumptive diagnosis of *Yersinia enterocolitica* O:9 infection and it could be used to identify suspected *Yersinia enterocolitica* O:9 isolates from selective media (stool culture).

## YERSINIOSIS

*Yersinia* is a genus of bacteria found naturally in the environment and in animals. The *Yersinia* bacterium can be transmitted by lack of hygiene when handling food throughout the food chain: on the farm or livestock farm, in the slaughterhouse, in industry and, of course, at home.

Yersiniosis was, in 2018, the 4th most detected alimentary zoonosis in Europe. It is transmitted by eating contaminated food, especially undercooked pork, unpasteurized dairy, or vegetables irrigated with contaminated water, especially roots (for example, carrots) stored cold for long periods of time.

Numerous biotypes and serotypes of the species *Yersinia enterocolitica* have been described. Eleven of these serotypes have been most frequently associated with infections in humans. Most *Yersinia enterocolitica* strains associated with human yersiniosis belong to bioserotypes 1B/O:8, 2/O:5,27, 2/O:9, 3/O:3, and 4/O:3.

biotical *Yersinia enterocolitica* O:9 Card is a qualitative immunochromatographic test for the detection of *Yersinia Enterocolitica* O:9 in stool samples and suspected colonies of *Yersinia* in stool culture.


## TABLE OF REFERENCES

In the following table will be able to find all the product references related to the acquired test.

REFERENCE	NAME	CONTENT
RTB25Y9	biotical <i>Yersinia enterocolitica</i> O:9 Card	25 tests without positive or negative controls.
RTB25Y9C	biotical <i>Yersinia enterocolitica</i> O:9 Card (+PC)	25 tests with positive control, without negative control.
RTB25Y9CN	biotical <i>Yersinia enterocolitica</i> O:9 Card (+PC +NC)	25 tests with positive and negative control.

## HOW DOES BIOTICAL *YERSINIA ENTEROCOLITICA* O:9 CARD WORK?

The card contains inside a test strip with a nitrocellulose membrane that has been previously precoated with:

- In the **GREEN** line (C): with rabbit polyclonal antibodies against a specific protein.
- In the **RED** line (T): with mouse monoclonal antibodies against *Yersinia enterocolitica* O:9.
- **Sample area (S)**: the pad where the sample is placed has been dispensed with a preparation of reagents (mouse monoclonal antibodies against *Yersinia enterocolitica* O:9) conjugated with latex of red polystyrene and with another preparation for the control line C (specific binding protein) conjugated to green polystyrene latex.

### TEST LINES:

**If the sample is positive:** the antigens of the patient's sample will react with the red coloured conjugate complex (anti-*Yersinia enterocolitica* O:9 monoclonal antibodies red latex microspheres for the line T) and advance by capillarity through the membrane. The anti-*Yersinia enterocolitica* O:9 antibodies present in the T line will capture the coloured conjugate complex resulting in a red line.

**If the sample is negative:** if the sample is negative and therefore does not contain antigens (or the antigen is present in a concentration below the detection limit), there is no reaction with red coloured complex. The anti-*Yersinia enterocolitica* O:9 antibodies present in the T line will capture the coloured antigen-complex (not formed) and the red line will not appear.

### CONTROL LINE (C):

**Regardless of whether the sample is positive or not**, to verify that the test has been performed correctly (sufficient sample volume, appropriate flow, internal control of the reagents, etc.) the mixture will continue moving through the membrane towards the immobilized antibodies against the specific protein located in the control line (C). These specific anti-protein antibodies present in the membrane will capture the control conjugate complex and the green control line will always appear.

## STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

## GENERAL PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The test and samples should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the test that are supplied. Do not use any other commercial kit component if compatibility is not specified.

- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

## MATERIALS

### MATERIALS PROVIDED

- Cassette or Test Card.
- Instructions for use.
- Stool collection tubes with diluent.
- Positive Control (only RTB25Y9C & RTB25Y9CN references).
- Negative Control (only RTB25Y9CN reference).

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container.
- Disposable gloves.
- Timer.
- CIN media (Agar)


The kit contains materials for the analysis of a maximum of 25 samples or controls.

## SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen (-20°C). In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

### Specimen preparation (see illustration):


- Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample **1**, to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
- Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion **2**. Use only fresh samples or fresh frozen samples without preservatives or fixatives.

Suspected *Yersinia enterocolitica* O:9 colonies in stool culture:

The CIN media (faecal culture) is the method of choice for the isolation of *Yersinia* spp. (aerobes, 48 hours/30°C). After incubation in CIN Agar (Cefsulodin-Irgasan-Novobiocin Agar) the typical *Yersinia enterocolitica* colonies appears.

- Examine CIN plates after 1 day incubation. Select small (1-2mm diameter) colonies having deep red centre with sharp border surrounded by clear colourless zone with entire edge. Take out the cap of the collection tube ①. Use the stick or an inoculating needle to pick up 3 or 4 suspected *Yersinia* colonies ② and add them to the collection tube.
- Close the collection tube with the diluent and suspected colonies. Shake the collection tube in order to assure good sample dispersion ③.

*Yersinia* en colonies in Agar CIN

① Collection tube

② *Yersinia* colonies in CIN Agar


③ Shake

Note: After 42 to 48 hours of incubation, the colonies are often completely pink.

## TEST PROCEDURE

Allow tests and samples to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

③ Cut the end of the cap.


④ Add 4 drops in the circular window (S).

- Proceed to **shake the stool collection tube** in order to assure good sample dispersion.
- Remove the test card from its sealed bag just before using it.
- Take the stool collection tube, **cut the end of the cap** ③.
- **Dispense 4 drops** from the testing tube in the circular window marked with the letter S ④. Avoid adding solid particles with the liquid.
- **Read the results at 10 minutes.** The results will not be valid after 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

## INTERPRETATION OF RESULTS

(Please refer to the illustration below)

	biotical Test	Interpretation of the results
1	C = GREEN LINE T = NO LINE	<b>NEGATIVE</b> Only the <b>GREEN</b> line appears across the results window in the site marked with the letter C (control line). No presence of <i>Yersinia enterocolitica</i> O:9 is detected in the sample.
2	C = GREEN LINE T = RED LINE	<b>POSITIVE</b> In addition to the <b>GREEN</b> control line, a <b>RED</b> line also appears in the site marked with the letter T (test line) in the results window. <i>Yersinia enterocolitica</i> O:9 presence is detected in the sample.
3-4	Any other result	<b>INVALID</b> Total absence of the control line ( <b>GREEN</b> ) regardless the appearance or not of the test line ( <b>RED</b> ) causes an invalid result. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the result still persists, discontinue using the test kit and contact your local distributor.

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window as an internal control, which confirms sufficient specimen volume and correct procedural technique.

The following sections (**POSITIVE CONTROL** & **NEGATIVE CONTROL**) only apply if the corresponding positive and negative controls are available (sold separately or with kits references ended with "C" or "CN").


## POSITIVE CONTROLS: INTENDED USE & COMPOSITION

The biotal Yersinia enterocolitica O:9 Positive control is an external quality control for biotal Yersinia enterocolitica O:9 Card. Use this control swabs to check that the sample reagents and the test card are working properly.

biotal Yersinia enterocolitica O:9 Positive control is manufactured with *Yersinia enterocolitica* O:9 extract dried in swab containing a preservative.

## POSITIVE CONTROLS: TESTING PROCEDURE

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- Add 15 drops of diluent and immediately introduce the positive control swab into the testing tube ①.
- Mix the solution by rotating the positive control swab forcefully against the side of the testing tube at least 30 seconds ②. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotal device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S ③.
- Read the result at 10 minutes.

## POSITIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = RED LINE</p>	<p><b>RESULT: CONTROL POSITIVE</b></p> <p>A GREEN line (C) and a RED line (T) appear across the results window during the test performance due to presence of the antigens in the positive control swab.</p>
2		<p><b>INVALID RESULT</b></p> <p>Total absence of the control line (GREEN), regardless the appearance or not of the test line (RED), or only appears the control line.</p> <p>Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

## NEGATIVE CONTROL: INTENDED USE & COMPOSITION

The biotical Negative Control is an external quality control for all immunochromatographic tests manufactured by biotical. Use this control swab to check that the extraction reagents and the test are working properly.

biotical Negative Control is manufactured with negative control buffer dried in swab containing a preservative.

## NEGATIVE CONTROL: TESTING PROCEDURE

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
GREEN line (Control)  
RED line (Test Yersinia enterocolitica O:9)

Sample window (S)

- 1 Add 15 drops of diluent and introduce the negative control swab into the testing tube.
- 2 Rotating 30 seconds and extract the liquid.
- 3 Add 4 drops in the circular window (S).

- Add 15 drops of diluent and immediately introduce the negative control swab into the testing tube 1.
- Mix the solution by rotating the negative control swab forcefully against the side of the testing tube at least 30 seconds 2. Extract as much liquid as possible from the swab, squeezing the sides of the testing tube as the swab is withdrawn. Discard the swab.
- Remove the biotical device from its sealed bag just before using it.
- Use a separate pipette and test for each control. Dispense 4 drops from the testing tube in the circular window marked with the letter S 3.
- Read the result at 10 minutes.

## NEGATIVE CONTROL: INTERPRETATION OF RESULTS

	biotical Test	Interpretation of the results
1	 <p>C = GREEN LINE T = NO LINE</p>	<p><b>RESULT: CONTROL NEGATIVE</b></p> <p>Only the control line (C) appears across the results window during the test performance.</p>
2		<p><b>RESULT: INVALID FOR NEGATIVE CONTROL</b></p> <p>Total absence of the control line, and/or the card shows the test line. Incorrect procedural techniques or deterioration of the reagents are the main reasons for control line failure.</p> <p>Review the procedure and repeat the control with a new test and a new negative control swab.</p>

## LIMITATIONS OF THIS TEST

- The test should be performed no later than 2 hours after opening the bag containing it.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of the test line can depend on whether the antigen concentration is close to the detection limit value of the test or not.
- The test has been validated with human faecal samples. Its use with other class of samples has not been tested so its use is not recommended. Proper sampling is essential to ensure the validity of the test results.
- A positive test result confirms the presence of *Yersinia enterocolitica* O:9 in the sample, nevertheless it must be confirmed together with other laboratory techniques (biochemical and serological methods or by PCR). The infection must be confirmed by a specialist, after evaluating the clinical tests and all the laboratory results.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity:

Cut-off value of biotal Yersinia enterocolitica O:9 Card is  $2 \times 10^7$  CFU/mL.

### Clinical sensitivity and specificity:

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (biotal Yersinia Enterocolitica O:9 Card, biotal) and an agglutination test (Yersinia enterocolitica Agglutination Kit, Progen).

The results were as follows:

		AG test: Yersinia enterocolitica Agglutination Kit		
		+	-	Total
IC Test: Biotal Yersinia enterocolitica O:9 Card	+	2	0	2
	-	0	61	61
	Total	2	61	63

<i>Biotal Yersinia enterocolitica O:9 Card vs Yersinia enterocolitica Agglutination Kit</i>	MIDDLE VALUE	95% CONFIDENCE INTERVAL
SENSITIVITY	>99.0%	15.8% - 100%
SPECIFICITY	>99.0%	94.1% - 100%
PPV	>99.0%	-
NPV	>99.0%	-
ACCURACY	>99.0%	94.3% - 100%

The results showed a high sensitivity and specificity to detect *Yersinia enterocolitica O:9* using biotal Yersinia enterocolitica O:9 Card.

## CROSS REACTIVITY

An evaluation was performed to determine the cross reactivity of biotical Yersinia enterocolitica O:9 Card; no cross reactivity against gastrointestinal pathogens occasionally present in faeces.

*Campylobacter jejuni*

*Campylobacter coli*

*Clostridium difficile*

*Escherichia coli* O157:H7

*Helicobacter pylori*

*Listeria monocytogenes*

*Salmonella enteritidis*

*Salmonella paratyphi*

*Salmonella typhi*

*Salmonella typhimurium*

*Shigella boydii*

*Shigella dysenteriae*

*Shigella flexneri*

*Shigella sonnei*

*Streptococcus aureus*

*Yersinia enterocolitica* O:3

## SYMBOLS FOR IVD COMPONENTS & REAGENTS

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## REFERENCES / BIBLIOGRAFÍA /

---

1. WEYNANTS V., et al. "Detection of Yersinia enterocolitica Serogroup O:3 by a PCR Method". Journal of clinical microbiology, May 1996. p. 1224-1227, Vol. 34, No.5.
2. FREDRIKSSON-AHOMAA M, et al. "Low occurrence of Pathogenic Yersinia enterocolitica in Clinical, Food, and Environmental Samples: a Methodological Problem". Clinical microbiology reviews, Apr. 2003. p. 22-229, Vol. 16, No. 2.
3. ROSNER B., et al. "Epidemiology of reported Yersinia enterocolitica infections in Germany, 2001-2008". BMC Public Health, 2010,10:337.

Sierra de Guadarrama, 1. Ed. Dublín  
28830 San Fernando de Henares. Madrid  
Tel.: + 34 91 677 43 08

[info@biotical.es](mailto:info@biotical.es)  
[www.biotical.es](http://www.biotical.es)

biotical health S.L.U.